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
PVC measured by DSC and TGA

Sample Unplasticized PVC, Trosiplast 3255 injection-molded sheet material

Conditions Measuring cell: DSC ,
TGA
Pan: DSC: Aluminum standard 40 µl, pierced lid
TGA: Alumina 70 µl without lid
Sample preparation: Disk cut off sheet
DSC measurement: Heating from 30 to 120 °C at 10 K/min, quench-cooled in the sample changer
Second run from 30 to 300 °C at 10 K/min
TGA measurement: Heating from 30 to 300 °C at 10 K/min
Atmosphere: Nitrogen, 50 cm³/min

Evaluation

Glass temperature, onset, first run in °C	68.8
Glass temperature, onset, second run in °C	72.5

Decomposition onset DSC in °C	250.0
Decomposition onset TGA in °C	281.0
Decomposition onset DTG in °C	242.0

All onsets are based on the inflectional tangents.

Interpretation

The first DSC curve shows typical discontinuities. They are caused by stress re-lief and enthalpy relaxation. Usually the second run gives glass transitions free of any disturbances.

Above about 200 °C, thermal degradation often begins in the DSC with an exo-thermal peak and always with weight loss due to the formation of HCl, measured by TGA. The HCl could attack the measuring cell. Purging the cell with nitrogen is a countermeasure to prevent damage. If only the beginning of decomposition is to be studied, the temperature program can be interrupted and the cell cooled down. The first decomposition peak is evaluated as onset for comparison purposes.

Part of the exothermal decomposition is caused by the chemical reaction of the HCl with the aluminum pan.

Conclusion

PVC is a polymer that is sensitive to heat. It decomposes with the formation of HCl.

The characteristic features are the glass transition temperature and the temperature at which degradation begins (which is a measure of its stability toward heat).