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Influences of Heating Rate on Decomposition, Metolazone

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Sample Metolazone

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Application Active ingredient (diuretic)

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Pan: Al 40 μ l, with pierced lid with a hole of 1 mm diameter.

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Measurement: Heating from 20 $^{\circ}$ C to 200 $^{\circ}$ C at 2 und 5 K/min. Both measurements are blank curve corrected.

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Figure 3.6 DSC curves of α -Metolazone measured at two different heating rates

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The DSC curves of α -Metolazone measured at two different heating rates are compared in Fig. 3.6. At 20 K/min the melting of the α modification at approx. 227 °C followed by recrystallization in the γ modification can be observed. This melts with accompanying decomposition at approx. 255 °C.

With a heating rate 5 K/min, a much smaller endothermic effect and a depression of the melting point by approx. 7 °C is observed. At 220 °C the hermetically sealed pan bursts (shown as a displacement of the base line).

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Only the endothermic peaks are evaluated. To integrate the peaks the baseline types 'horizontal left' or 'horizontal right' are selected.

升温速率	起始点 1	ΔH	起始点 2	ΔH
Heating rate	Onset 1	J/g	Onset 2	J/g
K/min	$^{\circ}\text{C}$		$^{\circ}\text{C}$	
5	218.9	53.2	-	-
20	227.0	83.1	255.1	31.0

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Substances which decompose at the melting point can only characterized to a limited extent by DSC. In such cases, a more reproducible characterization is often possible by using faster heating rates, since the decomposition reaction is moved to higher temperatures. TGA can also be used in the same way to gain additional information, especially when volatile decomposition products are formed.