

# SOP for Periodic Sensitivity Tests of Moisture Analyzers

**Title:** SOP for Periodic Sensitivity Tests of Moisture Analyzers

Document No.:

Author: METTLER TOLEDO Laboratory & Weighing Technologies

Number of Pages: 4

## Controls

Replacement Document: N/A

Reason for Revision: N/A

Release Date: 11. Oct. 2011

## Release for Routine Operation

Reviewed by:

Date:

Signature:

**Approved by:**

Date:

Signature:

## General

To test the sensitivity of the weighing unit of a moisture analyzer, a large test weight (between 50% and 100% of capacity) is recommended.

### Basic Rules for Handling Moisture Analyzers

- Before using a moisture analyzer, make sure the instrument was left on power for a sufficient period of time (mentioned in the user manual).
- Make sure the moisture analyzer is leveled.
- Minimize environmental influences, e.g. air draft, vibrations or direct sunlight.

### Basic Rules for Handling Weights

#### Important to know

- Only an external weight with calibration certificate can make the weighing unit of a moisture analyzer a traceable piece of equipment.
- Weights should always be placed gently on the sample pan holder and put back immediately in their storage place after use. They may be parked on a clean piece of paper.
- Weights (since they are also part of measuring equipment) need to be re-calibrated at specified intervals (ISO 9001).
- Any damage, which might have affected the value of the weight, should trigger an immediate re-calibration. METTLER TOLEDO's calibration services will give advice on this.

#### How to Store Weights

- Weights should be stored in their original box.
- Weights should be stored at room temperature, since temperature differences between weights and their surrounding could lead to measurement errors.
- Weights that have not been stored at the same temperature need acclimatization before use.

#### How to Move Weights

- Weights should only be handled with appropriate tools such as tweezers, forks, handles or gloves (see METTLER TOLEDO's accessories for weights).
- These tools should be exclusively used for transferring weights, due to possible contamination.

# Sensitivity Test

## Preparation

- Before the test is performed, the test weights must be acclimatized to the ambient temperature of the moisture analyzer.
- When the test is performed, the heating unit of the moisture analyzer should be at room temperature.

## Test Procedure HX and HS Models

- Start test mode and follow the instructions.

## Test Procedure other Moisture Analyzers

- Remove the sample pan.
- Tare the balance.
- Place the test weight(s) on the sample pan holder.
- Read the stable value from the display and note it.

## Evaluation

- Evaluate whether the noted value exceeds the defined "Warning Limit". <sup>1)</sup>
- Evaluate whether the noted value exceeds the defined "Control Limit". <sup>2)</sup>

## Deviation

### **Warning Limit** <sup>1)</sup> (where defined)

- If the warning limit is exceeded, repeat the test.
- If the warning limit is exceeded again, report that the warning limit was not met to the laboratory supervisor or the person responsible for the moisture analyzer.
- Level the moisture analyzer, perform adjustment with built-in (if applicable) or external weight and repeat the test.
- If the warning limit is still exceeded, report the problem to the laboratory supervisor or the person responsible of the moisture analyzer. Optionally, contact METTLER TOLEDO's service organization for advice.

### **Control Limit** <sup>2)</sup>

- If the control limit is exceeded, report the problem to the laboratory supervisor or the person responsible of the moisture analyzer.
- Mark the moisture analyzer as "out of control limits".
- Contact METTLER TOLEDO service organization for advice.

---

<sup>1)</sup> – Value within the warning limit: no action is necessary.  
– Value between the warning and control limit are within the tolerance but must be kept under surveillance.

<sup>2)</sup> – Value within the control limit, see <sup>1)</sup>  
– Value beyond the control limit show that weighing process is no longer under control and immediate action is therefore required.

## Recommended Test Weights, Warning and Control Limits

Moisture Analyzer	HX/HS	HR/HG	HB/MJ
Test weight	100 g	50 g	20 g
Warning Limit	± 25 mg	± 12.5 mg	± 5 mg
Control Limit	± 50 mg	± 25 mg	± 10 mg

Recommendation is based on 0.1% weighing accuracy and safety factor of 2. The absolute tolerances depend on the mass of the test weight.

Various weighing parameters can contribute to the measurement uncertainty of a weighing result. With the exception of the repeatability test, all control limits are set to 1/2 of the weighing accuracy in order to have a security reserve accounting for any other influences such as eccentricity and linearity on the result (control limit = test weight \* weighing accuracy / 2)

[www.mt.com/moisture](http://www.mt.com/moisture)

For more information

### **Mettler-Toledo AG**

Laboratory & Weighing Technologies

CH-8606 Greifensee

Tel. +41 44 944 22 11

Fax. +41 44 944 31 70

Subject to technical changes

© 11/2011 Mettler-Toledo AG

Printed in Switzerland 30031488

Global MarCom Switzerland