

METTLER TOLEDO

Table of Contents

| | | |
|----------|--|-----------|
| 1 | Introduction | 7 |
| 1.1 | Document purpose | 7 |
| 1.2 | Further documents and information | 7 |
| 1.3 | Explanation of conventions and symbols used | 7 |
| 1.4 | Acronyms and abbreviations | 8 |
| 1.5 | Product range | 8 |
| 1.5.1 | MX analytical balances | 8 |
| 1.5.2 | MX precision balances | 9 |
| 1.5.3 | MX large balances | 9 |
| 2 | Safety Information | 10 |
| 2.1 | Definition of signal words and warning symbols | 10 |
| 2.2 | Product-specific safety notes | 10 |
| 3 | Design and Function | 12 |
| 3.1 | Overview analytical balances | 12 |
| 3.2 | Overview precision balances, small | 13 |
| 3.2.1 | Balances with draft shield | 13 |
| 3.2.2 | Balances without draft shield | 14 |
| 3.3 | Overview precision balances, large | 15 |
| 3.4 | Overview terminal | 15 |
| 3.5 | Overview interface connections | 16 |
| 3.6 | Components description | 16 |
| 3.6.1 | Draft shield | 16 |
| 3.6.2 | Weighing pan | 17 |
| 3.6.3 | Draft-protection element | 17 |
| 3.6.4 | Drip tray | 17 |
| 3.6.5 | Door handle | 17 |
| 3.6.6 | ErgoDoor handle | 18 |
| 3.6.7 | Leveling feet | 18 |
| 3.6.8 | Terminal | 18 |
| 3.6.9 | QuickLock for draft shield | 18 |
| 3.6.10 | QuickLock for top door and front panel | 19 |
| 3.6.11 | QuickLock for side door | 19 |
| 3.6.12 | Release button for back panel | 19 |
| 3.7 | Overview type label | 19 |
| 3.8 | User interface | 20 |
| 3.8.1 | Main sections at a glance | 20 |
| 3.8.2 | Main weighing screen | 20 |
| 3.8.3 | Balance menu | 21 |
| 3.8.4 | Main configuration | 22 |
| 3.8.5 | Applications | 22 |
| 3.8.6 | Icons and symbols | 23 |
| 3.8.6.1 | System status icons | 23 |
| 3.8.6.2 | Weighing status icons | 23 |
| 3.8.6.3 | Process status icons | 23 |
| 4 | Installation and Putting into Operation | 25 |
| 4.1 | Selecting the location | 25 |
| 4.2 | Unpacking the balance | 25 |
| 4.3 | Installation | 26 |
| 4.3.1 | Balances with draft shield | 26 |
| 4.3.2 | Balances without draft shield | 29 |
| 4.3.2.1 | Assembling balances 1 mg | 29 |
| 4.3.2.2 | Assembling balances 0.01 g | 29 |

| | | |
|-----------|---|-----------|
| 4.3.2.3 | Assembling balances 0.1 g | 29 |
| 4.3.3 | Balances, large | 30 |
| 4.4 | Putting into operation..... | 30 |
| 4.4.1 | Connecting the balance | 30 |
| 4.4.2 | Switching on the balance..... | 31 |
| 4.4.3 | Leveling the balance | 31 |
| 4.4.4 | Performing an internal adjustment | 31 |
| 4.4.5 | Entering / Exiting standby mode | 32 |
| 4.4.6 | Entering / Exiting power-saving mode | 32 |
| 4.4.7 | Switching off the balance | 32 |
| 4.5 | Performing a simple weighing | 32 |
| 4.5.1 | Opening and closing draft shield doors | 32 |
| 4.5.2 | Zeroing the balance..... | 33 |
| 4.5.3 | Taring the balance | 33 |
| 4.5.4 | Performing a weighing | 34 |
| 4.6 | Transporting, packing, and storing..... | 34 |
| 4.6.1 | Transporting the balance over short distances | 34 |
| 4.6.2 | Transporting the balance over long distances..... | 34 |
| 4.6.3 | Packing and storing | 35 |
| 4.7 | Weighing below the balance..... | 35 |
| 5 | Operation | 36 |
| 5.1 | Touchscreen | 36 |
| 5.1.1 | Selecting or activating an item | 36 |
| 5.1.2 | Scrolling | 36 |
| 5.1.3 | Opening the fly-in panel..... | 36 |
| 5.1.4 | Entering characters and numbers..... | 36 |
| 5.1.5 | Changing values..... | 37 |
| 5.1.6 | Sliding..... | 37 |
| 5.2 | General balance settings | 37 |
| 5.2.1 | Date / Time / Language | 37 |
| 5.2.2 | Screen / StatusLight / Sound..... | 38 |
| 5.2.3 | Standby, Power-saving mode..... | 38 |
| 5.2.4 | Weighing / Quality | 39 |
| 5.2.4.1 | Warnings and reminders | 39 |
| 5.2.4.2 | Weighing profiles..... | 40 |
| 5.2.4.2.1 | Environment..... | 40 |
| 5.2.4.2.2 | Weighing mode..... | 40 |
| 5.2.4.2.3 | Value release | 40 |
| 5.2.4.3 | Test weights | 40 |
| 5.2.4.3.1 | Defining an individual test weight | 40 |
| 5.2.4.3.2 | Defining a combined test weight | 41 |
| 5.2.4.3.3 | Deleting a test weight | 41 |
| 5.3 | Weighing applications | 42 |
| 5.3.1 | Weighing applications overview | 42 |
| 5.3.2 | General settings for weighing applications | 42 |
| 5.3.2.1 | Defining target weight and tolerances | 42 |
| 5.3.2.2 | Defining a sample ID..... | 43 |
| 5.3.2.3 | Configuring a weighing application | 43 |
| 5.3.2.4 | Configuring a weighing series | 44 |
| 5.3.2.5 | Using automated functions..... | 44 |
| 5.3.2.6 | Configuring a report | 45 |
| 5.3.3 | Application "Weighing" | 45 |
| 5.3.4 | Application "Counting" | 46 |
| 5.3.5 | Application "Check weighing" | 46 |
| 5.3.6 | Application "Dynamic weighing" | 47 |
| 5.3.7 | Application "Formulation" | 48 |
| 5.3.8 | Application "Totaling" | 49 |

| | | |
|----------|--|----|
| 5.3.9 | Application "Back weighing" | 50 |
| 5.3.10 | Application "Density" | 50 |
| 5.3.11 | Application "Differential weighing" | 51 |
| 5.4 | Adjustments..... | 52 |
| 5.4.1 | Adjustment strategy | 52 |
| 5.4.2 | Editing an adjustment..... | 52 |
| 5.4.3 | Performing an internal adjustment | 53 |
| 5.4.4 | Performing an external adjustment..... | 53 |
| 5.5 | Tests..... | 54 |
| 5.5.1 | Editing a test | 54 |
| 5.5.2 | Performing a test | 54 |
| 5.5.2.1 | Sensitivity test..... | 54 |
| 5.5.2.2 | Repeatability test..... | 55 |
| 5.5.2.3 | Eccentricity test..... | 55 |
| 5.6 | Interfaces | 56 |
| 5.6.1 | Ethernet | 56 |
| 5.6.2 | Bluetooth | 56 |
| 5.7 | Devices / Printers | 57 |
| 5.7.1 | Printer | 57 |
| 5.7.1.1 | Installing a USB printer | 57 |
| 5.7.1.2 | Installing an RS232 printer | 58 |
| 5.7.1.3 | Installing a printer through Bluetooth | 58 |
| 5.7.2 | Barcode reader | 60 |
| 5.7.2.1 | Scanning a sample ID with a barcode reader | 60 |
| 5.7.3 | Foot switch | 60 |
| 5.7.4 | Keyboard | 61 |
| 5.7.5 | Adding and deleting a device | 61 |
| 5.7.6 | Editing the settings of a device | 62 |
| 5.8 | Services | 62 |
| 5.8.1 | Setting up services | 62 |
| 5.8.1.1 | MT-SICS service..... | 62 |
| 5.8.1.2 | EasyDirect Balance | 63 |
| 5.8.1.3 | Drop to cursor | 63 |
| 5.8.1.4 | File server | 64 |
| 5.8.2 | Transferring data to services..... | 64 |
| 5.8.2.1 | Transferring data: MT-SICS service | 64 |
| 5.8.2.2 | Transferring data: EasyDirect Balance | 65 |
| 5.8.2.3 | Transferring data: Drop to cursor | 66 |
| 5.8.2.4 | Transferring data: File server | 66 |
| 5.9 | Publishing..... | 67 |
| 5.9.1 | Printing data | 67 |
| 5.9.1.1 | Printing results manually through USB..... | 67 |
| 5.9.1.2 | Printing results automatically through Bluetooth | 68 |
| 5.9.2 | Exporting data to a USB storage device | 68 |
| 5.9.3 | Publishing options | 69 |
| 5.9.4 | Indicators for weighing results | 69 |
| 5.10 | User management | 69 |
| 5.10.1 | Activating / deactivating the user management | 69 |
| 5.10.2 | Managing users and user groups | 70 |
| 5.10.2.1 | Automatic logout..... | 70 |
| 5.10.2.2 | Creating a new user | 70 |
| 5.10.2.3 | Setting the user language..... | 70 |
| 5.10.2.4 | Deleting a user | 71 |
| 5.10.2.5 | Managing groups | 71 |
| 5.11 | Password protection | 71 |
| 5.11.1 | Logging in and logging out..... | 71 |
| 5.11.2 | Changing a password | 72 |

| | | |
|-----------|--|-----------|
| 5.11.3 | Resetting a password | 72 |
| 5.11.4 | Blocking and unblocking the balance | 72 |
| 6 | Software Description | 73 |
| 6.1 | Balance menu settings | 73 |
| 6.1.1 | Leveling aid | 73 |
| 6.1.2 | History | 73 |
| 6.1.2.1 | Adjustments history | 73 |
| 6.1.2.2 | Tests history | 74 |
| 6.1.2.3 | Service history | 74 |
| 6.1.2.4 | Activity log | 74 |
| 6.1.2.5 | Software update history | 74 |
| 6.1.2.6 | Error log | 75 |
| 6.1.3 | Information | 75 |
| 6.1.3.1 | Balance information | 75 |
| 6.1.3.2 | Service and support information | 75 |
| 6.1.4 | User management | 75 |
| 6.1.4.1 | User management – General | 76 |
| 6.1.4.2 | User management – Users | 76 |
| 6.1.4.3 | User management – Groups | 76 |
| 6.1.5 | Settings | 77 |
| 6.1.5.1 | Settings: Balance | 77 |
| 6.1.5.1.1 | Settings: Weighing / Quality | 77 |
| 6.1.5.1.2 | Settings: Publishing | 79 |
| 6.1.5.1.3 | Settings: Date / Time / Language | 80 |
| 6.1.5.1.4 | Settings: Screen / StatusLight / Sound | 81 |
| 6.1.5.1.5 | Settings: General | 82 |
| 6.1.5.2 | Settings: Interfaces | 82 |
| 6.1.5.2.1 | Settings: Ethernet | 82 |
| 6.1.5.2.2 | Settings: Bluetooth | 83 |
| 6.1.5.3 | Settings: Devices / Printers | 83 |
| 6.1.5.4 | Settings: Services | 83 |
| 6.1.5.4.1 | Settings: MT-SICS service | 84 |
| 6.1.5.4.2 | Settings: EasyDirect Balance | 84 |
| 6.1.5.4.3 | Settings: Drop to cursor | 84 |
| 6.1.5.4.4 | Settings: File server | 85 |
| 6.1.6 | Maintenance | 85 |
| 6.2 | Weighing applications settings | 85 |
| 6.2.1 | Settings: application "Weighing" | 85 |
| 6.2.1.1 | Main | 86 |
| 6.2.1.2 | ID format | 86 |
| 6.2.1.3 | Weighing | 87 |
| 6.2.1.4 | Automation | 88 |
| 6.2.1.5 | Report | 88 |
| 6.2.2 | Settings: application "Counting" | 88 |
| 6.2.2.1 | Main | 89 |
| 6.2.2.2 | ID format | 89 |
| 6.2.2.3 | Weighing | 89 |
| 6.2.2.4 | Automation | 90 |
| 6.2.2.5 | Report | 90 |
| 6.2.3 | Settings: application "Check weighing" | 91 |
| 6.2.3.1 | Main | 91 |
| 6.2.3.2 | ID format | 91 |
| 6.2.3.3 | Weighing | 92 |
| 6.2.3.4 | Automation | 93 |
| 6.2.3.5 | Report | 93 |
| 6.2.4 | Settings: application "Dynamic weighing" | 93 |
| 6.2.4.1 | Main | 94 |

| | | |
|----------|---|------------|
| 6.2.4.2 | ID format | 94 |
| 6.2.4.3 | Weighing | 94 |
| 6.2.4.4 | Automation | 95 |
| 6.2.4.5 | Report | 95 |
| 6.2.5 | Settings: application "Formulation" | 95 |
| 6.2.5.1 | Main | 96 |
| 6.2.5.2 | ID format | 96 |
| 6.2.5.3 | Weighing | 96 |
| 6.2.5.4 | Automation | 97 |
| 6.2.5.5 | Report | 97 |
| 6.2.6 | Settings: application "Totaling" | 97 |
| 6.2.6.1 | Main | 97 |
| 6.2.6.2 | ID format | 98 |
| 6.2.6.3 | Weighing | 98 |
| 6.2.6.4 | Automation | 98 |
| 6.2.6.5 | Report | 99 |
| 6.2.7 | Settings: application "Back weighing" | 99 |
| 6.2.7.1 | Main | 99 |
| 6.2.7.2 | ID format | 100 |
| 6.2.7.3 | Weighing | 100 |
| 6.2.7.4 | Automation | 101 |
| 6.2.7.5 | Report | 101 |
| 6.2.8 | Settings: application "Density" | 101 |
| 6.2.8.1 | Main | 102 |
| 6.2.8.2 | ID format | 102 |
| 6.2.8.3 | Weighing | 103 |
| 6.2.8.4 | Report | 103 |
| 6.2.9 | Settings: application "Differential weighing" | 103 |
| 6.2.9.1 | Main | 104 |
| 6.2.9.2 | ID format | 104 |
| 6.2.9.3 | Weighing | 105 |
| 6.2.9.4 | Automation | 105 |
| 6.2.9.5 | Report | 106 |
| 6.3 | Adjustment settings | 106 |
| 6.3.1 | Settings: Adjustment strategy | 106 |
| 6.3.2 | Settings: Internal adjustment | 106 |
| 6.3.3 | Settings: External adjustment | 108 |
| 6.4 | Test settings | 108 |
| 6.4.1 | Settings: Sensitivity test | 108 |
| 6.4.2 | Settings: Repeatability test | 110 |
| 6.4.3 | Settings: Eccentricity test | 111 |
| <hr/> | | |
| 7 | Maintenance | 113 |
| 7.1 | Maintenance tasks | 113 |
| 7.2 | Cleaning | 113 |
| 7.2.1 | Disassembling for cleaning | 113 |
| 7.2.1.1 | Balances with draft shield | 113 |
| 7.2.1.2 | Balances without draft shield | 116 |
| 7.2.1.3 | Balances, large | 116 |
| 7.2.2 | Cleaning agents | 117 |
| 7.2.3 | Cleaning the balance | 117 |
| 7.2.4 | Putting into operation after cleaning | 118 |
| 7.3 | Service | 118 |
| 7.4 | Software update | 118 |
| 7.4.1 | Updating the software | 119 |
| 7.4.2 | Putting into operation after software update | 119 |
| 7.5 | Resetting the balance | 119 |

| | | |
|-----------|--|------------|
| 8 | Troubleshooting | 120 |
| 8.1 | Error messages..... | 120 |
| 8.2 | Error symptoms | 121 |
| 8.3 | Saving a support file..... | 124 |
| 8.4 | Putting into operation after fixing an error | 124 |
| 9 | Technical Data | 125 |
| 9.1 | General data..... | 125 |
| 9.2 | Materials | 126 |
| 9.3 | Explanatory notes for the METTLER TOLEDO AC/DC adapter..... | 126 |
| 9.4 | Model-specific data | 128 |
| 9.4.1 | Analytical balances, readability 0.01 mg or 0.1 mg | 128 |
| 9.4.2 | Precision balances, readability 1 mg | 130 |
| 9.4.3 | Precision balances, readability 0.01 g or 0.1 g..... | 132 |
| 9.4.4 | Precision balances, large..... | 135 |
| 9.5 | Dimensions | 136 |
| 9.5.1 | Analytical balances, readability 0.01 mg | 136 |
| 9.5.2 | MX analytical balances, readability 0.1 mg | 137 |
| 9.5.3 | MX precision balances, readability 1 mg, with draft shield | 138 |
| 9.5.4 | MX precision balances, readability 1 mg, without draft shield..... | 139 |
| 9.5.5 | MX precision balances, readability 0.01 g | 140 |
| 9.5.6 | MX precision balances, readability 0.1 g | 141 |
| 9.5.7 | MX large balances, readability 0.1 g / 1 g | 142 |
| 10 | Accessories and Spare Parts | 143 |
| 10.1 | Accessories | 143 |
| 10.2 | Spare parts..... | 149 |
| 10.2.1 | MX analytical balances, readability 0.01 mg | 149 |
| 10.2.2 | MX analytical balances, readability 0.1 mg | 150 |
| 10.2.3 | MX precision balances, readability 1 mg, with draft shield | 151 |
| 10.2.4 | MX precision balances, readability 1 mg, without draft shield..... | 153 |
| 10.2.5 | MX precision balances, readability 0.01 g | 155 |
| 10.2.6 | MX precision balances, readability 0.1 g | 156 |
| 10.2.7 | MX large balances, readability 0.1 g / 1 g | 157 |
| 10.2.8 | AC/DC adapters | 158 |
| 10.2.8.1 | AC/DC adapter | 158 |
| 10.2.8.2 | AC/DC adapter, universal | 159 |
| 10.2.9 | Packaging | 160 |
| 10.2.9.1 | Balances with draft shield | 160 |
| 10.2.9.2 | Balances without draft shield..... | 161 |
| 10.2.9.3 | Balances, large | 161 |
| 11 | Disposal | 162 |
| 12 | Compliance Information | 163 |
| | Index | 165 |

1 Introduction

Thank you for choosing a METTLER TOLEDO balance. The balance combines high performance with ease of use.

This document is based on the software version V 1.2.

EULA

The software in this product is licensed under the METTLER TOLEDO End User License Agreement (EULA) for Software.

When using this product you agree to the terms of the EULA.

▶ www.mt.com/EULA

1.1 Document purpose

This Reference Manual provides detailed instructions on how to use the instrument.

1.2 Further documents and information

This document is available in other languages online.



▶ www.mt.com/MX-RM

Product page:

▶ www.mt.com/MX-balances

Instructions for cleaning a balance, "8 Steps to a Clean Balance":

▶ www.mt.com/lab-cleaning-guide

Search for software:

▶ www.mt.com/labweighing-software-download

Search for documents:

▶ www.mt.com/library

For further questions, please contact your authorized METTLER TOLEDO dealer or service representative.

▶ www.mt.com/contact

1.3 Explanation of conventions and symbols used

Conventions and symbols

Key and/or button designations and display texts are shown in graphic or bold text, e.g.,  **Publish**.

 **Note**

For useful information about the product.



Refers to an external document.

Elements of instructions

In this manual, step-by-step instructions are presented as follows. The action steps are numbered and can contain prerequisites, intermediate results and results, as shown in the example. Sequences with less than two steps are not numbered.

- Prerequisites that must be fulfilled before the individual steps can be executed.

1 Step 1

➔ Intermediate result

2 Step 2


➔ Result


1.4 Acronyms and abbreviations

| Original term | Explanation |
|---------------|---|
| AC | Alternating Current |
| ASTM | American Society for Testing and Materials |
| DC | Direct Current |
| EMC | Electromagnetic Compatibility |
| FCC | Federal Communications Commission |
| GWP | Good Weighing Practice |
| ID | Identification |
| IP | Ingress Protection |
| LAN | Local Area Network |
| LED | Light-Emitting Diode |
| LPS | Limited Power Source |
| MAC | Media Access Control |
| MT-SICS | METTLER TOLEDO Standard Interface Command Set |
| NA | Not Applicable |
| OIML | Organisation Internationale de Métrologie Légale (International Organization of Legal Metrology) |
| RM | Reference Manual |
| SOP | Standard Operating Procedure |
| TDNR | Type Definition Number |
| UM | User Manual |
| USB | Universal Serial Bus |
| USP | United States Pharmacopeia |




1.5 Product range

1.5.1 MX analytical balances


| Balance | Models designation |
|---|---|
|  | Readability: 0.01 mg <ul style="list-style-type: none">• MX105• MX105DU• MX205DU |

| Balance | Models designation |
|---|--|
|  | Readability: 0.1 mg <ul style="list-style-type: none"> • MX104 • MX204 • MX304 |

1.5.2 MX precision balances

| Balance | Models designation |
|---|---|
|  | Readability: 1 mg With draft shield: <ul style="list-style-type: none"> • MX303 • MX603 • MX1203 Without draft shield: <ul style="list-style-type: none"> • MX303N • MX603N • MX1203N |
|  | Readability: 0.01 g <ul style="list-style-type: none"> • MX2002 • MX4002 • MX6002 • MX6002DR • MX12002 |
|  | Readability: 0.1 g <ul style="list-style-type: none"> • MX6001 • MX8001 |

1.5.3 MX large balances

| Balance | Models designation |
|---|--|
|  | Readability: 0.1 g / 1 g <ul style="list-style-type: none"> • MX12001L • MX16001L • MX32001L • MX32000L |

2 Safety Information

Two documents named "User Manual" and "Reference Manual" are available for this instrument.

- The User Manual is available online in various languages.
- A printed version of the User Manual is delivered with the instrument.
- The Reference Manual is available online. This manual contains a full description of the instrument and its use.
- Keep both documents for future reference.
- Include both documents if you transfer the instrument to other parties.

Only use the instrument according to the User Manual and the Reference Manual. If you do not use the instrument according to these documents or if the instrument is modified, the safety of the instrument may be impaired and Mettler-Toledo GmbH assumes no liability.

2.1 Definition of signal words and warning symbols

Safety notes contain important information on safety issues. Ignoring the safety notes may lead to personal injury, damage to the instrument, malfunctions and false results. Safety notes are marked with the following signal words and warning symbols:

Signal words

| | |
|----------------|---|
| DANGER | A hazardous situation with high risk, resulting in death or severe injury if not avoided. |
| WARNING | A hazardous situation with medium risk, possibly resulting in death or severe injury if not avoided. |
| CAUTION | A hazardous situation with low risk, resulting in minor or moderate injury if not avoided. |
| NOTICE | A hazardous situation with low risk, resulting in damage to the instrument, other material damage, malfunctions and erroneous results, or loss of data. |

Warning symbols



General hazard



Notice

2.2 Product-specific safety notes

Intended use

This instrument is designed to be used by trained staff. The instrument is intended for weighing purposes.

Any other type of use and operation beyond the limits of use stated by Mettler-Toledo GmbH without consent from Mettler-Toledo GmbH is considered as not intended.

Responsibilities of the instrument owner

The instrument owner is the person holding the legal title to the instrument and who uses the instrument or authorizes any person to use it, or the person who is deemed by law to be the operator of the instrument. The instrument owner is responsible for the safety of all users of the instrument and third parties.

Mettler-Toledo GmbH assumes that the instrument owner trains users to safely use the instrument in their workplace and deal with potential hazards. Mettler-Toledo GmbH assumes that the instrument owner provides the necessary protective gear.

Safety notes



WARNING

Death or serious injury due to electric shock

Contact with parts that carry a live current can lead to death or injury.

- 1 Only use the METTLER TOLEDO power cable and AC/DC adapter designed for your instrument.
- 2 Connect the power cable to a grounded power outlet.
- 3 Keep all electrical cables and connections away from liquids and moisture.
- 4 Check the cables and the power plug for damage and replace them if damaged.



NOTICE


Damage to the instrument or malfunction due to the use of unsuitable parts

- Only use parts from METTLER TOLEDO that are intended to be used with your instrument.


3 Design and Function

3.1 Overview analytical balances

0.01 mg



0.1 mg




| | | | |
|----------|----------------------------|-----------|--------------------------------------|
| 1 | StatusLight | 9 | QuickLock, top door/panel |
| 2 | Terminal | 10 | Top door, draft shield |
| 3 | Draft-protection element | 11 | QuickLock, side door |
| 4 | Weighing pan | 12 | Side door, draft shield (right/left) |
| 5 | Drip tray | 13 | Handle, side door |
| 6 | Back panel *, draft shield | 14 | ErgoDoor handle |
| 7 | Front panel, draft shield | 15 | Leveling feet |
| 8 | Handle, top door | | |

* On balances with a readability of 0.01 mg, the back panel has a backlight.

3.2 Overview precision balances, small

3.2.1 Balances with draft shield


1 mg




| | | | |
|----------|---------------------------|-----------|--------------------------------------|
| 1 | StatusLight | 8 | QuickLock, top door/panel |
| 2 | Terminal | 9 | Top door, draft shield |
| 3 | SmartPan weighing pan | 10 | QuickLock, side door |
| 4 | Drip tray | 11 | Side door, draft shield (right/left) |
| 5 | Back panel, draft shield | 12 | Handle, side door |
| 6 | Front panel, draft shield | 13 | ErgoDoor handle |
| 7 | Handle, top door | 14 | Leveling feet |

3.2.2 Balances without draft shield


1 mg



0.01 g




0.1 g



| | | | |
|----------|--------------------------|----------|----------------------|
| 1 | StatusLight | 6 | Weighing pan support |
| 2 | Terminal | 7 | Weighing pan |
| 3 | Drip tray | 8 | Safety feet |
| 4 | Weighing pan support cap | 9 | Leveling feet |
| 5 | SmartPan weighing pan | | |


3.3 Overview precision balances, large







0.1 g / 1 g




| | | | |
|---|--------------------------|---|---------------|
| 1 | StatusLight | 4 | Weighing pan |
| 2 | Terminal | 5 | Leveling feet |
| 3 | Weighing pan support cap | | |

3.4 Overview terminal



| | Name | Description |
|---|------------------------------------|--|
|  | Standby / Power-saving mode | By tapping  , the balance goes into standby mode. By tapping and holding  , the balance goes into power-saving mode. To switch the balance completely off, it must be unplugged from the power supply. Note Do not disconnect the balance from the power supply unless the balance is not used for an extended period of time. After switching on the instrument, it must warm up before giving accurate results. |
|  | Tare | Tares the balance. This function is used when the weighing process involves containers. After taring the balance, the screen shows <i>Net</i> which indicates that all displayed weights are net. |
|  | Zero | Zeroes the balance. The balance must always be zeroed before starting the weighing process. After zeroing, the balance sets a new zero point. |
|  | Home | To return from any menu level to the main weighing screen. |

3.5 Overview interface connections




| | | | |
|----------|---------------------------|----------|--------------------------|
| 1 | Slot for anti-theft cable | 4 | USB-C port |
| 2 | Ethernet port (LAN) | 5 | Socket for AC/DC adapter |
| 3 | USB-A port | | |

3.6 Components description


3.6.1 Draft shield

The draft shield protects the weighing area against environmental impacts like drafts or moisture. The side doors and the top door can be opened manually.




3.6.2 Weighing pan

The weighing pan is the load receptor that serves to accommodate the weighing item.




3.6.3 Draft-protection element

The draft-protection element protects the weighing pan against draft. This element is only available for balances with a readability of 0.01 mg.




3.6.4 Drip tray

The drip tray is positioned below the weighing pan. The primary purpose of the drip tray is to ensure quick cleaning of the balance.




3.6.5 Door handle

The door handles are mounted on the draft shield doors. The handles are used to manually open the side doors and the top door of the draft shield.




3.6.6 ErgoDoor handle

The ErgoDoor handle is mounted on the platform. The ErgoDoor handle can be engaged with the handle of the side door. This allows to customize opening/closing of the side doors according to your needs.




3.6.7 Leveling feet

The balance stands on height-adjustable feet. These feet are used to level the balance.




3.6.8 Terminal

The balance terminal has a 7-inch touch-sensitive display. A StatusLight LED strip on the front side of the terminal indicates the current status of the balance. The terminal is protected by a replaceable cover.




3.6.9 QuickLock for draft shield

The QuickLock for the draft shield is used to secure the draft shield to the platform.




3.6.10 QuickLock for top door and front panel

Depending on the position, the QuickLock is used to lock/unlock the top door and the front panel of the draft shield.




3.6.11 QuickLock for side door

The QuickLock is used to lock/unlock the side door of the draft shield.



3.6.12 Release button for back panel

The release button is used to lock/unlock the back panel of the draft shield. This feature is only available for the draft shield of balances with a readability of 0.1 mg and 1 mg.



3.7 Overview type label

The information on the type label helps to identify the balance.




| | | | |
|---|---------------------|---|-----------------------|
| 1 | Balance model | 5 | Manufacturer |
| 2 | Year of manufacture | 6 | Balance serial number |
| 3 | Maximum capacity | 7 | Power consumption |
| 4 | Readability | | |

3.8 User interface

3.8.1 Main sections at a glance

The main weighing screen (1) is the central navigation point where all the menus and settings can be found. The sections **Balance menu** (2), **Main configuration** (3), and the applications section (4) open when tapping the corresponding icon or tab.



See also

- [Main weighing screen](#) ▶ Page 20
- [Balance menu](#) ▶ Page 21
- [Main configuration](#) ▶ Page 22
- [Applications](#) ▶ Page 22


3.8.2 Main weighing screen



| | Name | Description |
|---|------------------|--|
| 1 | Weighing results | Shows the results of the current weighing process. |


| | Name | Description |
|----|---------------------------|---|
| 2 | Level indicator | Indicates if the balance is leveled (green) or not (red). |
| 3 | Applications | Accesses available applications: Weighing, Adjustments, Tests. |
| 4 | Readability and capacity | Shows the readability and the capacity of the balance. |
| 5 | Additional information | Shows additional information about the current activity. Example: current weighing value in another unit |
| 6 | Information and warnings | Shows current information, warnings, and error messages. |
| 7 | Functions area | Shows the active functions according to the settings of the current weighing application. |
| 8 | Publish button | Publishes the results according to the settings of the current weighing application. Depending on the selected weighing application, the button can have different functions. |
| 9 | Action bar | Contains actions referring to the current weighing application. |
| 10 | Balance menu | Accesses the balance properties. |
| 11 | SmartTrac | Used as a weighing aid to define a target weight with upper and lower tolerances. |
| 12 | Main configuration | Accesses the configuration options for the current weighing application. |

3.8.3 Balance menu




| | Name | Description |
|---|------------------------|---|
| 1 | Leveling aid | Opens the leveling dialog. |
| 2 | History | Opens the history dialog. |
| 3 | Information | Opens the balance information. |
| 4 | User management | Opens the user management dialog (only visible if activated). |
| 5 | Settings | Opens the settings dialog. |
| 6 | Maintenance | Opens the maintenance dialog. |
| 7 | Logout, Block | Functions related to the menu topic User management . |

3.8.4 Main configuration



| | Name | Description |
|---|-------------------|-------------------------------------|
| 1 | Main | Opens the main configuration. |
| 2 | ID format | Opens the sample ID configuration. |
| 3 | Weighing | Opens the weighing configuration. |
| 4 | Automation | Opens the automation configuration. |
| 5 | Report | Opens the report configuration. |

3.8.5 Applications








| | Name | Description |
|---|--------------------|---|
| 1 | Weighing | Contains available weighing applications. |
| 2 | Adjustments | Contains available adjustments. |
| 3 | Tests | Contains available tests. |




3.8.6 Icons and symbols

3.8.6.1 System status icons





System messages appear due to a user action, a user input, or a system process. When tapping the icon, the corresponding system message is displayed.



| Icon | Name | Description |
|---|--------------------|--|
|  | Levelled | Indicates that the balance is correctly levelled. |
|  | Not levelled | Indicates that the balance is not levelled. |
|  | Information | Provides information related to the current action or process. |
|  | Warning | Provides information about an issue that requires attention. |
|  | Error | Provides information about a failed action or process. |

3.8.6.2 Weighing status icons

| Icon | Name | Description |
|---|--------------------------|--|
|  | Stability indicator | Indicates that the weighing process is ongoing. The weighing result is not yet stable. |
| Net | Net indicator | Appears when pressing the tare key, after the tare weight has been subtracted. |
|  | Calculated value | The current weight value is calculated. This symbol also appears when the balance is tared using the function Preset tare . |
|  | Minimum weight violation | The current weight value is smaller than the defined minimum weight. The weight must be larger than the minimum weight. |

3.8.6.3 Process status icons

| Icon | Name | Description |
|---|-----------------|--|
|  | Start | Starts the process. |
|  | Pause | Pauses the process |
|  | Continue | Continues a paused process. |
|  | Add | Adds the displayed result to a measurement series. |

| Icon | Name | Description |
|---|-----------------|------------------------|
|  | Complete | Completes the process. |
|  | Stop | Stops the process. |

4 Installation and Putting into Operation

4.1 Selecting the location

A balance is a sensitive precision instrument. The location where it is placed will have a profound effect on the accuracy of the weighing results.

Requirements of the location

Place indoors on stable table


Ensure sufficient spacing

Level the instrument


Provide adequate lighting




Avoid direct sunlight





Avoid vibrations



Avoid strong drafts



Avoid temperature fluctuations



Take into account the environmental conditions. See "Technical Data".

Sufficient spacing for balances: > 15 cm all around the instrument


4.2 Unpacking the balance

Check the package, the packaging elements and the delivered components for damages. If any components are damaged, please contact your METTLER TOLEDO service representative.


i Note

Depending on the balance model, the packaging elements and the components may look different.

- 1 Open the box and lift the package out using the lifting strap (1).



- 2 Open the lifting strap (1) and remove the User Manual (2).




- 3 Remove the upper part of the package and carefully unpack the platform (3).


NOTICE: Damage to the instrument

Do not touch the cone (4) protruding from the platform.

- 4 Remove the protective bag.
- 5 Keep the protective covers installed on the platform and on the terminal.



- 6 Open the lower part of the package (5).
- 7 Carefully lift out the draft shield (6), the draft shield doors (7), and all other items.
- 8 Store all parts of the packaging in a safe place for future use.
 - ➔ The balance is ready for assembling.



4.3 Installation

Note

Depending on the balance model, the components may look different.

4.3.1 Balances with draft shield




CAUTION

Injury due to sharp objects or broken glass


Instrument components, e.g., glass, can break and lead to injuries.

- Always proceed with focus and care.


- 1 Assemble the draft shield: Turn the QuickLock (1, right, left) and slide in the top door (2).




- 2 Attach the front panel (3), then turn the QuickLock (1, right, left) to hold the panel in place.




- 3 For balances without backlight: Attach the back panel (4) with the smooth side facing towards the front.




- 4 Slide in the side door (5) until the QuickLock (6) snaps in (right, left).
 - ➔ The draft shield is assembled.




- 5 Place the draft shield (7) on top of the platform (8).
 - i Note**
To protect your balance, keep the protective covers installed on the platform (8) and on the terminal (9).
- 6 Secure the draft shield (7) to the platform (8) by turning the QuickLock (10).



- 7 Fully open the side door (11).
- 8 Insert the drip tray (12).
 - i Note**
If you cannot install the drip tray properly, make sure the QuickLock (10) is correctly locked.




- 9 Only for balances with a readability of 0.01 mg: Place the draft-protection element (13) on top of the drip tray (12).
- 10 Install the weighing pan (14).
 - ➔ The balance is ready for use.



4.3.2 Balances without draft shield


4.3.2.1 Assembling balances 1 mg

- 1 Place the drip tray (1) on top of the platform (2).
Note
To protect your balance, keep the protective covers installed on the platform (2) and on the terminal (3).
- 2 Place the weighing pan (4) on top of the drip tray (1).
➔ The balance is ready for use.




4.3.2.2 Assembling balances 0.01 g

- 1 Place the drip tray (1) on top of the platform (2).
Note
To protect your balance, keep the protective covers installed on the platform (2) and on the terminal (3).
- 2 Attach the support caps (4).
- 3 Place the weighing pan (5) on top of the support caps (4).
➔ The balance is ready for use.



4.3.2.3 Assembling balances 0.1 g

- 1 Place the drip tray (1) on top of the platform (2).
Note
To protect your balance, keep the protective covers installed on the platform (2) and on the terminal (3).
- 2 Attach the support caps (4).
- 3 Place the weighing pan support (5) on top of the support caps (4).
- 4 Place the weighing pan (6) on top of the weighing pan support (5).
➔ The balance is ready for use.




4.3.3 Balances, large

- 1 Remove the transport locks (1).
- 2 Place the weighing pan (2) on top of the support caps (3).

Note

To protect your balance, keep the protective cover installed on the terminal (4).

➔ The balance is ready for use.



4.4 Putting into operation

4.4.1 Connecting the balance




WARNING

Death or serious injury due to electric shock

Contact with parts that carry a live current can lead to death or injury.

- 1 Only use the METTLER TOLEDO power cable and AC/DC adapter designed for your instrument.
- 2 Connect the power cable to a grounded power outlet.
- 3 Keep all electrical cables and connections away from liquids and moisture.
- 4 Check the cables and the power plug for damage and replace them if damaged.

- 1 Install the cables in such a way that they cannot be damaged or interfere with operation.
- 2 Insert the plug of the AC/DC adapter (1) into the power socket of the instrument.
- 3 Secure the plug by firmly tightening the knurled nut.
- 4 Insert the plug of the power cable into a grounded power outlet that is easily accessible.
➔ The balance automatically switches on.



Note

Do not connect the instrument to a power outlet controlled by a switch. After switching on the instrument, it must warm up before giving accurate results.

See also

[General data](#) ▶ Page 125

4.4.2 Switching on the balance

When connected to the power supply, the balance automatically switches on.

EULA (End User License Agreement)

When the balance is switched on the first time, the EULA (End User License Agreement) appears on the screen.

- 1 Read the conditions.
- 2 Tap **I accept the terms in the license agreement** and confirm with **✓ OK**.
 - ➔ The main weighing screen appears.

Acclimatization and warm up

Before the balance gives reliable results, it must:





- acclimatize to the room temperature
- warm up by being connected to the power supply

The acclimatization time and warm-up time for balances are available in "General data".

Note

When the balance is exiting standby, it is ready immediately.


See also

-  General data ▶ Page 125
-  Entering / Exiting standby mode ▶ Page 32
-  Entering / Exiting power-saving mode ▶ Page 32
-  Switching off the balance ▶ Page 32

4.4.3 Leveling the balance


Exact horizontal and stable positioning are essential for repeatable and accurate weighing results.

If the balance is out of level, the level indicator on the main screen turns red.

- 1 On the main weighing screen, tap .
 - ➔ The dialog **Leveling aid** opens.
- 2 Turn the leveling feet (**1**) as instructed on the display until the dot is in the center of the level indicator.

Alternative access to the dialog **Leveling aid**:

 **Navigation:**  **Balance menu** >  **Leveling aid**





4.4.4 Performing an internal adjustment

 **Navigation:**  >  **Applications** >  **Adjustments**



■ **Adjustments** is set to **Internal**.

- 1 Option 1: On the main weighing screen, tap  **Adjust**.
Option 2: Open the applications section, tap  **Adjustments**, select the adjustment, and tap  **Start**.
 - ➔ The adjustment is executed.
 - ➔ The adjustment results appear.
- 2 Tap **✓ Finish**.
 - ➔ The balance is ready.

4.4.5 Entering / Exiting standby mode

- 1 To enter standby mode, short press 
 - ➔ The display is blue. A QR code for more information about the balance is shown.
- 2 To exit standby mode, short press 
 - ➔ The display is turned on.

4.4.6 Entering / Exiting power-saving mode



- 1 To enter power-saving mode, long press  (more than two seconds).
 - ➔ The display is dark. The balance is in power-saving mode.
- 2 To exit power-saving mode, long press 
 - ➔ The balance is switched on.

Note

We recommend configuring power-saving times. When the balance exits power-saving mode automatically at the defined time, the balance is ready for use immediately.

If the power-saving mode is terminated manually, the balance must warm up before it can be used.

See also

-  General data ▶ Page 125
-  Standby, Power-saving mode ▶ Page 38




4.4.7 Switching off the balance

To completely switch off the balance, it must be disconnected from the power supply. By pressing , the balance goes only into standby mode or into power-saving mode.

Note

When the balance has been completely switched off for some time, it must warm up before it can be used.

See also

-  Switching on the balance ▶ Page 31
-  Entering / Exiting standby mode ▶ Page 32
-  Entering / Exiting power-saving mode ▶ Page 32

4.5 Performing a simple weighing

Note


A balance with draft shield is used to explain the procedure. For balances without a draft shield, skip the instructions steps concerning the draft shield.

4.5.1 Opening and closing draft shield doors


- 1 Open the door manually with the door handle (1).
- 2 Alternatively, use the ErgoDoor handle (2) to open the side door.

Note


The following instructions describe one use case for weighing in the sample from the right side.




- 3 Connect the ErgoDoor handle (3) to the door handle on the right side (4).



- 4 Move the ErgoDoor handle on the left side (5) to open and close the door on the right side.



See also

 ErgoDoor handle ▶ Page 18


4.5.2 Zeroing the balance

- 1 Open the draft shield.
- 2 Clear the weighing pan.
- 3 Close the draft shield.
- 4 Press **→0←** to zero the balance.
 - ➔ The balance is zeroed.


4.5.3 Taring the balance

If a sample vessel is used, the balance must be tared.

- The balance is zeroed.
- 1 Place the sample vessel (1) on the weighing pan (2).
 - 2 Press **→T←** to tare the balance.
 - ➔ The balance is tared. The icon ^{Net} appears.



4.5.4 Performing a weighing

- 1 Open the draft shield.
- 2 Place the weighing object into the sample vessel.
- 3 Close the draft shield.
 - ➔ The result is displayed.
- 4 Optional, if a printer is connected: Tap  to print the weighing result.

4.6 Transporting, packing, and storing

4.6.1 Transporting the balance over short distances



NOTICE


Flickering backlight due to disconnection of the draft shield (balances with readability of 0.01 mg)

Hold the balance by the platform. Never lift the balance by the draft shield.





- 1 Disconnect the AC/DC adapter and unplug all interface cables.
- 2 Hold the balance with both hands and carry it in horizontal position to the target location. Consider the requirements of the location.

If you want to put the balance into operation, proceed as follows:

- 1 Connect in reverse order.
- 2 Give the balance sufficient time to warm up.
- 3 Level the balance.
- 4 Perform an internal adjustment.



See also

-  [Selecting the location](#) ▶ Page 25
-  [Switching on the balance](#) ▶ Page 31
-  [Leveling the balance](#) ▶ Page 31
-  [Performing an internal adjustment](#) ▶ Page 31

4.6.2 Transporting the balance over long distances

METTLER TOLEDO recommends using the original packaging for transportation or shipment of the balance or balance components over long distances. The elements of the original packaging are developed specifically for the balance and its components and ensure maximum protection during transportation.

See also

-  [Unpacking the balance](#) ▶ Page 25

4.6.3 Packing and storing

Packing the balance

Store all parts of packaging in a safe place. The elements of the original packaging are developed specifically for the balance and its components, and ensures maximum protection during transportation and storage.

Storing the balance

Only store the balance under the following conditions:

- Indoor and in the original packaging
- According to the environmental conditions, see "Technical Data"

Note

When storing for longer than 6 months, the rechargeable battery may become empty (only date and time get lost).

See also

 Technical Data ▶ Page 125

4.7 Weighing below the balance

Your balance is equipped with a weighing hook for performing weighing operations below the work surface (weighing below the balance).

- A weighing table or workbench is available, through which the weighing hook can be accessed.

- 1 Disconnect the balance from the AC/DC adapter.
- 2 Disconnect all interface cables.
- 3 Carefully tilt the balance to its side.
- 4 Remove the weighing hook cover (1).


Note

The location of the weighing hook depends on the balance model.

➔ The hook (2) is accessible.

- 5 Carefully put the balance back on its feet.
- 6 Reconnect the AC/DC adapter and the interface cables.

➔ The weighing hook is accessible and can be used for below-the-balance weighing.



5 Operation

5.1 Touchscreen


5.1.1 Selecting or activating an item

1. Tap the item or function you want to select or activate.




5.1.2 Scrolling

1. Scroll up or down to view all items.



5.1.3 Opening the fly-in panel

1. Tap the tab, or pull the tab down, to open the fly-in panel.



5.1.4 Entering characters and numbers


Note

Tapping and holding a character gives access to special characters.




| | Name | Description |
|---|--------------------------------|---|
| 1 | Input field | Shows the entered characters and numbers. |
| 2 | Section title | Shows the icon and the title of the current section. |
| 3 | Discard | Closes the keyboard dialog. |
| 4 | Backspace | Tapping the backspace icon deletes the last character of the entry. Tapping and holding the backspace icon deletes the whole entry. |
| 5 | Numbers and special characters | Allows to enter special characters. |
| 6 | Confirm | Confirms the entered data. |
| 7 | Shift | Switches between lowercase and uppercase. |

5.1.5 Changing values



| | Name | Description |
|---|--------------|--------------------------|
| 1 | Plus button | Increases the value. |
| 2 | Value field | Shows the defined value. |
| 3 | Minus button | Decreases the value. |

5.1.6 Sliding



- Move the slider left or right to change the value.

5.2 General balance settings

5.2.1 Date / Time / Language

≡ **Navigation:** ≡ **Balance menu** > ⚙️ **Settings** > ⚖️ **Balance** > 🕒 **Date / Time / Language**

- The setting **Date / Time / Language** is open.
- 1 Optional: Tap the settings **Date format** and **Time format** to define how date and time are displayed.
- 2 Tap the setting **Date** to adjust the date.
- 3 Tap ✓ **OK**.

- 4 Tap the setting **Time** to adjust the time.
- 5 Tap ✓ **OK**.
- 6 Tap the setting **System language** and select the language of your choice.
- 7 Tap ✓ **OK**.
- 8 Tap ✓ **Save**.

Note

The internal clock can show a time deviation. Adjust the time if necessary.

Note

The system language applies to all users if the function **User management** is inactive. If the function **User management** is activated, users can set the system language individually under the menu topic **User management**.

See also

- [Setting the user language ▶ Page 70](#)
- [Settings: Date / Time / Language ▶ Page 80](#)

5.2.2 Screen / StatusLight / Sound

≡ **Navigation:** ≡ **Balance menu** > ⚙ **Settings** > ⚖ **Balance** > 📱 **Screen / StatusLight / Sound**

- The setting **Screen / StatusLight / Sound** is open.
- 1 Tap the setting **Screen brightness** to adjust the brightness of the display.
- 2 Tap ✓ **OK**.
- 3 For balances with a backlit draft shield: Tap the setting **Draft shield backlight brightness** to adjust the brightness.
 - Note**
This function can be deactivated.
- 4 Tap ✓ **OK**.
- 5 Tap the setting **Sound volume** to adjust the volume level.
- 6 Tap ✓ **OK**.
- 7 Tap the setting **StatusLight** to adjust the light strip on the terminal.
 - Note**
This function can be deactivated.
- 8 Tap ✓ **OK**.
- 9 If available: Tap the setting **StatusLight brightness** to adjust the brightness.
- 10 Tap ✓ **OK**.
- 11 Tap ✓ **Save**.

See also

- [Settings: Screen / StatusLight / Sound ▶ Page 81](#)

5.2.3 Standby, Power-saving mode

The function **Standby** helps to save power consumption during the working hours. Outside of the working hours, the function **Power-saving mode** serves to hibernate the balance.

When the balance is switched on from status **Standby**, it is ready for use immediately. When switched on from status **Power-saving mode**, the balance must warm up before use.

≡ **Navigation:** ≡ **Balance menu** > ⚙️ **Settings** > ⚖️ **Balance** > ⚙️ **General**

■ The setting **General** is open.

1 Tap the setting **Standby**.

Note

This function can be activated or deactivated.

2 Enter the time after which the balance goes into standby.

3 Tap ✓ **OK**.

4 Tap the setting **Power-saving mode**.

Note

This function can be activated or deactivated.

5 Define the settings **Start work** and **End work**.

Note

When the balance exits power-saving mode automatically at the defined time, the balance is ready for use immediately.

6 Select the working days.

Note

Between the defined settings **Start work** and **End work**, the balance does not go into power-saving mode.

7 Tap ✓ **OK**.

8 Tap ✓ **Save**.

See also

🔗 [Settings: General](#) ▶ Page 82

5.2.4 Weighing / Quality

5.2.4.1 Warnings and reminders

≡ **Navigation:** ≡ **Balance menu** > ⚙️ **Settings** > ⚖️ **Balance** > ⚖️ **Weighing / Quality**

Leveling warning

Typically, the leveling is done when required. If the option **Forced leveling** is selected, the balance must be leveled before it can be used.

■ The setting **Weighing / Quality** is open.

1 Tap the setting **Leveling warning**.

Note

This function can be activated or deactivated.

2 Select the option of your choice.

3 Tap ✓ **OK**.

Calibration reminder

If this function is activated, the balance reminds you when the balance or the test weights need to be calibrated.

– Tap the setting **Calibration reminder** to activate or deactivate the function.

Calibration expired

If this function is activated, the balance cannot be used when the calibration has expired.

– Tap the setting **Calibration expired** to activate or deactivate the function.

Service reminder

If this function is activated, the balance reminds you when a service is due.

– Tap the setting **Service reminder** to activate or deactivate the function.

See also

🔗 [Settings: Weighing / Quality](#) ▶ Page 77

5.2.4.2 Weighing profiles

≡ **Navigation:** ≡ **Balance menu** > ⚙️ **Settings** > ⚖️ **Balance** > ⚖️ **Weighing / Quality** > ⚖️ **Weighing profiles**

A weighing profile serves to adapt the balance to specific requirements. Up to three weighing profiles can be defined.

See also

🔗 [Settings: Weighing / Quality ▶ Page 77](#)

5.2.4.2.1 Environment

This setting serves to adapt the balance to the ambient conditions of a specific location.

- The setting **Weighing profiles** is open.
 - 1 Enter a name for the profile.
 - 2 Tap the setting **Environment**.
 - 3 Select the option that fits the environmental condition.
 - 4 Tap ✓ **OK**.

See also

🔗 [Settings: Weighing / Quality ▶ Page 77](#)

5.2.4.2.2 Weighing mode

This setting defines how weighing signals are filtered. For standard weighing applications, the option **Universal** is appropriate.

- The setting **Weighing profiles** is open.
 - 1 Tap the setting **Weighing mode**.
 - 2 Select the option that fits your needs.
 - 3 Tap ✓ **OK**.

5.2.4.2.3 Value release

This setting defines how quickly a weighing result is considered stable.

- The setting **Weighing profiles** is open.
 - 1 Tap the setting **Value release**.
 - 2 Select the option that fits your needs.
 - 3 Tap ✓ **OK**.

See also

🔗 [Settings: Weighing / Quality ▶ Page 77](#)

5.2.4.3 Test weights

See also

🔗 [Settings: Weighing / Quality ▶ Page 77](#)

5.2.4.3.1 Defining an individual test weight

The user should enter data related to each test weight based on the corresponding certificate. This enables each test weight to be clearly assigned to a specific certificate. Up to 10 test weights can be configured. These test weights can be used to carry out tests and adjustments.

≡ **Navigation:** ≡ **Balance menu** > ⚙️ **Settings** > ⚖️ **Balance** > ⚖️ **Weighing / Quality** > 🗑️ **Test weights**

Note

For each test weight, the actual weight has to be defined. Ideally, a test weight for an external adjustment corresponds to the balance capacity. Alternatively, use the maximum OIML weight recommended for the balance model.

- The section **Test weights** is open.
- 1 Tap **+**.
- 2 Select the option **Test weight**.
- 3 Tap **→ Next**.
- 4 Enter a name for the test weight.
- 5 Tap **→ Next**.
- 6 Enter the nominal weight of the test weight.
- 7 Tap **→ Next**.
- 8 Enter the actual weight of the test weight.
- 9 Tap **✓ OK**.
- 10 Tap **✓ Save**.
 - ➔ The new test weight is added to the list of available test weights.

5.2.4.3.2 Defining a combined test weight

The user can combine test weights to achieve a test weight capacity that is not available as a single standard weight. For example, a weight of 10 g and a weight of 20 g can be combined and used as a test weight of 30 g. Each combined test weight can include two or three test weights. The class of a specific combined weight can only be as good as the worst class of the individual test weights it contains. As for any other test weight, combined test weight can be used to carry out external tests and adjustments.

≡ **Navigation:** ≡ **Balance menu** > ⚙️ **Settings** > ⚖️ **Balance** > ⚖️ **Weighing / Quality** > 🗑️ **Test weights**

- The section **Test weights** is open.
- 1 Tap **+**.
- 2 Select the option **Combined weight**.
- 3 Tap **→ Next**.
- 4 Enter a name for the combined weight.
- 5 Tap **→ Next**.
- 6 Select an appropriate weight class.
- 7 Tap **→ Next**.
- 8 Select the weights you want to combine.
- 9 Tap **✓ OK**.
 - ➔ The new test weight is added to the list of available test weights.
 - ➔ The nominal weight of the combined weight is calculated automatically.
- 10 Tap **✓ Save**.

5.2.4.3.3 Deleting a test weight

≡ **Navigation:** ≡ **Balance menu** > ⚙️ **Settings** > ⚖️ **Balance** > ⚖️ **Weighing / Quality** > 🗑️ **Test weights**

- The section **Test weights** is open.
- 1 Tap **🗑️**.

- 2 Select the weight you want to delete.
- 3 Tap ✓ **OK**.
 - ➔ The test weight is removed from the list.
- 4 Tap ✓ **Save**.

5.3 Weighing applications










A weighing application serves to carry out specific weighing tasks. The balance offers various weighing applications with default parameters.

5.3.1 Weighing applications overview

The section **Weighing** provides an overview of weighing applications available on the balance. This section serves to select a weighing application for a specific weighing procedure.

☰ **Navigation:** ▼ >  **Applications** > 

The following weighing applications are available:

-  **Weighing**
-  **Counting**
-  **Check weighing**
-  **Dynamic weighing**
-  **Formulation**
-  **Totaling**
-  **Back weighing**
-  **Density**
-  **Differential weighing**








5.3.2 General settings for weighing applications

5.3.2.1 Defining target weight and tolerances

Some weighing applications provide the option to define a target weight. You can also define a tolerance range for the weighing result. Instead of a \pm tolerance range, you can define an upper (+) and/or a lower (-) tolerance limit. If the weighing result is out of range, this is indicated on the main weighing screen.

☰ **Navigation:** ▼ >  >  **Weighing** > 

This example shows how to define a target weight and a tolerance range for the application **Weighing**. The procedure for other weighing applications is similar.

- The section  **Main** is open.
- 1 Tap the setting **Target and tolerances**.
 -  **Note**
This function can be activated or deactivated.
 - ➔ The section  **Target** is open.
- 2 Enter a target weight.
Alternatively, tap  to define the target value with an actual weight.
- 3 Tap  **± Tol..**
 -  **Note**
This function can be activated or deactivated.
- 4 Enter a tolerance range [% or g].
 -  **Note**
Tap on the corresponding icon to switch between % and gram.
- 5 Tap ✓ **OK**.
- 6 Tap ✓ **Save**.

- ➔ The target weight and the tolerance range are shown on the main weighing screen.

See also

[Main](#) ▶ Page 86

5.3.2.2 Defining a sample ID

☰ **Navigation:** ▼ >  >  **Weighing** > 

This example shows how to define a sample ID for the application **Weighing**. The procedure for other weighing applications is similar.

- The settings section of the weighing application is open.

1 Tap  **ID format**.

2 Tap **Sample ID**.

 **Note**

This function can be activated or deactivated.

3 Tap **Default value** and enter a value.

4 Tap  **OK**.

Adding a description

You can add up to three descriptions to a sample.

- The section  **ID format** is open.

1 Tap **Description 1**.

 **Note**

This function can be activated or deactivated.

2 Tap **Type** and select the option **Sample**.

3 Tap **Label** to enter a description.

4 Tap  **OK**.

5 Tap **Default value** to enter a value.

6 Tap  **OK**.

7 Tap **Input prompt**. If this option is activated, you are prompted to enter a value for the sample ID.

8 Tap  **OK**.

9 Tap  **Save**.

See also

[ID format](#) ▶ Page 86

5.3.2.3 Configuring a weighing application

☰ **Navigation:** ▼ >  >  **Weighing** > 

This example shows how to configure the application **Weighing**. The procedure for other weighing applications is similar.

- The settings section of the weighing application is open.

1 Tap  **Weighing**.

2 Tap **Info weight** and select a unit for the secondary weight to be displayed on the main weighing screen.

 **Note**

This function can be activated or deactivated.

3 Tap **Weighing profile** and select the option of your choice.

4 Tap  **OK**.

5 Tap **Weight capture mode** and select the option of your choice.

6 Tap  **Save**.

See also









 Weighing ► Page 89

5.3.2.4 Configuring a weighing series

Some weighing applications provide the option to define a weighing series. If this option is activated, the balance can also provide statistical calculations.

☰ **Navigation:** ▼ >  >  **Weighing** > 

This example shows how to configure a weighing series for the application **Weighing**. The procedure for other weighing applications is similar.

- The settings section of the weighing application is open.
- 1 Optional, if activated: Tap  **ID format**, then tap **Description**.
- 2 Tap **Type** and select the option **Series**.
- 3 Tap  **OK**.
- 4 Tap **Automatic** value if you want to activate an automatic timestamp.
 **Note**
If activated, this setting deactivates the options **Default value** and **Input prompt**.
- 5 Tap  **OK**.
- 6 Tap  **Weighing**.
- 7 Tap **Measurement series** to activate this function.
- 8 Optional: Tap **Statistical calculations** to activate this function.
- 9 Optional: Tap **Acceptance range** and enter a value.
 **Note**
This function can be activated or deactivated.
- 10 Tap  **OK**.
- 11 Tap  **Save**.

See also








 Weighing ► Page 89

5.3.2.5 Using automated functions

Most weighing applications provide the option to automate specific functions. For example, with the option **Automatic tare**, the balance automatically stores the first stable weight as the tare weight.


☰ **Navigation:** ▼ >  >  **Weighing** > 

This example shows how to select automated functions for the application **Weighing**. All functions can be activated or deactivated individually. The procedure for other weighing applications is similar.

- The section  **Automation** is open.
- 1 Tap **Automatic zero** and enter a threshold below which the balance is automatically zeroed.
 **Note**
Select the unit of your choice.
- 2 Tap  **OK**.
- 3 Tap **Automatic tare** to activate or deactivate this function.
- 4 Tap **Preset tare** and enter a fixed tare weight.
Alternatively, tap  to define the value using an actual tare container.
 **Note**
Select the unit of your choice.
- 5 Tap  **OK**.
- 6 Tap **Recall weight** and select the option of your choice.
- 7 Tap  **OK**.

- 8 Tap **✓ Save**.
 - ➔ The main weighing screen opens.
- 9 Tap **→PT←** to set a pretare value.
- 10 Enter a preset tare weight.
- 11 Tap **✓ Apply**.
 - ➔ The preset tare weight is shown on the main weighing screen.

See also

 Automation ▶ Page 88


5.3.2.6 Configuring a report

By default, only the weighing result and the weight unit are published. The report can be configured to show more information. The report defines the content for the following publishing strategies:

- printing data on a printer
- exporting data to a file on a USB storage device
- transferring data to the software **EasyDirect Balance**

☰ **Navigation:** ▼ >  >  **Weighing** > 

This example shows how to configure a report for the application **Weighing**. The procedure for other weighing applications is similar.

- The settings section of the weighing application is open.
 - 1 Tap  **Report**.
 - 2 Tap **Header and Footer**.
 - 3 Tap the items you want to have in the report.
 - 4 Tap **Title** to enter a title name.
 - 5 Tap **✓ OK**.
 - 6 Tap **Empty lines** and enter a number.
 - 7 Tap **✓ OK**.
 - 8 Tap **➤** to get to the next section of the report configuration.
 - 9 Tap the items you want to have in the report.
 - 10 Continue until you have reached the last section of the report configuration.
 - 11 Tap **✓ Save**.

See also

 Report ▶ Page 88

5.3.3 Application "Weighing"



The application **Weighing** offers basic weighing functions. This application is used for simple weighing tasks, or to perform a measurement series.


The settings of the weighing item, such as target weight and tolerances, can be specified.


☰ **Navigation:** ▼ >  >  **Weighing**

Example procedure

This example shows how to weigh a sample. We use a balance with a draft shield.

- 1 Open the applications section.
 - ➔ The section  is selected.
- 2 Tap  **Weighing**.
 - ➔ The corresponding weighing application opens.
- 3 Press **→0←** to zero the balance.

- 4 Open the draft shield door (if applicable).
- 5 Place the weighing object (1) on the weighing pan (2).
- 6 Close the draft shield door (if applicable).
- 7 Wait until the weight stabilizes.
 - ➔ The result is displayed.
- 8 Optional, depending on the settings: Tap  **Publish** to print or export the weighing result.



See also

 Settings: application "Weighing" ▶ Page 85







5.3.4 Application "Counting"

The application **Counting** is used to determine the number of pieces put on the weighing pan. It is advantageous if all pieces are of approximately equal weight, since the unit quantity is determined on the basis of the average weight of a defined number of reference pieces.

 **Navigation:** ▼ >  >  **Counting**

Example procedure

This example shows how to weigh pieces in a sample vessel.

- 1 Open the applications section.
 - ➔ The section  is selected.
- 2 Tap  **Counting**.
 - ➔ The corresponding weighing application opens.
- 3 Tap the title section **Reference**.
Alternatively, tap  to access this setting.
 - ➔ The screen to define the reference weight opens.
- 4 Tap on the left title section to define the number of reference pieces. For example, enter 5.
- 5 Tap  **OK**.
- 6 Place an empty sample vessel on the weighing pan.
- 7 Press  to tare the balance.
- 8 Put the five reference pieces into the sample vessel.
 - ➔ The total weight of the reference pieces is displayed.
- 9 Tap  **OK**.
 - ➔ The number of reference pieces is indicated.
 - ➔ The weight of one reference piece is displayed in the left title section.
- 10 Add pieces to the sample vessel.
 - ➔ The total number of pieces in the sample vessel is displayed.

See also

 Settings: application "Counting" ▶ Page 88

5.3.5 Application "Check weighing"

The application **Check weighing** checks the deviation of a sample weight within a tolerance limit against a reference target weight.

Navigation: Check weighing


Example procedure

This example shows how to check a sample against a target weight. We use a \pm tolerance range.

- 1 Open the applications section.
 - ➔ The section is selected.
- 2 Tap **Check weighing**.
 - ➔ The corresponding weighing application opens.
- 3 Tap .
- ➔ The section **Main configuration** opens.
- 4 Tap **Target weight** and enter a value for the reference sample.
 - Note**
Alternatively, tap to weigh the reference sample.
- 5 Tap and enter a value for the tolerances.
- 6 Tap **OK**.
 - ➔ The section **Main configuration** opens.
- 7 Tap **Check threshold** and enter a value.
- 8 Tap **OK**.
 - ➔ The section **Main configuration** opens.
- 9 Tap **Save**.
 - ➔ The main weighing screen opens.
- 10 Place the sample on the weighing pan.
 - ➔ The result is displayed.

Note

If the result is within the tolerance range, the background appears green.



If the result is outside the tolerance range, the background appears red.



See also

Settings: application "Check weighing" ▶ Page 91










5.3.6 Application "Dynamic weighing"

The application **Dynamic weighing** determines the weight of unstable samples. It also allows weighing under unstable ambient conditions. The calculated weight is an average of several weighings over a defined time.

Navigation: > > **Dynamic weighing**

Example procedure


This example shows how to manually start a dynamic weighing in a sample vessel.

- 1 Open the applications section.
 - ➔ The section  is selected.
- 2 Tap  **Dynamic weighing**.
 - ➔ The corresponding weighing application opens.
- 3 Tap the title section to define the measuring duration in seconds. For example, enter 5.
 - Note**
Alternatively, tap  to access this setting.
- 4 Tap  **OK**.
- 5 Tap **Start mode**.
- 6 Select **Manual**.
- 7 Tap  **OK**.
- 8 Tap  **Save**.
 - ➔ The main weighing screen opens.
- 9 Place an empty sample vessel on the weighing pan.
- 10 Press  to tare the balance.
 - ➔ Net is displayed.
- 11 Place the sample into the sample vessel.
 - ➔ The result is displayed.
- 12 Tap  **Start**.
 - ➔ The balance is capturing the dynamic weight for the defined measuring duration.
 - ➔ The result is displayed on a blue background.
- 13 Tap  **Finish**.

See also

 Settings: application "Dynamic weighing" ▶ Page 93






5.3.7 Application "Formulation"


The application **Formulation** is used to weigh in multiple components one after the other. The balance displays the total weight of the added components. The function  **Fill up** serves to add a component to reach a defined target weight.

Navigation: > > **Formulation**

Example procedure Formulation





This example shows how to add components to a sample vessel.

- 1 Open the applications section.
 - ➔ The section  is selected.
- 2 Tap  **Formulation**.
 - ➔ The corresponding weighing application opens.
- 3 Press  to zero the balance.
- 4 Tap  **Start**.
- 5 Place an empty sample vessel on the weighing pan.
- 6 Press  to tare the balance.
 - ➔ Net is displayed.
- 7 Add the first component to the sample vessel.

- 8 Tap **+ Add**.
- 9 Add the second component to the sample vessel.
- 10 Tap **+ Add**.
- 11 Tap  **Complete**.
 - ➔ The result is displayed.

Example procedure Fill up

This example shows how to add a liquid to samples to reach a defined target weight.

- 1 Open the applications section.
 - ➔ The section  is selected.
- 2 Tap  **Formulation**.
 - ➔ The corresponding weighing application opens.
- 3 Press **→0←** to zero the balance.
- 4 Tap **▶ Start**.
- 5 Place an empty sample vessel on the weighing pan.
- 6 Press **→T←** to tare the balance.
 - ➔ Net is displayed.
- 7 Place the sample into the sample vessel.
- 8 Tap **+ Add**.
- 9 Place another sample into the sample vessel.
- 10 Tap **+ Add**.
 - ➔ The total weight of the samples is displayed in the title bar.
- 11 Repeat the procedure with all samples.
- 12 Tap  **Fill up**.
 - ➔ The total weight of the samples is displayed.
- 13 Add liquid to the sample vessel until the desired target weight is displayed.
 - ➔ The weight of the added liquid is displayed in the title section.
- 14 Tap **✓ OK**.
- 15 Tap  **Complete**.
 - ➔ The number of samples and the total weight of the samples is displayed.

See also

 Settings: application "Formulation" ▶ Page 95



5.3.8 Application "Totaling"

The application **Totaling** is used to separately weigh different samples. The balance automatically calculates the sum of the weighings.

 **Navigation:**      **Totaling**

Example procedure

This example shows how to automatically calculate the total weight of separately weighed samples.

- 1 Open the applications section.
 - ➔ The section  is selected.
- 2 Tap  **Totaling**.
 - ➔ The corresponding weighing application opens.
- 3 Press **→0←** to zero the balance.
- 4 Tap **▶ Start**.
- 5 Place the first sample on the weighing pan.

- 6 Wait until the weight stabilizes.
- 7 Tap **+ Add**.
- 8 Remove the sample from the weighing pan.
- 9 Place another sample on the weighing pan.
- 10 Wait until the weight stabilizes.
- 11 Tap **+ Add**.
 - ➔ The total weight of both samples is displayed in the title section.
- 12 Remove the sample from the weighing pan.
- 13 Repeat the procedure for all samples.
- 14 Tap **☰ Complete**.
 - ➔ The number of samples and their total weight is displayed.
- 15 Tap **✓ Complete**.

See also

[Settings: application "Totaling" ▶ Page 97](#)



5.3.9 Application "Back weighing"

The application **Back weighing** is used to calculate the difference of two weighing values.

☰ **Navigation:** ▼ >  >  **Back weighing**

Example procedure

This example shows how to calculate how much sample remains in a sample vessel after the sample vessel has been emptied.

- 1 Open the applications section.
 - ➔ The section  is selected.
- 2 Tap  **Back weighing**.
 - ➔ The corresponding weighing application opens.
- 3 Tap **▶ Start**.
- 4 Place an empty sample vessel on the weighing pan.
 - ➔ The balance is taring.
- 5 Place the sample into the sample vessel.
 - ➔ **Initial weight:** The weight of the sample is displayed.
- 6 Remove the sample vessel from the weighing pan and remove the sample.
- 7 Place the sample vessel on the weighing pan.
 - ➔ **Final weight:** The weight of the remaining sample is displayed.
 - ➔ δ : The weight of the removed sample is displayed in the title section.
- 8 Tap **☰ Complete**.
 - ➔ The result is displayed.
- 9 Tap **✓ Finish**.

See also

[Settings: application "Back weighing" ▶ Page 99](#)

5.3.10 Application "Density"

The application **Density** is used to determine the density of solids. Density determination is carried out based on **Archimedes' principle**: A body immersed in a fluid undergoes an apparent loss in weight that is equal to the weight of the fluid it displaces.

≡ **Navigation:** ▼ >  >  **Density**

Example procedure

This example describes how to determine the density of a solid with the help of a density kit. Distilled water is used as auxiliary liquid.

■ A density kit is available for the balance.

1 Open the applications section.

➔ The section  is selected.


2 Tap  **Density**.

➔ The corresponding weighing application opens.

➔ The title section shows the auxiliary liquid **Distilled water**.

Note

For custom auxiliary liquids, tap the title section and select the option **Custom**. Then define the density of the custom auxiliary liquid.

Alternatively, tap  to access these settings.

3 Tap ► **Start**.

4 Place the density kit with the auxiliary liquid on the weighing pan.

5 Tap ✓ **OK**.

➔ Net is displayed.

6 Place the sample on the weighing pan.

7 Tap ✓ **OK**.

8 Place the sample into the auxiliary liquid.

9 Tap ✓ **OK**.

➔ The result is displayed.

10 Tap ✓ **Finish**.

See also

 Settings: application "Density" ► Page 101

5.3.11 Application "Differential weighing"

The application **Differential weighing** is used to calculate the difference of two or more weighing values. It is similar to the application **Back weighing**, but it provides more options.

≡ **Navigation:** ▼ >  >  **Differential weighing**

Example procedure

This example uses the weighing sequence **Initial weights first**.

1 Open the applications section.


➔ The section  is selected.

2 Tap  **Differential weighing**

➔ The corresponding weighing application opens.

3 Tap the title section to define the number of samples and the weighing sequence.

Note

Alternatively, tap  to access these settings.

4 Tap ✓ **Save**.


5 Tap ► **Start**.

6 Place an empty sample vessel on the weighing pan.

➔ The balance is taring.

➔ Net is displayed.

7 Place the sample into the sample vessel.

- ➔ **Initial weight:** The weight of the sample is displayed.
- 8 Tap ✓ **OK**.
- 9 Remove the sample vessel from the weighing pan.
- 10 Tap ➔ **Next**.
- 11 Repeat steps 6 - 10 for each sample.
- 12 Place the sample vessel with the treated sample 1 on the weighing pan.
 - ➔ **Final weight:** The weight of the sample is displayed.
- 13 Tap ✓ **OK**.
 - ➔ δ : The weight difference is displayed in the title section.
- 14 Remove the sample vessel from the weighing pan.
- 15 Tap ➔ **Next**.
- 16 Repeat steps 12 - 15 for each sample.
- 17 Tap ✓ **OK**.
- 18 Tap  **Complete**.
 - ➔ The result is published.

See also

 Settings: application "Differential weighing" ▶ Page 103

5.4 Adjustments

This section describes how to set up and perform an internal or an external adjustment.

The internal adjustment uses the built-in weights to adjust the balance. Typically, the balance is set to automatically perform an internal adjustment after a certain event.

The external adjustment requires separate weights to adjust the balance. Typically, an external adjustment is only performed when required by the customer's SOP.




 **Navigation:** ▼ >  **Applications** >  **Adjustments**

5.4.1 Adjustment strategy

This setting defines which adjustment type is performed whenever you start an adjustment.

Example procedure

This example shows how to change the adjustment type **Internal** to the adjustment type **External**.

- 1 Open the applications section.
- 2 Tap  **Adjustments**.
- 3 Tap  **External (OFF)**.
 - ➔ The section **Adjustment strategy** opens.
- 4 Select the option **External adjustment**.
- 5 Tap ✓ **Save**.
- 6 Tap  to return to the main weighing screen.




See also

 Settings: Adjustment strategy ▶ Page 106

5.4.2 Editing an adjustment

This example shows how to edit the adjustment type **Internal**. The procedure to edit the adjustment type **External** is similar.

- 1 Open the applications section.
- 2 Tap  **Adjustments**.

- 3 Tap  **Internal**.
 - ➔ The main weighing screen opens.
- 4 Tap .
 - ➔ The section with the settings opens.
- 5 Change the settings according to your needs.
- 6 Tap  **Save**.
 - ➔ The main weighing screen opens.

See also

 Adjustment settings ▶ Page 106

5.4.3 Performing an internal adjustment

 **Navigation:** ▼ >  **Applications** >  **Adjustments**

■ **Adjustments** is set to **Internal**.

- 1 Option 1: On the main weighing screen, tap  **Adjust**.
 Option 2: Open the applications section, tap  **Adjustments**, select the adjustment, and tap  **Start**.
 - ➔ The adjustment is executed.
 - ➔ The adjustment results appear.
- 2 Tap  **Finish**.
 - ➔ The balance is ready.







5.4.4 Performing an external adjustment

An external test weight for an external adjustment has to weigh at least 10% of the balance capacity. External test weights under 10% of the balance capacity are not displayed on the balance.

 **Navigation:** ▼ >  **Applications** >  **Adjustments**

This example shows how to define a test weight and how to perform an external adjustment.

■ **Adjustments** is set to **External**.

- 1 On the main weighing screen, tap  **Adjust**.
 - ➔ If no suitable test weight has been defined, you are prompted to define a test weight.
- 2 Tap .
 - ➔ The section **Test weights** opens.
- 3 Tap the weight of your choice.
- 4 Tap **Actual weight** and enter a value.
- 5 Tap  **OK** twice.
- 6 Tap  **Save**.
 - ➔ The main weighing screen opens.
- 7 Tap  **Adjust**.
- 8 Place the test weight on the weighing pan.
 - ➔ The adjustment is executed.
- 9 When prompted, remove the weight from the weighing pan.
 - ➔ The adjustment results appear.
- 10 Tap  **Finish**.
 - ➔ The balance is ready.

See also

 Settings: External adjustment ▶ Page 108

5.5 Tests





Routine tests ensure accurate weighing results according to GWP® or other quality-management systems. The tests should be performed in regular intervals, and the results should be documented in a traceable way.

METTLER TOLEDO can help you to define the routine tests to be performed based on your process requirements. Please contact your local METTLER TOLEDO representative for additional information.




≡ **Navigation:** ▼ >  **Applications** >  **Tests**

5.5.1 Editing a test

This example shows how to edit the sensitivity test. The procedure to edit other routine tests is similar.

- 1 Open the applications section.
- 2 Tap  **Tests**.
- 3 Tap  **Sensitivity**.
 - ➔ The main weighing screen opens.
- 4 Tap .
- ➔ The section with the settings opens.
- 5 Change the settings according to your needs.
- 6 If required: Tap  to access the list of available test weights.
- 7 Tap ✓ **OK**.
- 8 Tap ✓ **Save**.
 - ➔ The main weighing screen opens.

See also

-  [Settings: Sensitivity test ▶ Page 108](#)
-  [Settings: Repeatability test ▶ Page 110](#)
-  [Settings: Eccentricity test ▶ Page 111](#)

5.5.2 Performing a test



NOTICE

Incorrect weighing results due to wrong handling of the test weights.

- Only handle test weights with gloves, tweezers, weight forks, or weight handles.

5.5.2.1 Sensitivity test

The sensitivity of the balance defines the deviation between the balance reading and the actual load. The sensitivity test allows you to measure the sensitivity using one or two test points.

≡ **Navigation:** ▼ >  **Applications** >  **Tests** >  **Sensitivity**

This example shows how to perform a sensitivity test with one test point. The procedure with two test points or a tare container is similar, but additional test weights and test containers are required.

- **Test point 1** is defined.
 - A test weight and an appropriate tool to handle the test weight are available.
- 1 Open the application **Sensitivity test**.
 - 2 Tap ► **Start**.
 - ➔ The balance is zeroing.
 - 3 Place the test weight on the weighing pan.
 - 4 When prompted, remove the test weight from the weighing pan.

➔ The result is displayed.

5 Tap ✓ **Finish**.

See also

 Settings: Sensitivity test ▶ Page 108

5.5.2.2 Repeatability test

The repeatability test calculates the standard deviation of a series of measurements with a single test weight in order to determine the repeatability of the balance.

Repeatability is highly affected by the ambient conditions (drafts, temperature fluctuations and vibrations) and also by the skill of the person performing the weighing. Therefore, the series of measurements must be carried out by the same operator, in the same location, under constant ambient conditions, and without interruption.

☰ **Navigation:** ▼ >  **Applications** >  **Tests** >  **Repeatability**

This example shows how to perform a repeatability test.

■ A test weight and an appropriate tool to handle the test weight are available.


1 Open the application **Repeatability test**.

2 If required: Tap the left title section to define the nominal weight of the test weight.

3 Tap ✓ **OK**.

4 If required, tap the right title section to define the number of repetitions.

Note

Alternatively, tap  to access these settings.

5 Tap ✓ **Save**.

➔ The main weighing screen opens.

6 Tap ▶ **Start**.

7 Place the test weight on the weighing pan.

8 When prompted, remove the test weight.

➔ The balance is zeroing.

9 Repeat this procedure as many times as defined.

10 When prompted, remove the test weight.

➔ The result is displayed.

11 Tap ✓ **Finish**.

See also

 Settings: Repeatability test ▶ Page 110

5.5.2.3 Eccentricity test

The eccentricity test checks if every eccentric load deviation (corner load deviation) is within the user SOP tolerances. The corner load is the deviation of the measurement value through off-center (eccentric) loading. The corner load increases with the weight of the load and its distance from the center of the weighing pan (1). If the display remains consistent, even when the same load is placed on different parts of the weighing pan, the balance does not have corner load deviation.

The result corresponds to the highest of the four determined eccentric load deviations (2 to 5).

☰ **Navigation:** ▼ >  **Applications** >  **Tests** >  **Eccentricity**

This example shows how to perform an eccentricity test.

■ A test weight and an appropriate tool to handle the test weight are available.


1 Open the application **Eccentricity test**.

2 If required: Tap the title section to define the nominal weight of the test weight.

Note

Alternatively, tap  to access this setting.

- 3 Tap ► **Start**.
- 4 When prompted, place the test weight on the appropriate positions of the weighing pan.
- 5 When prompted, remove the test weight.
 - ➔ The result is displayed.
- 6 Tap ✓ **Finish**.



See also

🔗 Settings: Eccentricity test ► Page 111

5.6 Interfaces

5.6.1 Ethernet

☰ **Navigation:** ☰ **Balance menu** > ⚙️ **Settings** > 📶 **Interfaces**

This example shows how to configure the balance such that it can communicate with a peripheral device or a service through Ethernet.



NOTICE

Possible electromagnetic interference with other devices

If the Ethernet cable is longer than 30 meters, electromagnetic interference with other devices may be possible.

- Use an Ethernet cable that is shorter than 30 meters.

- The section **Interfaces** is open.

- 1 Tap the function **Ethernet**.

Note

This function can be activated or deactivated.

- 2 Tap the parameter **Host name** to change the name.
- 3 Tap ✓ **OK**.
- 4 Tap the parameter **Network configuration**.
- 5 Select the option of your choice.
- 6 Tap ✓ **OK**.
- 7 If you selected the option **Manual**: Change the other parameters if needed, for example, **IP address**.
- 8 Tap ✓ **OK**.
- 9 Tap ✓ **Save**.

See also

🔗 Settings: Interfaces ► Page 82

5.6.2 Bluetooth

☰ **Navigation:** ☰ **Balance menu** > ⚙️ **Settings** > 📶 **Interfaces**

This example shows how to configure the balance such that it can communicate with a printer through Bluetooth.

Note

This function is only available if a Bluetooth adapter is connected to the balance.

- A Bluetooth adapter is connected to the balance.
- The section **Interfaces** is open.
- 1 Tap the function **Bluetooth**.

Note

This function can be activated or deactivated.

- 2 Tap the parameter **Bluetooth identification** to change the name.
- 3 Tap **✓ OK**.
- 4 Tap **✓ Save**.

See also

[Settings: Interfaces](#) ▶ Page 82

5.7 Devices / Printers

Navigation: **Balance menu** > **Settings** > **Devices / Printers**

See also

[Settings: Devices / Printers](#) ▶ Page 83

[Accessories](#) ▶ Page 143

5.7.1 Printer

Printers serve to document your processes and results. Each weighing application offers the possibility to trigger the printing process manually. The balance can also be configured such that the results are automatically printed.



NOTICE

Damage to the device due to inappropriate use

- Consult the User Manual of the device before using it.

5.7.1.1 Installing a USB printer

Installing and connecting the printer

This example describes how to install a USB printer and connect it to the balance with a USB cable.

Note

A suitable cable from METTLER TOLEDO must be used to ensure proper function.

Navigation: **Balance menu** > **Settings** > **Devices / Printers**

- The USB printer is switched on.
- A suitable cable to connect the printer to the balance is available.
- On the balance, the main weighing screen is open.
- 1 Connect the cable to the USB printer.
- 2 Connect the cable to the USB-A port of the balance.
 - ➔ The printer appears on the list **Devices / Printers**.
 - ➔ The printer is ready for use.

Printing a test page

≡ **Navigation:** ≡ **Balance menu** > ⚙️ **Settings** > 🖨️ **Devices / Printers**

- The printer is connected to the balance.
- The main weighing screen is open.
- 1 Navigate to the section 🖨️ **Devices / Printers**.
- 2 Tap the appropriate printer.
- 3 Tap ✓🖨️.
 - ➔ The printer prints a short text.
- 4 Tap ✓ **OK**.

5.7.1.2 Installing an RS232 printer

Installing and connecting the printer

This example describes how to install an RS232 printer and connect it to the balance with a USB cable.

Note

A suitable cable from METTLER TOLEDO must be used to ensure proper function.

≡ **Navigation:** ≡ **Balance menu** > ⚙️ **Settings** > 🖨️ **Devices / Printers**

- The RS232 printer is switched on.
- A suitable cable to connect the printer to the balance is available.
- On the balance, the main weighing screen is open.
- 1 Connect the cable to the RS232 printer.
- 2 Connect the cable to the USB-A port of the balance.
- 3 Navigate to the section **Devices / Printers**.
- 4 Tap +.
- 5 Select the option **USB-RS232 converter**.
- 6 Tap → **Next**.
- 7 Configure the printer.
- 8 Tap ✓ **Save**.
 - ➔ The printer appears on the list 🖨️ **Devices / Printers**.
 - ➔ The printer is ready for use.

Printing a test page

≡ **Navigation:** ≡ **Balance menu** > ⚙️ **Settings** > 🖨️ **Devices / Printers**

- The printer is connected to the balance.
- The main weighing screen is open.
- 1 Navigate to the section 🖨️ **Devices / Printers**.
- 2 Tap the appropriate printer.
- 3 Tap ✓🖨️.
 - ➔ The printer prints a short text.
- 4 Tap ✓ **OK**.

5.7.1.3 Installing a printer through Bluetooth

This example describes how to install a printer and connect it to the balance through Bluetooth.



For more information about how to install your Bluetooth adapter, consult the Installation Instructions provided with it.


Connecting the printer to the balance

Navigation: ≡ Balance menu > ⚙ Settings > 📶 Interfaces > Bluetooth

Navigation: ≡ Balance menu > ⚙ Settings > 🖨 Devices / Printers

- The printer is switched on.
- A Bluetooth RS adapter (to connect to the printer) and a Bluetooth USB adapter (to connect to the balance) are available.
- The switch on the Bluetooth RS adapter is in the position DCE.
- You have identified the MAC address (unique device address) on the Bluetooth RS adapter.
- The main weighing screen is open.

1 Connect the Bluetooth USB adapter (1) to the USB-A port (2) of the balance.



2 Connect the Bluetooth RS adapter (3) to the printer (4).
➔ The light on the Bluetooth RS adapter starts blinking.

3 Navigate to the section **Bluetooth** and activate the function.

4 Tap ✓ **Save**.

5 Tap 🖨 **Devices / Printers**.

6 Tap +.

7 Select the option **Bluetooth**.

8 Tap → **Next**.

➔ The balance is searching for devices.

9 Select the MAC address of the Bluetooth RS adapter (3).


10 Tap → **Next**.

➔ The balance is pairing the Bluetooth USB adapter (1) with the Bluetooth RS adapter (3) from the printer.

11 Tap → **Next**.

➔ The balance is connecting to the printer.

12 Tap ✓ **Finish**.



Printing a test page

Navigation: ≡ Balance menu > ⚙ Settings > 🖨 Devices / Printers

- The printer is connected to the balance.
 - The main weighing screen is open.
- 1 Navigate to the section 🖨 **Devices / Printers**.
 - 2 Tap the appropriate printer.
 - 3 Tap ✓ 🖨.
➔ The printer prints a short text.
 - 4 Tap ✓ **OK**.

See also

🔗 Settings: Bluetooth ▶ Page 83

5.7.2 Barcode reader

The barcode reader can be used to enter text or numbers in any character-input field on the display. The format of the field must be compatible with the scanned code.



NOTICE


Damage to the device due to inappropriate use

- Consult the User Manual of the device before using it.





5.7.2.1 Scanning a sample ID with a barcode reader

This example shows how to scan a sample ID with a barcode reader in the application **Weighing**.

Installing the barcode reader

- A barcode reader is available.
- The main weighing screen is open.
- Connect the USB cable of the barcode reader to the appropriate USB port of the balance.
 - ➔ The balance automatically detects the barcode reader.
 - ➔ The barcode reader appears on the list  **Devices / Printers**.
 - ➔ The barcode reader is ready for use.

Scanning a sample ID with the barcode reader

- The barcode reader is configured: **End of line character** is set to "Enter".
- The barcode reader is connected to the balance.
- The application **Weighing** is open.
- 1 Tap  to open the settings.
- 2 Tap  **ID format**.
- 3 Tap **Sample ID**.
- 4 Tap **Default value**.
- 5 Scan the code of the sample ID with the barcode reader.
 - ➔ The scanned sample ID appears in the corresponding field.
- 6 Optional: Tap **Default value** again to manually change the scanned sample ID.
- 7 Tap  **OK**.
- 8 Tap  **Save**.

5.7.3 Foot switch

The foot switch can be used to perform certain operations on your balance without having to use the terminal.



NOTICE


Damage to the device due to inappropriate use

- Consult the User Manual of the device before using it.

This example shows how to install and use a foot switch through USB.

Installing and configuring the foot switch

- A foot switch is available.
- The main weighing screen is open.
- 1 Connect the USB cable of the foot switch to the appropriate USB port of the balance.

- ➔ The balance automatically detects the foot switch.
 - ➔ The foot switch appears on the list  **Devices / Printers**.
- 2 Tap on the foot switch.
 - 3 Tap **Function** and configure how the foot switch is to be used.
 - 4 Tap **✓ OK**.
 - 5 Tap **✓ Save**.
 - ➔ The foot switch is ready for use.

5.7.4 Keyboard

A keyboard can be used to perform certain operations on your balance without having to use the terminal.




NOTICE

Damage to the device due to inappropriate use

- Consult the User Manual of the device before using it.

This example shows how to install and use a keyboard through USB.

Installing and configuring the keyboard

- A keyboard with USB cable is available.
 - The main weighing screen is open.
- 1 Connect the USB cable of the keyboard to the appropriate USB port of the balance.
 - ➔ The balance automatically detects the keyboard.
 - ➔ The keyboard appears on the list  **Devices / Printers**.
 - 2 Tap **✓ OK**.
 - ➔ The keyboard is ready for use.

5.7.5 Adding and deleting a device


≡ **Navigation:** ≡ **Balance menu** > ⚙ **Settings** >  **Devices / Printers**

This example shows how to add and delete a printer with a USB interface.

Adding a device

- 1 Open the settings section **Devices / Printers**.
- 2 Tap **+**.
- 3 Select the option **USB**.
- 4 Tap **➔ Next**.
- 5 When prompted, connect the device to the balance.
 - ➔ The device is automatically detected.
- 6 Tap **✓ Save**.
 - ➔ The device appears on the list **Devices / Printers**.

Deleting a device

- 1 Open the settings section **Devices / Printers**.
- 2 Select the device you want to delete.
- 3 Tap .
- 4 Tap **✓ OK**.

5.7.6 Editing the settings of a device

≡ **Navigation:** ≡ **Balance menu** > ⚙ **Settings** > 🖨 **Devices / Printers**

- 1 Open the settings section **Devices / Printers**.
 - ➔ A list of available devices is displayed.
- 2 Adjust the settings if needed.

5.8 Services

The balance offers several ways to control the balance remotely or to manage data.

≡ **Navigation:** ≡ **Balance menu** > ⚙ **Settings** > ⚙ **Services**

See also

🔗 Settings: Services ▶ Page 83

5.8.1 Setting up services

5.8.1.1 MT-SICS service

MT-SICS is a service that allows you to operate the balance by sending commands from a computer. This enables you to integrate your balances into your systems.

The full documentation related to MT-SICS for MX and MR balances is available online.

▶ www.mt.com/labweighing-software-download

This example shows how to establish a connection between your balance and a computer through USB. Other connection options work in a similar way. The computer can then be used to control the balance and receive data using the commands of MT-SICS.

Configuring the balance

≡ **Navigation:** ≡ **Balance menu** > ⚙ **Settings** > ⚙ **Services**

- The section **Services** is open.
- 1 Tap the function **MT-SICS service**.
 - Note**
This function can be activated or deactivated.
 - 2 Tap the setting **Interface**.
 - 3 Select the option **USB**.
 - 4 Tap ✓ **OK**.
 - 5 Tap the setting **Command set**.
 - 6 Select the option **MT-SICS**.
 - 7 Tap ✓ **OK**.
 - 8 Tap ✓ **Save**.

Connecting the balance to the computer

When connecting MT-SICS through USB, a USB driver must be installed on your computer. This creates a COM port for communication with the balance.

The USB driver is available online:

▶ www.mt.com/labweighing-software-download

- The USB driver is installed on the computer.
 - A terminal program is installed and running on the computer.
 - A suitable cable from METTLER TOLEDO is available.
- 1 Provide the necessary connection settings to the terminal program.

- 2 Test the connection by sending a command to the balance, for example, **s** to retrieve the stable weight from the balance.
 - ➔ If a string is received by the terminal program with the weight, date, and time, the connection has been successfully established.
 - ➔ If no response is received by the terminal program, check the connection settings.

See also

- 🔗 Settings: MT-SICS service ▶ Page 84
- 🔗 Transferring data: MT-SICS service ▶ Page 64

5.8.1.2 EasyDirect Balance

This example shows how to establish a connection between your balance and a computer through USB. The computer can then be used to control the balance and receive data using the software **EasyDirect Balance**.

Configuring the balance

☰ **Navigation:** ☰ **Balance menu** > ⚙️ **Settings** > ⚙️ **Services**

- The section **Services** is open.
- 1 Tap the function **EasyDirect Balance**.
 - 📘 **Note**
This function can be activated or deactivated.
- 2 Tap the setting **Interface**.
- 3 Select the option **USB**.
- 4 Tap ✓ **OK**.
- 5 Tap ✓ **Save**.

Connecting the balance to the computer

The software **EasyDirect Balance** must be installed on your computer. A trial version of the software is available online:

▶ www.mt.com/EasyDirectBalance

- A suitable cable from METTLER TOLEDO to connect the balance to the computer is available.
- 1 Install the software **EasyDirect Balance** on your computer.
- 2 Follow the instructions to establish a connection with the balance.

See also

- 🔗 Settings: EasyDirect Balance ▶ Page 84
- 🔗 Transferring data: EasyDirect Balance ▶ Page 65

5.8.1.3 Drop to cursor

This example shows how to configure the balance such that data can be transferred to a computer using the service **Drop to cursor**.

📘 **Note**

The use of special characters is limited when using the service **Drop to cursor**.

Configuring the balance

☰ **Navigation:** ☰ **Balance menu** > ⚙️ **Settings** > ⚙️ **Services**

- The section **Services** is open.
- 1 Tap the service **Drop to cursor**.
 - 📘 **Note**
This function can be activated or deactivated.
- 2 Tap the setting **Interface**.

- 3 Select the option **USB**.
- 4 Tap ✓ **OK**.

Connecting the balance to the computer

- A suitable cable from METTLER TOLEDO is available.
- Connect the balance to the computer.

See also

- 🔗 Settings: Drop to cursor ▶ Page 84
- 🔗 Transferring data: Drop to cursor ▶ Page 66

5.8.1.4 File server

This example shows how to configure the balance such that data can be transferred using the service **File server**.

Configuring the balance

☰ **Navigation:** ☰ **Balance menu** > ⚙️ **Settings** > ⚙️ **Services**

- The interface **Ethernet** is activated.
- The section **Services** is open.
- 1 Tap the service **File server**.
 - Note**
This function can be activated or deactivated.
- 2 Enter the name of the target server.
- 3 Tap ✓ **OK**.
- 4 Enter the name of the shared folder.
- 5 Tap ✓ **OK**.
- 6 Optional: Activate the option **Credentials** to define a user name and a password.
- 7 Tap ✓ **Save**.
 - ➔ The connection to the file server is verified.

See also

- 🔗 Ethernet ▶ Page 56
- 🔗 Settings: File server ▶ Page 85
- 🔗 Transferring data: File server ▶ Page 66

5.8.2 Transferring data to services

This setting serves to define what kind of data is transferred to a target service.

See also

- 🔗 Settings: Publishing ▶ Page 79

5.8.2.1 Transferring data: MT-SICS service

All MX balances can be integrated into a network. The balance can be configured to communicate with a computer. The service MT-SICS (METTLER TOLEDO Standard Interface Command Set) serves to send commands to operate the balance.

For further information, contact your METTLER TOLEDO representative.

The full documentation related to MT-SICS for MX and MR balances is available online.

▶ www.mt.com/labweighing-software-download

Configuring the balance

☰ **Navigation:** ☰ **Balance menu** > ⚙️ **Settings** > ⚖️ **Balance** > 📄 **Publishing**

- A computer is connected to the balance.
- The service **MT-SICS** is activated and configured.
- The section **Publishing** is open.
- 1 Tap the function **Transfer data**.
 - 📘 **Note**
This function can be activated or deactivated.
- 2 Tap the setting **Transfer to**.
- 3 Select the option **MT-SICS service**.
- 4 Tap ✓ **OK**.
- 5 Tap the setting **Output mode** and select the option of your choice.
- 6 Tap ✓ **OK**.
- 7 Tap ✓ **Save**.

Transferring data

In this example, weighing data is transferred to the service **MT-SICS**. The data format is defined in **MT-SICS**.

- A computer with **MT-SICS** is connected to the balance.
- The balance is configured as described above.
- Perform a weighing and tap **Publish**.
 - ➔ The weighing data is sent to the **MT-SICS** client.

See also

🔗 [MT-SICS service](#) ▶ Page 62

5.8.2.2 Transferring data: EasyDirect Balance

EasyDirect Balance is a software to collect, analyze, store, and export measurement results and balance details from up to ten balances.

Configuring the balance

☰ **Navigation:** ☰ **Balance menu** > ⚙️ **Settings** > ⚖️ **Balance** > 📄 **Publishing**

- A computer with the software **EasyDirect Balance** is connected to the balance.
- The service **EasyDirect Balance** is activated and configured.
- The section **Publishing** is open.
- 1 Tap the function **Transfer data**.
 - 📘 **Note**
This function can be activated or deactivated.
- 2 Tap the setting **Transfer to**.
- 3 Select the option **EasyDirect Balance**.
- 4 Tap ✓ **OK**.
- 5 Tap ✓ **Save**.

Transferring data

In this example, weighing data is transferred to the software **EasyDirect Balance**. Which data is transferred is defined in the application-specific section **Report**.

- A computer with the software **EasyDirect Balance** is connected to the balance.
- The balance is configured as described above.
- 1 On the computer, open the software **EasyDirect Balance** and select the balance.

- 2 Perform a weighing and tap **Publish**.
 - ➔ The weighing data is sent to the software **EasyDirect Balance**.

See also

- [EasyDirect Balance ▶ Page 63](#)
- [Configuring a report ▶ Page 45](#)

5.8.2.3 Transferring data: Drop to cursor

The balance offers the option to send weighing results to a computer. This feature can be used, for example, to send results to an Excel sheet or to a text file. With the service **Drop to cursor**, the result is sent to the computer where the cursor is located, as if it were a keyboard input.

Configuring the balance

≡ **Navigation:** ≡ **Balance menu** > ⚙ **Settings** > ⚖ **Balance** > 📄 **Publishing**

- A computer is connected to the balance.
 - The service **Drop to cursor** is activated and configured.
 - The section **Publishing** is open.
- 1 Tap the function **Transfer data**.

Note

This function can be activated or deactivated.

- 2 Tap the setting **Transfer to**.
- 3 Select the option **Drop to cursor**.
- 4 Tap ✓ **OK**.
- 5 In the section **Type of data**, select the type of data you want to transfer.
- 6 In the section **Field configuration**, define the layout of the transferred data.
- 7 Tap ✓ **OK**.
- 8 Tap ✓ **Save**.

Transferring data

In this example, weighing data is transferred to Excel through the service **Drop to cursor**. Which data is transferred is defined here:

≡ **Navigation:** ≡ **Balance menu** > ⚙ **Settings** > ⚙ **Services** > 🖨 **Drop to cursor**

- A computer is connected to the balance.
 - The balance is configured as described above.
- 1 On the computer, open Excel and select a target cell.
 - 2 Perform a weighing and tap **Publish**.
 - ➔ The weighing data is added to the target cell in Excel.
 - 3 The next cell is automatically selected for the next weighing data.

See also

- [Drop to cursor ▶ Page 63](#)

5.8.2.4 Transferring data: File server

This setting serves to define a storage location and a file format for exported data.

☰ **Navigation:** ☰ **Balance menu** > ⚙️ **Settings** > ⚖️ **Balance** > 📄 **Publishing**

- A file server is connected to the balance.
- The section **Report** of the weighing application is configured.
- The section **Publishing** is open.

1 Tap 📄 **Export file**.

Note

This function can be activated or deactivated.

- 2 Tap the option **Export to** and select the option **File server**.
- 3 Tap ✓ **OK**.
- 4 Tap the option **File type** and select a format.
- 5 Tap ✓ **OK**.
- 6 Tap ✓ **Save**.

See also

🔗 File server ▶ Page 64

5.9 Publishing

The balance offers various ways to publish results or transfer data to another device or service. The settings in this section apply to the devices that are defined for the balance. For publishing to services, refer to [Transferring data to services ▶ Page 64].

5.9.1 Printing data

This setting serves to define a target printer and a print format for printed data.

☰ **Navigation:** ☰ **Balance menu** > ⚙️ **Settings** > ⚖️ **Balance** > 📄 **Publishing**

- A printer is connected to the balance.
- The section **Publishing** is open.

1 Tap the function **Printout**.

Note

This function can be activated or deactivated.

- 2 Tap ✓ **OK**.
- 3 Tap ✓ **Save**.

See also

🔗 Settings: Publishing ▶ Page 79

5.9.1.1 Printing results manually through USB



This example shows how to manually print results on a printer that is connected to the balance through USB.

Note

A suitable cable from METTLER TOLEDO must be used to ensure proper function.





- A printer is connected to the balance through USB.
 - A weighing application of your choice is open.
 - The section **Report** of the weighing application is configured.
- 1 Place the sample on the weighing pan.
 - ➔ The result is displayed.
 - 2 Tap 📄 **Publish**.
 - ➔ The result is printed according to the report configuration.

See also

-  Installing an RS232 printer ▶ Page 58
-  Configuring a report ▶ Page 45

5.9.1.2 Printing results automatically through Bluetooth

This example shows how to automatically print results on a printer that is connected to the balance through Bluetooth.

- A printer is connected to the balance through Bluetooth.
 - A weighing application of your choice is open, for example, **Weighing**.
 - The section **Report** of the weighing application is configured.
- 1 Navigate to the settings section of the weighing application, for example, .
 - 2 Tap  **Weighing**.
 - 3 Tap **Weight capture mode**.
 - 4 Select the option **Automatic, stable (zero excluded)** or **Automatic, stable (zero included)**.
 - 5 Tap  **OK**.
 - 6 Tap  **Save**.
 - ➔ The main weighing screen is open.
 - 7 Place a sample on the weighing pan.
 - ➔ The result is printed automatically.

See also

-  Installing a printer through Bluetooth ▶ Page 58
-  Configuring a report ▶ Page 45

5.9.2 Exporting data to a USB storage device

This setting serves to define a storage location and a file format for exported data.

Note

The export takes at least 15 seconds. Do not remove the USB storage device during data export.




 **Navigation:**  **Balance menu** >  **Settings** >  **Balance** >  **Publishing**

- A USB storage device is connected to the balance.
- The section **Report** of the weighing application is configured.
- The section **Publishing** is open.



- 1 Tap  **Export file**.

Note

This function can be activated or deactivated.

- 2 Tap the option **Export to** and select a USB storage device.
- 3 Tap  **OK**.
- 4 Tap the option **File type** and select a format.
- 5 Tap  **OK**.
- 6 Tap  **Save**.

See also

-  Settings: Publishing ▶ Page 79
-  Configuring a report ▶ Page 45

5.9.3 Publishing options

These settings serve to define how a specific type of result is published. A result type can be, for example, test results.

≡ **Navigation:** ≡ **Balance menu** > ⚙ **Settings** > ⚖ **Balance** > 📄 **Publishing**

- The section **Publishing** is open.
- 1 Tap **Single results**.
 - ➔ The information is shown that the behavior is defined in the application-specific setting **Weight capture mode**.
- 2 Tap ✓ **OK**.
- 3 Tap **Workflow results**, **Adjustment results**, and/or **Test results**.
- 4 Select an option.
- 5 Tap ✓ **OK**.
- 6 Tap ✓ **Save**.

See also

🔗 [Settings: Publishing](#) ▶ Page 79

5.9.4 Indicators for weighing results

When published, weighing results might be marked with indicators.

| Indicator | Main weighing screen | Published |
|--------------------|----------------------|-----------|
| Net weight | Net | N |
| Tare weight | — | T |
| Preset tare weight | — | PT |
| Gross weight | — | G |
| Calculated weight | * | * |
| Unstable weight | ○ | D |

5.10 User management

5.10.1 Activating / deactivating the user management

≡ **Navigation:** ≡ **Balance menu** > ⚙ **Settings** > ⚖ **Balance** > ⚙ **General**

- The setting **General** is open.
- 1 Tap the setting **User management**.
- 2 Select the option **Active** or **Inactive**.
- 3 Tap ✓ **OK**.
 - ➔ The current user is logged in as administrator.
 - ➔ When the setting **User management** is deactivated, the current user is logged out automatically.

5.10.2 Managing users and user groups

≡ **Navigation:** ≡ **Balance menu** > 👤 **User management**

5.10.2.1 Automatic logout

≡ **Navigation:** ≡ **Balance menu** > 👤 **User management** > 👤 **User management – General**

- The section **User management – General** is open.
- 1 Tap the setting **Automatic logout**.
 - Note**
This function can be activated or deactivated.
- 2 Define a waiting time before automatic logout.
 - ➔ When the balance is not used, the current user is automatically logged out after the defined waiting time.
- 3 Tap **✓ Save**.

See also

🔗 [User management – General](#) ▶ Page 76

5.10.2.2 Creating a new user

≡ **Navigation:** ≡ **Balance menu** > 👤 **User management** > 👤 **User management – Users**

- The section **User management – Users** is open.
- 1 Tap **+**.
- 2 Enter a user name.
- 3 Tap **➔ Next**.
- 4 Assign a group.
- 5 Tap **➔ Next**.
- 6 Optional: Enter the user's first name and last name.
- 7 Select whether the user is currently active or not active.
- 8 Select a language.
- 9 Optional: Set a password.
- 10 Tap **✓ Save**.
 - ➔ The new user appears on the list of users.

See also

🔗 [User management – Users](#) ▶ Page 76

5.10.2.3 Setting the user language

If the function **User management** is activated, users can set their preferred system language individually.

≡ **Navigation:** ≡ **Balance menu** > 👤 **User management** > 👤 **User management – Users**

- The user is logged in.
- The section **User management – Users** is open.
- 1 Tap the user name.
- 2 Tap **User language**.
- 3 Select a system language.
 - ➔ The system language for this specific user changes to the selected language.
- 4 Tap **✓ Save**.
- 5 Tap **✓ OK**.

See also



[Date / Time / Language ▶ Page 37](#)

[User management – Users ▶ Page 76](#)

5.10.2.4 Deleting a user

☰ **Navigation:** ☰ **Balance menu** > 👤 **User management** > 👤 **User management – Users**

■ The section **User management – Users** is open.



- 1 Tap the user you want to delete.
 - ➔ The user details open.
- 2 Tap .
- 3 Tap  **OK**.
 - ➔ The user is removed from the list of users.

5.10.2.5 Managing groups

This example shows how to manage permissions for a user group. Whether you are allowed to change these settings depends on your permissions.

☰ **Navigation:** ☰ **Balance menu** > 👤 **User management** > 👤 **User management – Groups**

■ The section **User management – Groups** is open.

- 1 Tap a group.
- 2 Tap **Group name** to change the name.
- 3 Tap **Run applications** to select the applications this group is allowed to run.
- 4 Tap  **OK**.
- 5 Tap the other settings to activate or deactivate the corresponding permission.
- 6 Tap  **Save**.

See also

[User management – Groups ▶ Page 76](#)

5.11 Password protection



If the function **User management** is activated, each user has an individual password.

- Users can define and change their own password.
- Users with the permission to configure user management can change the password of any user.
- If users forgot their password, they can request a reset.



5.11.1 Logging in and logging out

If the function **User management** is activated, users need to log in to use the balance.

Logging in

- The login dialog is open.
- 1 Select a user and enter the password.
 - 2 Tap  **OK**.
 - 3 Tap  **Login**.

Logging out

- The user is logged in.
- 1 Tap  **Menu**.
 - 2 Tap  **Logout**.

5.11.2 Changing a password

≡ **Navigation:** ≡ **Balance menu** > 👤 **User management** > 👤 **User management – Users**

- The user is logged in.
- The section **User management – Users** is open.
- 1 Tap the appropriate user.
- 2 Tap 🔑 **Password**.
- 3 Enter the old password.
- 4 Tap ✓ **OK**.
- 5 Enter the new password twice.
- 6 Tap ✓ **OK**.
- 7 Tap ✓ **Save**.

5.11.3 Resetting a password

If users with the permission to configure the function **User management** have forgotten their password, a password reset can be requested.

- The login dialog is open.
- 1 Tap ⋮ **More**.
- 2 Tap ↻ **Request reset password**.
- 3 Enter the user name.
- 4 Tap ✓ **OK**.
- 5 Note the service code and tap ✉ **Service request**.
 - ➔ Information about your METTLER TOLEDO service representative appears.
- 6 Contact your METTLER TOLEDO service representative via phone or email.
 - ➔ You get a temporary password with which you can log in once.
- 7 Log in with your temporary password and select a new password.

5.11.4 Blocking and unblocking the balance

If the function **User management** is activated, the balance can be blocked and unblocked. The balance can only be blocked/unblocked by users with corresponding rights.

Blocking the balance


- A user with access right **Quality management** is present.
- 1 Tap ≡ **Menu**.
- 2 Tap 🔒 **Block**.
- 3 Tap ▶ **Block** to confirm.

Unblocking the balance





- The balance is blocked.
- A user with access right **Quality management** is present.
- 1 Log in to the balance.
 - ➔ The dialog to unblock the balance appears.
- 2 Tap ▶ **Unblock**.
 - ➔ The balance is ready for use.

6 Software Description

6.1 Balance menu settings

The section **Balance menu** contains general settings and information. To open the section **Balance menu**, tap the symbol  on the main screen.

The section **Balance menu** is divided into the following topics:

-  **Leveling aid**
-  **History**
-  **Information**
-  **User management**
-  **Settings**
-  **Maintenance**

6.1.1 Leveling aid

Exact horizontal positioning and stable installation is essential for repeatable and accurate weighing results. The menu topic **Leveling aid** serves to level the balance.

 **Navigation:**  **Balance menu** >  **Leveling aid**

Note

After leveling the balance, an internal adjustment must be performed.

See also







 [Leveling the balance](#) ▶ Page 31

6.1.2 History

The balance records the tests and adjustments that are performed in the menu topic **History**.

 **Navigation:**  **Balance menu** >  **History**



The menu topic **History** is divided into the following sections.

-  **Adjustments history**
-  **Tests history**
-  **Service history**
-  **Activity log**
-  **Software update history**
-  **Error log**

6.1.2.1 Adjustments history

 **Navigation:**  **Balance menu** >  **History** >  **Adjustments history**



A maximum of 500 entries can be stored. If this value is exceeded, the oldest entry is overwritten.

| Button | Name | Description |
|---|----------------|---|
|  | Filter | Tap to filter the list: <ul style="list-style-type: none">• Filter by date• Filter by user |
|  | Publish | Tap to publish or print the displayed entries. |

6.1.2.2 Tests history

≡ Navigation: ≡ Balance menu > 📄 History > 📄 Tests history



A maximum of 500 entries can be stored. If this value is exceeded, the oldest entry is overwritten.

| Button | Name | Description |
|---|---------|---|
|  | Filter | Tap to filter the list: <ul style="list-style-type: none">• Filter by date• Filter by user |
|  | Publish | Tap to publish or print the displayed entries. |

6.1.2.3 Service history

≡ Navigation: ≡ Balance menu > 📄 History > 📄 Service history



A maximum of 100 entries can be stored. If this value is exceeded, the oldest entry is overwritten.

| Button | Name | Description |
|---|---------|---|
|  | Filter | Tap to filter the list: <ul style="list-style-type: none">• Filter by date• Filter by user |
|  | Publish | Tap to publish or print the displayed entries. |

6.1.2.4 Activity log

≡ Navigation: ≡ Balance menu > 📄 History > 📄 Activity log


A maximum of 500 entries can be stored. If this value is exceeded, the oldest entry is overwritten.

| Button | Name | Description |
|---|---------|---|
|  | Filter | Tap to filter the list: <ul style="list-style-type: none">• Filter by date• Filter by user |
|  | Publish | Tap to publish or print the displayed entries. |

6.1.2.5 Software update history

≡ Navigation: ≡ Balance menu > 📄 History > 📄 Software update history



A maximum of 100 entries can be stored. If this value is exceeded, the oldest entry is overwritten.

| Button | Name | Description |
|---|--------|---|
|  | Filter | Tap to filter the list: <ul style="list-style-type: none">• Filter by date• Filter by user |

6.1.2.6 Error log

≡ Navigation: ≡ Balance menu > 📅 History > 📄 Error log

A maximum of 500 entries can be stored. If this value is exceeded, the oldest entry is overwritten.

| Button | Name | Description |
|---|---------|---|
|  | Filter | Tap to filter the list: <ul style="list-style-type: none">• Filter by date• Filter by user |
|  | Publish | Tap to publish or print the displayed entries. |

6.1.3 Information

≡ Navigation: ≡ Balance menu > ⓘ Information

The menu topic **Information** is divided into the following sections:

- 📄 **Balance information**
- 📧 **Service and support information**

6.1.3.1 Balance information

≡ Navigation: ≡ Balance menu > ⓘ Information > 📄 Balance information

The section **Balance information** provides information about the following topics:

- **Balance identification**
- **Logged in user** (if **User management** activated)
- **Software**
- **Hardware**
- **Network**
- **End user license agreement**

6.1.3.2 Service and support information

≡ Navigation: ≡ Balance menu > ⓘ Information > 📧 Service and support information

The section **Service and support information** provides information about the following topics:

- **Service information**
- **Service support contact**

6.1.4 User management

In the menu topic **User management**, rights for users and user groups can be defined. Users can be assigned to user groups.

The menu topic **User management** is only visible if it is activated under the menu topic **Settings**. As a consequence, a login dialog opens at every system start.

A maximum of 20 users can be created. A user is always part of a user group and has the permissions of the corresponding group. Which user has which permissions can be defined or changed by users with the appropriate permission rights.

📘 Note

The settings related to screen brightness and sound can be edited by all users and changes are applied to all users. Any user can set a user-specific language for the balance interface without influencing the settings of other users.

Navigation: Balance menu > User management

The menu topic **User management** is divided into the following sections:

- **User management – General**: settings for all users
- **User management – Users**: settings for individual users
- **User management – Groups**: settings for user groups

See also

[User management](#) ▶ Page 69

6.1.4.1 User management – General

Navigation: Balance menu > User management > User management – General

| Parameter | Description | Values |
|------------------|---|-------------------------------|
| Automatic logout | Defines whether the user is automatically logged out after a predefined waiting time. | Active Inactive* Numeric |

* Factory setting

6.1.4.2 User management – Users

Navigation: Balance menu > User management > User management – Users

| Parameter | Description | Values |
|----------------|--|---------------------|
| User name | Defines a unique identifier for the user. Once the user profile is defined, the value User name cannot be changed anymore. | Text |
| First name | Defines the first name of the user. | Text |
| Last name | Defines the last name of the user. | Text |
| Active | Activates or deactivates the current user. | Active* Inactive |
| Assigned group | Assigns the user to user groups. | Defined groups |
| User language | Defines the language of the user profile. | Available languages |
| Password | Allows the user to set a password. | Active Inactive* |

* Factory setting

6.1.4.3 User management – Groups

Navigation: Balance menu > User management > User management – Groups

Note

This section is only accessible to users with corresponding rights.

| Parameter | Description | Values |
|------------|--------------------------------|--------------------------|
| Group name | Defines the name of the group. | Text (1...22 characters) |

Activity permissions

| Parameter | Description | Values |
|------------------------|---|--|
| Run applications | Defines which applications the group is allowed to run. | Active (All)* Active (number/total number) |
| Run adjustments | Defines whether the group is allowed to perform adjustments. | Active (All)* Inactive |
| Run tests | Defines whether the group is allowed to perform routine tests. | Active (All)* Inactive |
| Configure applications | Defines whether the group is allowed to configure applications. | Active Inactive |
| Cancel results | Defines whether the group is allowed to cancel results. | Active* Inactive |

| | | |
|--------------|--|-------------------|
| Show history | Defines whether the group is allowed to view the menu topic History . | Active Inactive |
|--------------|--|-------------------|

* Factory setting

General configuration permissions

| Parameter | Description | Values |
|--------------------|---|-------------------|
| Quality management | Defines whether the group is allowed to configure the balance settings Weighing / Quality . | Active Inactive |
| User management | Defines whether the group is allowed to configure settings of the menu topic User management . | Active Inactive |
| General | Defines whether the group is allowed to configure the balance settings General . | Active Inactive |

6.1.5 Settings

This section describes the settings of the balance that can be changed to suit specific requirements. The balance settings apply to the entire weighing system and to all users.

≡ **Navigation:** ≡ **Balance menu** > ⚙️ **Settings**

The menu topic **Settings** is divided into the following sections:

- 🏠 **Balance**
- 🖨️ **Interfaces**
- 🖨️ **Devices / Printers**
- ⚙️ **Services**

6.1.5.1 Settings: Balance

≡ **Navigation:** ≡ **Balance menu** > ⚙️ **Settings** > 🏠 **Balance**

The section **Balance** is divided into the following subsections:

- ⚙️ **Weighing / Quality**
- 🖨️ **Publishing**
- 🕒 **Date / Time / Language**
- 🔊 **Screen / StatusLight / Sound**
- ⚙️ **General**

6.1.5.1.1 Settings: Weighing / Quality

≡ **Navigation:** ≡ **Balance menu** > ⚙️ **Settings** > 🏠 **Balance** > ⚙️ **Weighing / Quality**

| Parameter | Description | Values |
|----------------------|--|---|
| Leveling warning | Defines the action when the balance is out of level. For approved balances, the default setting is Forced leveling . | Active* Inactive Optional leveling* Forced leveling |
| Weighing profiles | A weighing profile stores the balance settings needed for a certain weighing application. It is possible to create separate weighing profiles for different weighing applications. Detailed settings are described in the table Weighing profiles below. | Weighing profile 2, Weighing profile 3: Active Inactive |
| Test weights | Allows to define test weights. Detailed settings are described in the table Test weights below. | – |
| Calibration reminder | Defines whether the user is reminded about the upcoming expiry date of the calibration. | Active* Inactive |

| | | |
|---------------------|--|--------------------|
| Calibration expired | Defines whether the balance is blocked if the calibration has expired. | Active Inactive* |
| Service reminder | Defines whether the user is reminded about the upcoming due date of the service. | Active* Inactive |

* Factory setting

Weighing profiles

Settings related to weighing performance and data from balance calibration can be stored in a weighing profile.

| Parameter | Description | Values |
|-------------------------|--|--|
| Profile name | Defines the name of the profile. | Text (1...22 characters) |
| Indicator | Defines the color and the text of the indicator icon. | Active Inactive* Color Text (1...3 characters) |
| Calibration certificate | Defines the ID, the creation date, and the expiry date of the certificate. New certificates can only be created by a service technician based on a performed balance calibration. | Active Inactive* ID (1...32 characters) Date Next date |
| Environment | Defines the environmental conditions of the balance. Stable: For a working environment that is practically free from drafts and vibrations. Standard: For an average working environment subject to moderate variations in the ambient conditions. Unstable: For a working environment subject to variations in the ambient conditions. Very unstable: For a working environment subject to strong variations in the ambient conditions. | Stable Standard* Unstable Very unstable |
| Weighing mode | Defines the filter settings of the balance. Universal: For all standard weighing applications. Sensor mode: Depending on the setting of the ambient conditions, this setting delivers a filtered weighing signal of varying strength. The filter has a linear characteristic in relation to time (not adaptive) and is suitable for continuous measured value processing. | Universal* Sensor mode |
| Value release | Defines the speed at which the balance regards the measured value as stable and available for capture. Very fast: Recommended if you require fast results and repeatability is not very important. Very reliable: Provides very good repeatability of measurement results, but increases stabilization time. | Very fast Fast Fast and reliable* Reliable Very reliable |
| Display readability | Determines the readability d of the balance display. 1d: maximum resolution 2d: 2 times smaller resolution 5d: 5 times smaller resolution 10d: 10 times smaller resolution 100d: 100 times smaller resolution 1000d: 1000 times smaller resolution For approved balances, the values available for this setting depend on the balance model. | 1d* 2d 5d 10d 100d 1000d |
| Zero drift compensation | Performs ongoing corrections of deviations from zero. This setting is not available for approved balances. | Active* Inactive |

| | | |
|---------------------|--------------------------------------|-------------------------------|
| Smallest net weight | Defines the smallest net weight [g]. | Active Inactive* Numeric |
|---------------------|--------------------------------------|-------------------------------|

* Factory setting

Test weights

| Parameter | Description | Values |
|-----------------------|--|--|
| Test weight name | Defines the name of the test weight. | Text (1...22 characters) |
| Test weight ID | Defines the ID of the test weight. | Text (0...22 characters) |
| Nominal weight | Defines the approximate, rounded value of the test weight. | Numeric The available units depend on the balance model. |
| Weight class | Defines the weight class according to OIML or ASTM. Alternatively, create a customized tolerance class with the parameter Own . | E1 E2 F1 F2* M1 M2 M3 ASTM000 ASTM00 ASTM0 ASTM1 ASTM2 ASTM3 ASTM4 ASTM5 ASTM6 ASTM7 Own |
| Actual weight | Defines the actual weight. The actual weight is a specific weight with a specific Conventional Mass Value (CMV) from the weight calibration certificate. | Numeric |
| Next calibration date | Defines the next calibration date. | Active Inactive* Date |
| Weight set ID | Defines the ID of the weight set. | Text (0...22 characters) |

* Factory setting

6.1.5.1.2 Settings: Publishing

Navigation: Balance menu > Settings > Balance > Publishing

| Parameter | Description | Values |
|---------------|---|---|
| Printout | Print to: Defines on which printer the results are printed. Printout type: Defines how the results are printed. | Active* Inactive |
| Export file | Export to: Defines where the results are exported to. File type: Defines the type of the export file. <ul style="list-style-type: none"> USB storage: csv, txt, xml, pdf File server: xml, pdf | Active Inactive* USB storage File server |
| Transfer data | Transfer to: Defines where the results are transferred to when being published. Detailed settings are described in the tables Type of data and Field configuration below. | Active Inactive* Drop to cursor MT-SICS service EasyDirect Balance |

* Factory setting

Publishing options

These settings apply to all available publishing options.

| Parameter | Description | Values |
|------------------|--|----------------------|
| Single results | Single results are published as defined in the setting Weight capture mode . | Application specific |
| Workflow results | Defines whether the workflow results are published immediately after the result is calculated. | Automatic* Manual |

| | | |
|--------------------|--|---|
| Adjustment results | Defines whether the adjustment results are published immediately after the result is calculated. | Automatic Manual* |
| Test results | Defines whether the test results are published immediately after the result is calculated. | Automatic Manual* |
| Approval brackets | Indicates uncertified digits. Relevant for approved balances only. | Active Inactive* <ul style="list-style-type: none"> [] : first decimal place [] : first decimal place for dual range balances |

* Factory setting

Type of data

| Parameter | Description | Values |
|---|--|--------------------|
| Sample ID, Description ID 1, Description ID 2, Description ID 3, Date, Time | Defines whether the corresponding field is included in the output. | Active Inactive* |

* Factory setting

Field configuration

| Parameter | Description | Values |
|-----------------------|--|---|
| Weight status | Defines whether the weight status is included in the output. | Active Inactive* |
| Sign | Defines whether the weighing results are published with a plus sign or a minus sign to indicate positive or negative values. | Always Only negative values* |
| Decimal delimiter | Defines the character used to separate decimal values. | , (comma) . (full stop)* |
| Net indicator | Defines whether the net weights are specially marked in the output. | Active Inactive* |
| Unit | Defines whether the weighing results are published with a unit. | Active* Inactive |
| Field delimiter | Defines the character used to separate data fields. | None TAB* , (comma) ; (semicolon) SPACE |
| End of line character | Defines the character used at the end of a line. | TAB Enter* None |

* Factory setting

See also

[Publishing](#) ▶ Page 67

6.1.5.1.3 Settings: Date / Time / Language

≡ Navigation: ≡ Balance menu > ⚙ Settings > ⚖ Balance > 🕒 Date / Time / Language

| Parameter | Description | Values |
|-----------|---------------------------|--------|
| Date | Defines the current date. | Date |
| Time | Defines the current time. | Time |

| | | |
|----------------------|---|---|
| System language | Defines the language of the interface navigation. This applies to all users if the function User management is inactive. | English* Deutsch Français Español Italiano Polski Český Magyar Nederlands Português Türkçe 中文 日本語 한국어 |
| Show date / time | Shows the current date and time on the screen, in the defined format. | Active* Inactive |
| Time zone | Selects a time zone. When the time zone is set, the balance changes automatically between summer and winter time. | Active Inactive* |
| Time synchronization | Enables synchronization with an NTP server in the network. This setting is only available if the parameter Time zone is activated. | Active Inactive* Text (1...32 characters) |
| Date format | Defines the date format. | DD.MM.YYYY* MM/DD/YYYY YYYY-MM-DD YYYY/MM/DD |
| Time format | Defines the time format. | 24:MM* 12:MM 24.MM 12.MM |

* Factory setting

6.1.5.1.4 Settings: Screen / StatusLight / Sound

≡ Navigation: ≡ Balance menu > ⚙ Settings > ⚖ Balance > *⏪ Screen / StatusLight / Sound

| Parameter | Description | Values |
|-----------------------------------|--|--|
| Screen brightness | Defines the brightness of the display. | 10% ... 100% 60%* |
| Draft shield backlight brightness | Defines the brightness of the draft shield backlight (if applicable). | Active Inactive* 10% ... 100% |
| Sound volume | Defines the volume of the sound. | Inactive Low Medium* High |
| Sound on key press | Defines whether a sound is audible when a key is pressed. | Active* Inactive |
| Sound on feedback | Defines whether a sound is audible when an information appears on the display. | Active* Inactive |
| Sound on stability | Defines whether a sound is audible when the weight value becomes stable. | Active* Inactive |
| StatusLight | <p>Defines whether the light strip on the terminal is used to indicate the status of the balance.</p> <p>Active (without green light): The status of the balance is monitored, but the light strip only glows in red or yellow. The green light is not used.</p> <ul style="list-style-type: none"> Red light: Error. The balance must not be used until the error is corrected. Yellow light: Warning. The balance can still be used. Example: The yellow light glows if you are operating the balance between the date of the calibration reminder and the scheduled date of the next calibration. Green light / no light: No issues detected. The balance is ready for use. | Active* Inactive Active* Active (without green light) |

| | | |
|------------------------|---|--------------|
| StatusLight brightness | Defines the brightness of the light bar on the terminal (StatusLight). This setting is only available if the parameter StatusLight is activated. | 10% ... 100% |
|------------------------|---|--------------|

* Factory setting

6.1.5.1.5 Settings: General

Navigation: Balance menu > Settings > Balance > General

| Parameter | Description | Values |
|-------------------|---|-------------------------------|
| Balance ID | Defines the ID of the balance. This name can be used to communicate with the balance over a network. No space or special characters are allowed. | Text (1...24 characters) |
| Standby | Defines the time before the balance goes into standby mode when it is not in use. | Active* Inactive Numeric |
| Power-saving mode | Defines the working hours and working days. Outside of the defined times, the balance goes into power-saving mode. The setting Start work defines when the balance is ready for use. | Active Inactive* |
| Communication | Defines whether the balance interfaces are open or blocked for communication with connected devices. | Active* Blocked |
| User management | Activates or deactivates the menu topic User management . | Active* Inactive |

* Factory setting

6.1.5.2 Settings: Interfaces

Navigation: Balance menu > Settings > Interfaces

The section **Interfaces** is divided into the following subsections:

- Ethernet
- Bluetooth

| Parameter | Description | Values |
|-----------|--|--------------------|
| Ethernet | With the option Ethernet , the balance can communicate with peripheral devices, such as a printer. | Active Inactive* |
| Bluetooth | With the option Bluetooth , the balance can communicate with peripheral devices, such as a printer. | Active Inactive* |

* Factory setting

See also

Interfaces ▶ Page 56

6.1.5.2.1 Settings: Ethernet

The interface **Ethernet** allows to connect the balance to a network and perform the following actions:

- store weighing results as XML or PDF files
- communicate remotely with the balance using the MT-SICS communication protocol

Navigation: Balance menu > Settings > Interfaces > Ethernet

| Parameter | Description | Values |
|-------------|--|--------------------------|
| Host name | Defines the host name of the balance. | Text (1...24 characters) |
| MAC address | Information on the MAC address that is used to uniquely identify the balance in the network. | not editable |

| | | |
|-----------------------|--|---------------------------------------|
| Network configuration | DHCP: The settings of the Ethernet connection will be automatically set. Manual: The settings of the Ethernet connection must be set manually by the user. If this option is selected, the following parameters are editable. | DHCP* Manual |
| IP address | Defines the IP address of the balance. | 000.000.000.000... 255.255.255.255 |
| Subnet mask | Defines the subnet mask that is used by the TCP/IP protocol to determine whether a host is on the local subnet or on a remote network. | 000.000.000.000... 255.255.255.255 |
| DNS server | Defines the address of the DNS (domain name system) server. | 000.000.000.000... 255.255.255.255 |
| Standard gateway | Defines the address of the standard gateway that links the subnet of the host to other networks. | 000.000.000.000... 255.255.255.255 |

* Factory setting

6.1.5.2.2 Settings: Bluetooth

Navigation: Balance menu > Settings > Interfaces > Bluetooth

| Parameter | Description | Values |
|--------------------------|--|--------------------------|
| Bluetooth identification | Serves to identify the balance when the option Bluetooth is used. | Text (1...24 characters) |

6.1.5.3 Settings: Devices / Printers

Navigation: Balance menu > Settings > Devices / Printers

| Parameter | Description | Values |
|---------------------|--|--|
| Physical connection | Defines the type of physical connection between the balance and a peripheral device. | USB* USB-RS232 converter Network Bluetooth |





* Factory setting

6.1.5.4 Settings: Services

Several services are available to communicate with the balance. Note that only one service can be enabled at any given time.

Navigation: Balance menu > Settings > Services

The section **Services** is divided into the following subsections:

-  **MT-SICS service**
-  **EasyDirect Balance**
-  **Drop to cursor**
-  **File server**

See also

[Services](#) ▶ Page 62

6.1.5.4.1 Settings: MT-SICS service

Navigation: Balance menu > Settings > Services >  MT-SICS service

| Parameter | Description | Values |
|-------------|--|--|
| Interface | If the option MT-SICS service is activated, a corresponding port is opened. | USB-C USB-RS232 converter* Network Bluetooth |
| Command set | Available set of commands to communicate with the balance. | MT-SICS* Sartorius commands 22 Sartorius commands 16 |
| Baudrate | Defines the speed of data transmission. | 600 bps 1200 bps 2400 bps 4800 bps 9600 bps * 19200 bps 38400 bps 57600 bps 115200 bps |
| Bits/Parity | Number of data bits / Checksum for error detection during data transmission | 8/No* 7/No 7/Even 7/Odd |
| Data flow | Also known as "handshake". Defines the synchronization for data transmission. | Xon/Xoff* RTS/CTS None |
| Stop bit | Marks the end of data transmission. | 1 bit* 2 bit |
| Line end | Defines the character at the end of a line. | <CR><LF>* <CR> <LF> <TAB> |

* Factory setting

See also

 MT-SICS service ▶ Page 62

6.1.5.4.2 Settings: EasyDirect Balance

Navigation: Balance menu > Settings > Services >  EasyDirect Balance

| Parameter | Description | Values |
|-----------|---|--|
| Interface | Defines how the service EasyDirect Balance communicates with the balance. If the option Network is selected, the parameter Port can be defined. | USB-C* Network Port: 1024...65535 |

* Factory setting

See also

 EasyDirect Balance ▶ Page 63

6.1.5.4.3 Settings: Drop to cursor

Navigation: Balance menu > Settings > Services >  Drop to cursor

| Parameter | Description | Values |
|-----------|--|--------|
| Interface | The option Drop to cursor can only be used through USB. | USB-C* |

* Factory setting

See also

 Drop to cursor ▶ Page 63

6.1.5.4.4 Settings: File server

Navigation: Balance menu > Settings > Services > File server

| Parameter | Description | Values |
|-------------|--|---------------------------|
| Server name | Defines the name of the target server. | Text (1...63 characters) |
| Share name | Defines the name of the shared folder. | Text (1...140 characters) |

Credentials

| Parameter | Description | Values |
|-------------|--|--------------------------|
| Domain name | Defines the name of the user domain. | Text (1...15 characters) |
| User name | Defines a user name to access the file server. | Text (1...22 characters) |
| Password | Defines a password to access the file server. | Text (1...22 characters) |

See also

[File server](#) ▶ Page 64

6.1.6 Maintenance

Navigation: Balance menu > Maintenance

Note

This section is only accessible to users with corresponding rights.
The menu topic **Maintenance** is divided into the following sections:

- Import / Export
- Software update
- Reset
- Save support file
- Level center adjustment
- Service tool connection

See also

[Exporting data to a USB storage device](#) ▶ Page 68

[Updating the software](#) ▶ Page 119


[Resetting the balance](#) ▶ Page 119

[Saving a support file](#) ▶ Page 124

6.2 Weighing applications settings

6.2.1 Settings: application "Weighing"

This section describes the settings of the application **Weighing**.



The settings for this weighing application are grouped as follows:

- **Main**
- **ID format**
- **Weighing**
- **Automation**
- **Report**

See also

Application "Weighing" ▶ Page 45

6.2.1.1 Main

| Parameter | Description | Values |
|-----------------------|--|--|
| Unit | Defines the unit of the weighing result. | The available units depend on the balance model. |
| Target and tolerances | The target weight can be added manually or by weighing. The definition of tolerances is optional. Depending on the settings, the target weight and the tolerance limits appear on the main weighing screen. The section SmartTrac indicates whether the current weighing result is within the tolerance limits. | Active Inactive* Numeric Tolerances: % g |

* Factory setting

6.2.1.2 ID format

Sample ID

| Parameter | Description | Values |
|-----------------|--|---------------------------|
| Sample ID | Defines a sample identification. | Active Inactive* |
| Default value | Defines a default value for the sample description. | Text (1...200 characters) |
| Automatic value | Defines whether an automatic value is generated for the sample description. | not editable |
| Input prompt | Defines whether you are prompted to enter a value. If a default value is defined, this parameter is not editable. | Active Inactive* |

* Factory setting

Description

| Parameter | Description | Values |
|-----------------|---|---------------------------|
| Description | Allows to define a sample description. | Active Inactive* |
| Type | Defines the sample type. | Sample* Series |
| Label | Describes the sample. | Text (1...25 characters) |
| Default value | Defines a default value for the sample description. | Text (1...200 characters) |
| Automatic value | Defines whether an automatic value is generated for the sample description. | Active Inactive* |
| Input prompt | Defines whether you are prompted to enter a value. | Active Inactive* |

* Factory setting

6.2.1.3 Weighing

| Parameter | Description | Values |
|---------------------|---|--|
| Info weight | A secondary weight is displayed on the main weighing screen. | Active Inactive* The available units depend on the balance model. |
| Weighing profile | Defines the weighing profile. | General* 10d |
| Weight capture mode | Defines the behavior when the button to add the result was tapped, or when the add result was triggered by the automatic weighing result creation. Stable: The balance waits for a stable weight. Immediate: The balance does not wait for a stable weight. Automatic, stable (zero excluded): The results are published as soon as the weight is stable. Values of 0 g are not published. Automatic, stable (zero included): The results are published as soon as the weight is stable. Values of 0 g are also published. Continuous: The results are published at the defined interval. | Stable* Immediate Automatic, stable (zero excluded) Automatic, stable (zero included) Continuous |

* Factory setting

Series / Statistics

| Parameter | Description | Values |
|--------------------------|--|-----------------------------------|
| Measurement series | A measurement series can be performed. | Active Inactive* |
| Statistical calculations | Statistical information is provided. This setting is only available if the parameter Measurement series is activated. | Active Inactive* |
| Acceptance range | Defines the acceptance range for the statistical calculations. This setting is only available if the parameter Statistical calculations is activated. | Active Inactive* Numeric (%) |

* Factory setting

See also

 Weighing profiles ▶ Page 40

6.2.1.4 Automation

| Parameter | Description | Values |
|----------------|---|---|
| Automatic zero | The balance is automatically zeroed when the weight falls below a predefined threshold. This setting is not available for approved balances. | Active Inactive* Numeric The available units depend on the balance model. |
| Automatic tare | The balance automatically stores the first stable weight as the tare weight. | Active Inactive* |
| Preset tare | A fixed tare weight can be defined manually or by weighing. | Active Inactive* Numeric The available units depend on the balance model. |
| Recall weight | Displays the last weighing result. | Active Inactive* Automatic Manual* |

* Factory setting

6.2.1.5 Report






| Parameter | Description | Values |
|---------------------------|---|--|
| Header and Footer | Defines the header and/or the footer to be published. | Title Date/time User Signature Empty lines |
| Balance information | Defines which information about the balance is published. | Balance type Balance ID Balance serial number Software version |
| Quality information | Defines which quality information is published. | Weighing profile Adjustment date/time Routine test information Result state Level state MinWeigh state Tolerance state |
| Task information | Defines which information about the task is published. | Application settings |
| Result detail information | Defines which information related to the measurement result is published. | Tare / Gross weight Info weight Date/time |

6.2.2 Settings: application "Counting"

This section describes the settings of the application **Counting**.

Navigation:  **Counting** > 

The settings for this weighing application are grouped as follows:

-  **Main**
-  **ID format**
-  **Weighing**
-  **Automation**
-  **Report**

See also

 Application "Counting" ▶ Page 46

6.2.2.1 Main

| Parameter | Description | Values |
|--------------------------|--|--|
| Reference PCS | Defines the number of items used to determine the average weight per item. | Numeric |
| Reference average weight | Defines the average weight for one piece. The average weight of one piece serves as the basis for piece counting. During task execution, the balance calculates the actual number of pieces on the weighing pan based on the measured weight and the average weight of one piece. | Numeric |
| Target and tolerances | The target weight can be added manually or by weighing. The definition of tolerances is optional. Depending on the settings, the target weight and the tolerance limits appear on the main weighing screen. The section SmartTrac indicates whether the current weighing result is within the tolerance limits. | Active Inactive* Numeric Tolerances: PCS % |

* Factory setting

6.2.2.2 ID format

Sample ID

| Parameter | Description | Values |
|-----------------|--|---------------------------|
| Sample ID | Defines a sample identification. | Active Inactive* |
| Default value | Defines a default value for the sample description. | Text (1...200 characters) |
| Automatic value | Defines whether an automatic value is generated for the sample description. | not editable |
| Input prompt | Defines whether you are prompted to enter a value. If a default value is defined, this parameter is not editable. | Active Inactive* |

* Factory setting

Description

| Parameter | Description | Values |
|-----------------|---|---------------------------|
| Description | Allows to define a sample description. | Active Inactive* |
| Type | Defines the sample type. | Sample* Series |
| Label | Describes the sample. | Text (1...25 characters) |
| Default value | Defines a default value for the sample description. | Text (1...200 characters) |
| Automatic value | Defines whether an automatic value is generated for the sample description. | Active Inactive* |
| Input prompt | Defines whether you are prompted to enter a value. | Active Inactive* |

* Factory setting

6.2.2.3 Weighing

| Parameter | Description | Values |
|------------------|--|--|
| Info weight | A secondary weight is displayed on the main weighing screen. | Active Inactive* The available units depend on the balance model. |
| Weighing profile | Defines the weighing profile. | General* 10d |

| | | |
|---------------------|---|---|
| Weight capture mode | <p>Defines the behavior when the button to add the result was tapped, or when the add result was triggered by the automatic weighing result creation.</p> <p>Stable: The balance waits for a stable weight.</p> <p>Immediate: The balance does not wait for a stable weight.</p> <p>Automatic, stable (zero excluded): The results are published as soon as the weight is stable. Values of 0 g are not published.</p> <p>Automatic, stable (zero included): The results are published as soon as the weight is stable. Values of 0 g are also published.</p> | Stable* Immediate Automatic, stable (zero excluded) Automatic, stable (zero included) |
|---------------------|---|---|

* Factory setting

Series / Statistics

| Parameter | Description | Values |
|--------------------------|---|-----------------------------------|
| Measurement series | A measurement series can be performed. | Active Inactive* |
| Statistical calculations | <p>Statistical information is provided.</p> <p>This setting is only available if the parameter Measurement series is activated.</p> | Active Inactive* |
| Acceptance range | <p>Defines the acceptance range for the statistical calculations.</p> <p>This setting is only available if the parameter Statistical calculations is activated.</p> | Active Inactive* Numeric (%) |

* Factory setting

See also

[🔗 Weighing profiles](#) ▶ Page 40

6.2.2.4 Automation

| Parameter | Description | Values |
|----------------|--|---|
| Automatic zero | <p>The balance is automatically zeroed when the weight falls below a predefined threshold.</p> <p>This setting is not available for approved balances.</p> | Active Inactive* Numeric The available units depend on the balance model. |
| Automatic tare | The balance automatically stores the first stable weight as the tare weight. | Active Inactive* |
| Preset tare | A fixed tare weight can be defined manually or by weighing. | Active Inactive* Numeric The available units depend on the balance model. |
| Recall weight | Displays the last weighing result. | Active Inactive* Automatic Manual* |

* Factory setting

6.2.2.5 Report

| Parameter | Description | Values |
|-------------------|---|--|
| Header and Footer | Defines the header and/or the footer to be published. | Title Date/time User Signature Empty lines |






| | | |
|---------------------------|---|--|
| Balance information | Defines which information about the balance is published. | Balance type Balance ID Balance serial number Software version |
| Quality information | Defines which quality information is published. | Weighing profile Adjustment date/time Routine test information Result state Level state MinWeigh state Tolerance state |
| Task information | Defines which information about the task is published. | Application settings |
| Result detail information | Defines which information related to the measurement result is published. | Tare / Gross weight Info weight Date/time |

6.2.3 Settings: application "Check weighing"

This section describes the settings of the application **Check weighing**.

☰ **Navigation:** ▼ >  >  **Check weighing** > 

The settings for this weighing application are grouped as follows:

-  **Main**
-  **ID format**
-  **Weighing**
-  **Automation**
-  **Report**

See also

 Application "Check weighing" ▶ Page 46

6.2.3.1 Main

| Parameter | Description | Values |
|-------------------------------|--|--|
| Unit | Defines the unit of the weighing result. | The available units depend on the balance model. |
| Target weight ± Tolerances | The target weight can be added manually or by weighing. The definition of tolerances is optional. Depending on the settings, the target weight and the tolerance limits appear on the main weighing screen. The section SmartTrac indicates whether the current weighing result is within the tolerance limits. | Numeric Tolerances: Active* Inactive % g |
| Check threshold | Defines the target threshold. Values below the defined threshold are not checked. | Active* Inactive Numeric (%) |

* Factory setting

6.2.3.2 ID format

Sample ID

| Parameter | Description | Values |
|---------------|---|---------------------------|
| Sample ID | Defines a sample identification. | Active Inactive* |
| Default value | Defines a default value for the sample description. | Text (1...200 characters) |

| | | |
|-----------------|--|--------------------|
| Automatic value | Defines whether an automatic value is generated for the sample description. | not editable |
| Input prompt | Defines whether you are prompted to enter a value. If a default value is defined, this parameter is not editable. | Active Inactive* |

* Factory setting

Description

| Parameter | Description | Values |
|-----------------|---|---------------------------|
| Description | Allows to define a sample description. | Active Inactive* |
| Type | Defines the sample type. | Sample* Series |
| Label | Describes the sample. | Text (1...25 characters) |
| Default value | Defines a default value for the sample description. | Text (1...200 characters) |
| Automatic value | Defines whether an automatic value is generated for the sample description. | Active Inactive* |
| Input prompt | Defines whether you are prompted to enter a value. | Active Inactive* |

* Factory setting

6.2.3.3 Weighing

| Parameter | Description | Values |
|---------------------|--|---|
| Info weight | A secondary weight is displayed on the main weighing screen. | Active Inactive* The available units depend on the balance model. |
| Weighing profile | Defines the weighing profile. | General* 10d |
| Weight capture mode | Defines the behavior when the button to add the result was tapped, or when the add result was triggered by the automatic weighing result creation. Stable: The balance waits for a stable weight. Immediate: The balance does not wait for a stable weight. Automatic, stable (zero excluded): The results are published as soon as the weight is stable. Values of 0 g are not published. Automatic, stable (zero included): The results are published as soon as the weight is stable. Values of 0 g are also published. | Stable* Immediate Automatic, stable (zero excluded) Automatic, stable (zero included) |

* Factory setting

Series / Statistics

| Parameter | Description | Values |
|--------------------------|--|-----------------------------------|
| Measurement series | A measurement series can be performed. | Active Inactive* |
| Statistical calculations | Statistical information is provided. This setting is only available if the parameter Measurement series is activated. | Active Inactive* |
| Acceptance range | Defines the acceptance range for the statistical calculations. This setting is only available if the parameter Statistical calculations is activated. | Active Inactive* Numeric (%) |

* Factory setting

See also

[🔗 Weighing profiles](#) ▶ Page 40

6.2.3.4 Automation

| Parameter | Description | Values |
|----------------|---|---|
| Automatic zero | The balance is automatically zeroed when the weight falls below a predefined threshold. This setting is not available for approved balances. | Active Inactive* Numeric The available units depend on the balance model. |
| Automatic tare | The balance automatically stores the first stable weight as the tare weight. | Active Inactive* |
| Preset tare | A fixed tare weight can be defined manually or by weighing. | Active Inactive* Numeric The available units depend on the balance model. |
| Recall weight | Displays the last weighing result. | Active Inactive* Automatic Manual* |



* Factory setting

6.2.3.5 Report






| Parameter | Description | Values |
|---------------------------|---|--|
| Header and Footer | Defines the header and/or the footer to be published. | Title Date/time User Signature Empty lines |
| Balance information | Defines which information about the balance is published. | Balance type Balance ID Balance serial number Software version |
| Quality information | Defines which quality information is published. | Weighing profile Adjustment date/time Routine test information Result state Level state MinWeigh state Tolerance state |
| Task information | Defines which information about the task is published. | Application settings |
| Result detail information | Defines which information related to the measurement result is published. | Tare / Gross weight Info weight Date/time |

6.2.4 Settings: application "Dynamic weighing"

This section describes the settings of the application **Dynamic weighing**.

☰ **Navigation:** ▼ >  >  **Dynamic weighing** > 

The settings for this weighing application are grouped as follows:

-  **Main**
-  **ID format**
-  **Weighing**
-  **Automation**
-  **Report**

See also

 Application "Dynamic weighing" ▶ Page 47

6.2.4.1 Main

| Parameter | Description | Values |
|----------------------|--|--|
| Measurement duration | Defines the measuring duration in seconds. | Numeric |
| Start mode | Defines how the measurement is started. | Manual Automatic - After 3 seconds* |
| Unit | Defines the unit of the weighing result. | The available units depend on the balance model. |

* Factory setting

6.2.4.2 ID format

Sample ID

| Parameter | Description | Values |
|-----------------|--|---------------------------|
| Sample ID | Defines a sample identification. | Active Inactive* |
| Default value | Defines a default value for the sample description. | Text (1...200 characters) |
| Automatic value | Defines whether an automatic value is generated for the sample description. | not editable |
| Input prompt | Defines whether you are prompted to enter a value. If a default value is defined, this parameter is not editable. | Active Inactive* |

* Factory setting

Description

| Parameter | Description | Values |
|-----------------|---|---------------------------|
| Description | Allows to define a sample description. | Active Inactive* |
| Type | Defines the sample type. | Sample* Series |
| Label | Describes the sample. | Text (1...25 characters) |
| Default value | Defines a default value for the sample description. | Text (1...200 characters) |
| Automatic value | Defines whether an automatic value is generated for the sample description. | Active Inactive* |
| Input prompt | Defines whether you are prompted to enter a value. | Active Inactive* |

* Factory setting

6.2.4.3 Weighing

| Parameter | Description | Values |
|------------------|--|--|
| Info weight | A secondary weight is displayed on the main weighing screen. | Active Inactive* The available units depend on the balance model. |
| Weighing profile | Defines the weighing profile. | General* 10d |

* Factory setting

Series / Statistics

| Parameter | Description | Values |
|--------------------|--|--------------------|
| Measurement series | A measurement series can be performed. | Active Inactive* |

| | | |
|--------------------------|--|-----------------------------------|
| Statistical calculations | Statistical information is provided. This setting is only available if the parameter Measurement series is activated. | Active Inactive* |
| Acceptance range | Defines the acceptance range for the statistical calculations. This setting is only available if the parameter Statistical calculations is activated. | Active Inactive* Numeric (%) |

* Factory setting

6.2.4.4 Automation

| Parameter | Description | Values |
|----------------|---|---|
| Sample tare | After the result has been calculated, the balance is automatically tared when the sample is removed from the weighing pan. | Active Inactive |
| Automatic zero | The balance is automatically zeroed when the weight falls below a predefined threshold. This setting is not available for approved balances. | Active Inactive* Numeric The available units depend on the balance model. |
| Automatic tare | The balance automatically stores the first stable weight as the tare weight. | Active Inactive* |
| Preset tare | A fixed tare weight can be defined manually or by weighing. | Active Inactive* Numeric The available units depend on the balance model. |

* Factory setting

6.2.4.5 Report

| Parameter | Description | Values |
|---------------------------|---|--|
| Header and Footer | Defines the header and/or the footer to be published. | Title Date/time User Signature Empty lines |
| Balance information | Defines which information about the balance is published. | Balance type Balance ID Balance serial number Software version |
| Quality information | Defines which quality information is published. | Weighing profile Adjustment date/time Routine test information Result state Level state MinWeigh state |
| Task information | Defines which information about the task is published. | Application settings |
| Result detail information | Defines which information related to the measurement result is published. | Tare / Gross weight Info weight Date/time |

6.2.5 Settings: application "Formulation"

This section describes the settings of the application **Formulation**.

☰ **Navigation:** ▼ >  >  **Formulation** > 

The settings for this weighing application are grouped as follows:

-  **Main**

-  ID format
-  Weighing
-  Automation
-  Report

See also

 Application "Formulation" ▶ Page 48

6.2.5.1 Main

| Parameter | Description | Values |
|-----------|--|--|
| Unit | Defines the unit of the weighing result. | The available units depend on the balance model. |

6.2.5.2 ID format

Sample ID

| Parameter | Description | Values |
|-----------------|--|---------------------------|
| Sample ID | Defines a sample identification. | Active Inactive* |
| Default value | Defines a default value for the sample description. | Text (1...200 characters) |
| Automatic value | Defines whether an automatic value is generated for the sample description. | not editable |
| Input prompt | Defines whether you are prompted to enter a value. If a default value is defined, this parameter is not editable. | Active Inactive* |

* Factory setting

Description

| Parameter | Description | Values |
|-----------------|---|---------------------------|
| Description | Allows to define a sample description. | Active Inactive* |
| Type | Defines the sample type. | Sample* Series |
| Label | Describes the sample. | Text (1...25 characters) |
| Default value | Defines a default value for the sample description. | Text (1...200 characters) |
| Automatic value | Defines whether an automatic value is generated for the sample description. | Active Inactive* |
| Input prompt | Defines whether you are prompted to enter a value. | Active Inactive* |

* Factory setting

6.2.5.3 Weighing

| Parameter | Description | Values |
|---------------------|---|---------------------|
| Weighing profile | Defines the weighing profile. | General* 10d |
| Weight capture mode | Defines the behavior when the button to add the result was tapped, or when the add result was triggered by the automatic weighing result creation. Stable: The balance waits for a stable weight. Immediate: The balance does not wait for a stable weight. | Stable* Immediate |

* Factory setting

6.2.5.4 Automation

| Parameter | Description | Values |
|----------------|---|---|
| Automatic zero | The balance is automatically zeroed when the weight falls below a predefined threshold. This setting is not available for approved balances. | Active Inactive* Numeric The available units depend on the balance model. |

* Factory setting

6.2.5.5 Report






| Parameter | Description | Values |
|---------------------------|---|--|
| Header and Footer | Defines the header and/or the footer to be published. | Title Date/time User Signature Empty lines |
| Balance information | Defines which information about the balance is published. | Balance type Balance ID Balance serial number Software version |
| Quality information | Defines which quality information is published. | Weighing profile Adjustment date/time Routine test information Result state Level state MinWeigh state |
| Result detail information | Defines which information related to the measurement result is published. | Tare / Gross weight Date/time |

6.2.6 Settings: application "Totaling"

This section describes the settings of the application **Totaling**.

☰ **Navigation:** ▼ >  > Σ **Totaling** > Σ *

The settings for this weighing application are grouped as follows:

-  **Main**
-  **ID format**
-  **Weighing**
-  **Automation**
-  **Report**

See also

 Application "Totaling" ▶ Page 49

6.2.6.1 Main

| Parameter | Description | Values |
|-----------|--|--|
| Unit | Defines the unit of the weighing result. | The available units depend on the balance model. |

6.2.6.2 ID format

Sample ID

| Parameter | Description | Values |
|-----------------|--|---------------------------|
| Sample ID | Defines a sample identification. | Active Inactive* |
| Default value | Defines a default value for the sample description. | Text (1...200 characters) |
| Automatic value | Defines whether an automatic value is generated for the sample description. | not editable |
| Input prompt | Defines whether you are prompted to enter a value. If a default value is defined, this parameter is not editable. | Active Inactive* |

* Factory setting

Description

| Parameter | Description | Values |
|-----------------|---|---------------------------|
| Description | Allows to define a sample description. | Active Inactive* |
| Type | Defines the sample type. | Sample* Series |
| Label | Describes the sample. | Text (1...25 characters) |
| Default value | Defines a default value for the sample description. | Text (1...200 characters) |
| Automatic value | Defines whether an automatic value is generated for the sample description. | Active Inactive* |
| Input prompt | Defines whether you are prompted to enter a value. | Active Inactive* |

* Factory setting

6.2.6.3 Weighing

| Parameter | Description | Values |
|---------------------|--|---|
| Weighing profile | Defines the weighing profile. | General* 10d |
| Weight capture mode | Defines the behavior when the button to add the result was tapped, or when the add result was triggered by the automatic weighing result creation. Stable: The balance waits for a stable weight. Immediate: The balance does not wait for a stable weight. Automatic, stable (zero excluded): The results are published as soon as the weight is stable. Values of 0 g are not published. Automatic, stable (zero included): The results are published as soon as the weight is stable. Values of 0 g are also published. | Stable* Immediate Automatic, stable (zero excluded) Automatic, stable (zero included) |

* Factory setting

6.2.6.4 Automation

| Parameter | Description | Values |
|----------------|---|---|
| Automatic zero | The balance is automatically zeroed when the weight falls below a predefined threshold. This setting is not available for approved balances. | Active Inactive* Numeric The available units depend on the balance model. |
| Automatic tare | The balance automatically stores the first stable weight as the tare weight. | Active Inactive* |

| | | |
|-------------|---|---|
| Preset tare | A fixed tare weight can be defined manually or by weighing. | Active Inactive* Numeric The available units depend on the balance model. |
|-------------|---|---|

* Factory setting

6.2.6.5 Report






| Parameter | Description | Values |
|---------------------------|---|--|
| Header and Footer | Defines the header and/or the footer to be published. | Title Date/time User Signature Empty lines |
| Balance information | Defines which information about the balance is published. | Balance type Balance ID Balance serial number Software version |
| Quality information | Defines which quality information is published. | Weighing profile Adjustment date/time Routine test information Result state Level state MinWeigh state |
| Result detail information | Defines which information related to the measurement result is published. | Tare / Gross weight Date/time |

6.2.7 Settings: application "Back weighing"

This section describes the settings of the application **Back weighing**.

☰ **Navigation:** ▼ >  >  **Back weighing** > 

The settings for this weighing application are grouped as follows:

-  **Main**
-  **ID format**
-  **Weighing**
-  **Automation**
-  **Report**

See also

 Application "Back weighing" ▶ Page 50

6.2.7.1 Main

| Parameter | Description | Values |
|-----------------|---|--|
| Tare container | Defines whether a tare container is used. | Active* Inactive |
| Difference unit | Selects the result view for the calculated difference. Percentage (%): Reports the difference between back-weighing and initial weighing as a percentage of the initial weight. Absolute percentage (% Abs.): Reports back-weighing as a percentage of the initial weight. ATRO moisture content (%AM): Reports the moisture content of the sample as a percentage of the dry weight. ATRO dry content (%AD): Reports the wet weight of the sample as a percentage of the dry weight. | Weight* Percentage (%) Absolute percentage (% Abs.) ATRO moisture content (%AM) ATRO dry content (%AD) |

| | | |
|------------------|---|-------------------------------------|
| Difference value | Displays the calculated difference in work area and result view. Unsigned (absolute value): Displays the absolute value. Signed: Displays the value by means of algebraic sign. | Unsigned (absolute value)* Signed |
|------------------|---|-------------------------------------|

* Factory setting

Initial values for weighing

| Parameter | Description | Values |
|-----------|--|--|
| Unit | Defines the unit of the weighing result. | The available units depend on the balance model. |

* Factory setting

6.2.7.2 ID format

Sample ID

| Parameter | Description | Values |
|-----------------|--|---------------------------|
| Sample ID | Defines a sample identification. | Active Inactive* |
| Default value | Defines a default value for the sample description. | Text (1...200 characters) |
| Automatic value | Defines whether an automatic value is generated for the sample description. | not editable |
| Input prompt | Defines whether you are prompted to enter a value. If a default value is defined, this parameter is not editable. | Active Inactive* |

* Factory setting

Description

| Parameter | Description | Values |
|-----------------|---|---------------------------|
| Description | Allows to define a sample description. | Active Inactive* |
| Type | Defines the sample type. | Sample* Series |
| Label | Describes the sample. | Text (1...25 characters) |
| Default value | Defines a default value for the sample description. | Text (1...200 characters) |
| Automatic value | Defines whether an automatic value is generated for the sample description. | Active Inactive* |
| Input prompt | Defines whether you are prompted to enter a value. | Active Inactive* |

* Factory setting

6.2.7.3 Weighing

| Parameter | Description | Values |
|---------------------|--|---|
| Weighing profile | Defines the weighing profile. | General* 10d |
| Weight capture mode | Defines the behavior when the button to add the result was tapped, or when the add result was triggered by the automatic weighing result creation. Stable: The balance waits for a stable weight. Immediate: The balance does not wait for a stable weight. Automatic, stable (zero excluded): The results are published as soon as the weight is stable. Values of 0 g are not published. Automatic, stable (zero included): The results are published as soon as the weight is stable. Values of 0 g are also published. | Stable Immediate Automatic, stable (zero excluded)* Automatic, stable (zero included) |

* Factory setting

Series / Statistics

| Parameter | Description | Values |
|--------------------------|--|-----------------------------------|
| Measurement series | A measurement series can be performed. | Active Inactive* |
| Statistical calculations | Statistical information is provided. This setting is only available if the parameter Measurement series is activated. | Active Inactive* |
| Acceptance range | Defines the acceptance range for the statistical calculations. This setting is only available if the parameter Statistical calculations is activated. | Active Inactive* Numeric (%) |

* Factory setting

6.2.7.4 Automation

| Parameter | Description | Values |
|----------------|---|---|
| Automatic zero | The balance is automatically zeroed when the weight falls below a predefined threshold. This setting is not available for approved balances. | Active Inactive* Numeric The available units depend on the balance model. |
| Preset tare | A fixed tare weight can be defined manually or by weighing. | Active Inactive* Numeric The available units depend on the balance model. |

* Factory setting

6.2.7.5 Report

| Parameter | Description | Values |
|---------------------------|---|--|
| Header and Footer | Defines the header and/or the footer to be published. | Title Date/time User Signature Empty lines |
| Balance information | Defines which information about the balance is published. | Balance type Balance ID Balance serial number Software version |
| Quality information | Defines which quality information is published. | Weighing profile Adjustment date/time Routine test information Result state Level state MinWeigh state |
| Task information | Defines which information about the task is published. | Application settings Measurement details |
| Result detail information | Defines which information related to the measurement result is published. | Tare / Gross weight Date/time |

6.2.8 Settings: application "Density"

This section describes the settings of the application **Density**.

Navigation: > > Density >

The settings for this weighing application are grouped as follows:

- Main
- ID format
- Weighing
- Report

See also

Application "Density" ▶ Page 50

6.2.8.1 Main

| Parameter | Description | Values |
|--------------------|---|---|
| Determination type | Defines the type of density measurement. | Solid* |
| Density result | Defines the number of decimal places of the result value. | 1 decimal place 2 decimal places 3 decimal places 4 decimal places 5 decimal places |

* Factory setting

Initial values for weighing

| Parameter | Description | Values |
|--------------------------|---|--|
| Unit | Defines the unit of the weighing result. | The available units depend on the balance model. |
| Auxiliary liquid | Defines the type of auxiliary liquid used. | Distilled water* Custom |
| Temperature | Defines the temperature of the auxiliary liquid. | Numeric (°C) |
| Auxiliary liquid density | Defines the density of the auxiliary liquid. For distilled water, the value is predefined. | Numeric (g/cm ³) |

* Factory setting

6.2.8.2 ID format

Sample ID

| Parameter | Description | Values |
|-----------------|--|---------------------------|
| Sample ID | Defines a sample identification. | Active Inactive* |
| Default value | Defines a default value for the sample description. | Text (1...200 characters) |
| Automatic value | Defines whether an automatic value is generated for the sample description. | not editable |
| Input prompt | Defines whether you are prompted to enter a value. If a default value is defined, this parameter is not editable. | Active Inactive* |

* Factory setting

Description

| Parameter | Description | Values |
|-------------|--|--------------------------|
| Description | Allows to define a sample description. | Active Inactive* |
| Type | Defines the sample type. | Sample* Series |
| Label | Describes the sample. | Text (1...25 characters) |

| | | |
|-----------------|---|---------------------------|
| Default value | Defines a default value for the sample description. | Text (1...200 characters) |
| Automatic value | Defines whether an automatic value is generated for the sample description. | Active Inactive* |
| Input prompt | Defines whether you are prompted to enter a value. | Active Inactive* |

* Factory setting

6.2.8.3 Weighing

| Parameter | Description | Values |
|---------------------|---|---------------------|
| Weighing profile | Defines the weighing profile. | General* 10d |
| Weight capture mode | Defines the behavior when the button to add the result was tapped, or when the add result was triggered by the automatic weighing result creation. Stable: The balance waits for a stable weight. Immediate: The balance does not wait for a stable weight. | Stable* Immediate |

* Factory setting

Series / Statistics

| Parameter | Description | Values |
|--------------------------|--|-----------------------------------|
| Measurement series | A measurement series can be performed. | Active Inactive* |
| Statistical calculations | Statistical information is provided. This setting is only available if the parameter Measurement series is activated. | Active Inactive* |
| Acceptance range | Defines the acceptance range for the statistical calculations. This setting is only available if the parameter Statistical calculations is activated. | Active Inactive* Numeric (%) |

* Factory setting

6.2.8.4 Report

| Parameter | Description | Values |
|---------------------------|---|--|
| Header and Footer | Defines the header and/or the footer to be published. | Title Date/time User Signature Empty lines |
| Balance information | Defines which information about the balance is published. | Balance type Balance ID Balance serial number Software version |
| Quality information | Defines which quality information is published. | Weighing profile Adjustment date/time Routine test information Result state Level state MinWeigh state |
| Task information | Defines which information about the task is published. | Application settings Measurement details |
| Result detail information | Defines which information related to the measurement result is published. | Tare / Gross weight Date/time |

6.2.9 Settings: application "Differential weighing"

This section describes the settings of the application **Differential weighing**.

Navigation: > > **Differential weighing** >

The settings for this weighing application are grouped as follows:

- **Main**
- **ID format**
- **Weighing**
- **Automation**
- **Report**

See also

Application "Differential weighing" ▶ Page 51

6.2.9.1 Main

| Parameter | Description | Values |
|-------------------|---|--|
| Samples | Defines the number of samples. | Numeric (10* 1...200) |
| Back-weighings | Defines the number of back-weighings per sample. | 1 per sample* 2 per sample 3 per sample |
| Weighing sequence | Defines the sequence of the weighings. | Initial weights first* Sample by sample |
| Spill correction | Corrects the weighing result if sample was spilled. | Active Inactive* |
| Tare container | Defines whether a tare container is used. | Active* Inactive |
| Difference unit | Selects the result view for the calculated difference. Percentage (%): Reports the difference between back-weighing and initial weighing as a percentage of the initial weight. Absolute percentage (% Abs.): Reports back-weighing as a percentage of the initial weight. ATRO moisture content (%AM): Reports the moisture content of the sample as a percentage of the dry weight. ATRO dry content (%AD): Reports the wet weight of the sample as a percentage of the dry weight. | Weight* Percentage (%) Absolute percentage (% Abs.) ATRO moisture content (%AM) ATRO dry content (%AD) |
| Difference value | Displays the calculated difference in work area and result view. Unsigned (absolute value): Displays the absolute value. Signed: Displays the value by means of algebraic sign. | Unsigned (absolute value)* Signed |

* Factory setting

Initial values for weighing

| Parameter | Description | Values |
|-----------|--|--|
| Unit | Defines the unit of the weighing result. | The available units depend on the balance model. |

6.2.9.2 ID format

Sample ID

| Parameter | Description | Values |
|-----------------|---|---------------------------|
| Sample ID | Defines a sample identification. | Active Inactive* |
| Default value | Defines a default value for the sample description. | Text (1...200 characters) |
| Automatic value | Defines whether an automatic value is generated for the sample description. | not editable |

| | | |
|--------------|--|--------------------|
| Input prompt | Defines whether you are prompted to enter a value. If a default value is defined, this parameter is not editable. | Active Inactive* |
|--------------|--|--------------------|

* Factory setting

Description

| Parameter | Description | Values |
|-----------------|---|---------------------------|
| Description | Allows to define a sample description. | Active Inactive* |
| Type | Defines the sample type. | Sample* Series |
| Label | Describes the sample. | Text (1...25 characters) |
| Default value | Defines a default value for the sample description. | Text (1...200 characters) |
| Automatic value | Defines whether an automatic value is generated for the sample description. | Active Inactive* |
| Input prompt | Defines whether you are prompted to enter a value. | Active Inactive* |

* Factory setting

6.2.9.3 Weighing

| Parameter | Description | Values |
|---------------------|--|---|
| Weighing profile | Defines the weighing profile. | General* 10d |
| Weight capture mode | Defines the behavior when the button to add the result was tapped, or when the add result was triggered by the automatic weighing result creation. Stable: The balance waits for a stable weight. Immediate: The balance does not wait for a stable weight. Automatic, stable (zero excluded): The results are published as soon as the weight is stable. Values of 0 g are not published. Automatic, stable (zero included): The results are published as soon as the weight is stable. Values of 0 g are also published. | Stable Immediate Automatic, stable (zero excluded)* Automatic, stable (zero included) |

* Factory setting

Series / Statistics

| Parameter | Description | Values |
|--------------------------|--|-----------------------------------|
| Statistical calculations | Statistical information is provided. This setting is only available if the parameter Measurement series is activated. | Active Inactive* |
| Acceptance range | Defines the acceptance range for the statistical calculations. This setting is only available if the parameter Statistical calculations is activated. | Active Inactive* Numeric (%) |

* Factory setting

6.2.9.4 Automation

| Parameter | Description | Values |
|----------------|---|---|
| Automatic zero | The balance is automatically zeroed when the weight falls below a predefined threshold. This setting is not available for approved balances. | Active Inactive* Numeric The available units depend on the balance model. |

| | | |
|-------------|---|---|
| Preset tare | A fixed tare weight can be defined manually or by weighing. | Active Inactive* Numeric The available units depend on the balance model. |
|-------------|---|---|

* Factory setting

6.2.9.5 Report

| Parameter | Description | Values |
|---------------------------|---|--|
| Header and Footer | Defines the header and/or the footer to be published. | Title Date/time User Signature Empty lines |
| Balance information | Defines which information about the balance is published. | Balance type Balance ID Balance serial number Software version |
| Quality information | Defines which quality information is published. | Weighing profile Adjustment date/time Routine test information Result state Level state MinWeigh state |
| Task information | Defines which information about the task is published. | Application settings Measurement details |
| Result detail information | Defines which information related to the measurement result is published. | Tare / Gross weight Date/time |

6.3 Adjustment settings

See also

[Editing an adjustment](#) ▶ Page 52

6.3.1 Settings: Adjustment strategy

Navigation: Applications > Adjustments > inactive adjustment

| Parameter | Description | Values |
|---------------------|--|--|
| Adjustment strategy | Defines the type of adjustment to be performed. For approved balances, this setting is not available. | No adjustment Internal adjustment* External adjustment |

* Factory setting

See also

[Adjustment strategy](#) ▶ Page 52

6.3.2 Settings: Internal adjustment

Navigation: Applications > Adjustments > Internal > *

The settings are divided into the following subsections:

- Specification**
- Management**
- Report**

Specification

| Parameter | Description | Values |
|------------------|---|---|
| Weighing profile | Defines the weighing profile. | General* 10d |
| 'As found' test | At the start of the adjustment sequence, an internal sensitivity test is automatically performed to evaluate the current status. The test results are displayed. | Active Inactive* |
| 'As left' test | When the adjustment is complete, an internal sensitivity test is automatically performed. The test results are displayed. | Active Inactive* |
| Control limit | Defines the error tolerance of a process with respect to its set value. Exceeding the value Control limit is a violation of quality requirements and therefore requires a correction of the process. If the value Control limit is exceeded: The adjustment failed, the balance is out of specification. This setting is only available if the settings ' As found ' test or ' As left ' test are active. | Numeric (0.1%* 0.001...100%) |
| Warning limit | Defines the upper or lower limit that, if exceeded or not reached, makes more stringent process monitoring necessary. The value Warning limit must be smaller than the value Control limit . Result if the value Warning limit is exceeded: The adjustment is passed, but the difference is higher than expected. This setting is only available if the settings ' As found ' test or ' As left ' test are active. | Active Inactive* Numeric (0.001...100%) |

* Factory setting

Management

| Parameter | Description | Values |
|-------------------------|--|---|
| Error management | Blocks the balance automatically when the adjustment fails. | Active Inactive* |
| Planning - Start events | Serves to plan after which event an adjustment is automatically executed. Multiple selections are allowed. | Active* Inactive <ul style="list-style-type: none"> Start after temperature change Start after leveling Start after power-on |
| Planning - Schedule | Serves to plan at what time and on which weekday an adjustment is automatically executed. <ul style="list-style-type: none"> Start time: Up to three start times can be defined. Preferred days: Monday, Tuesday, ... Sunday | Active* Inactive Numeric |

* Factory setting

Report

| Parameter | Description | Values |
|---------------------|--|--|
| Header and Footer | Defines the header and/or the footer to be published. | Title Date/time User Signature Empty lines |
| Balance information | Defines which information about the balance is published. | Balance type Balance ID Balance serial number Software version |
| Quality information | Defines whether the leveling status of the balance is published. | Level state |

6.3.3 Settings: External adjustment

Navigation: > Applications > Adjustments > External >

The settings are divided into the following subsections:

- Specification
- Report

Specification

| Parameter | Description | Values |
|------------------|-------------------------------|----------------|
| Weighing profile | Defines the weighing profile. | General* 10d |

* Factory setting

Report

| Parameter | Description | Values |
|---------------------|--|--|
| Header and Footer | Defines the header and/or the footer to be published. | Title Date/time User Signature Empty lines |
| Balance information | Defines which information about the balance is published. | Balance type Balance ID Balance serial number Software version |
| Quality information | Defines whether the leveling status of the balance is published. | Level state |

6.4 Test settings

Navigation: > Applications > Tests

The section **Tests** is divided into the following subsections:

- Sensitivity
- Repeatability
- Eccentricity

6.4.1 Settings: Sensitivity test

Navigation: > Applications > Tests > Sensitivity >

The settings are divided into the following subsections:

- Specification
- Management
- Report

Specification

| Parameter | Description | Values |
|------------------|---|--------------------|
| Weighing profile | Defines the weighing profile. | General* 10d |
| Tare container | Defines whether a tare container is used. | Active Inactive* |

* Factory setting

Test Point

Up to two test points can be defined.

| Parameter | Description | Values |
|----------------|---|--|
| Nominal weight | Defines the nominal value of the weight that is used for the test. | Numeric The available units depend on the balance model. |
| Weight class | Defines the weight class according to OIML or ASTM. Alternatively, create a customized tolerance class with the parameter Own . | E1 E2 F1 F2* M1 M2 M3 ASTM000 ASTM00 ASTM0 ASTM1 ASTM2 ASTM3 ASTM4 ASTM5 ASTM6 ASTM7 Own |
| Control limit | Defines the error tolerance of a process with respect to its set value. Exceeding the value Control limit is a violation of quality requirements and therefore requires a correction of the process. If the value Control limit is exceeded: The test failed, the balance is out of specification. | Numeric The available units depend on the balance model. |
| Warning limit | Defines the upper or lower limit that, if exceeded or not reached, makes more stringent process monitoring necessary. The value Warning limit must be smaller than the value Control limit . Result if the value Warning limit is exceeded: The test is passed, but the difference is higher than expected. | Active Inactive* Numeric |

* Factory setting

Management

| Parameter | Description | Values |
|------------------|---|--------------------|
| Error management | Blocks the balance automatically when the adjustment fails. | Active Inactive* |
| Planning | Serves to plan when a test is automatically executed. | Active Inactive* |

If the parameter **Planning** is activated, the following options are available.

| Parameter | Description | Values |
|--------------|---|--|
| Frequency | Serves to plan how often a test is automatically executed. | Daily* Weekly Biweekly Monthly Bimonthly Quarterly Twice a year Yearly |
| Time | Serves to plan at what time a test is automatically executed. | Numeric |
| Notification | Defines how early you will be notified about a planned test (in hours). | Numeric |

* Factory setting

Report

| Parameter | Description | Values |
|---------------------|---|--|
| Header and Footer | Defines the header and/or the footer to be published. | Title Date/time User Signature Empty lines |
| Balance information | Defines which information about the balance is published. | Balance type Balance ID Balance serial number Software version |

| | | |
|---------------------|--|-------------|
| Quality information | Defines whether the leveling status of the balance is published. | Level state |
|---------------------|--|-------------|

See also

 Sensitivity test ▶ Page 54

6.4.2 Settings: Repeatability test

≡ Navigation: ▼ >  Applications >  Tests >  Repeatability > .

The settings are divided into the following subsections:

-  **Specification**
-  **Management**
-  **Report**

Specification

| Parameter | Description | Values |
|-----------------------|--|------------------------|
| Weighing profile | Defines the weighing profile. | General* 10d |
| Number of repetitions | Defines the number of weight measurements of a series. | Numeric (10* 4...20) |
| Tare container | Defines whether a tare container is used. | Active Inactive* |

* Factory setting

Test Point

| Parameter | Description | Values |
|----------------|---|--|
| Nominal weight | Defines the nominal value of the weight that is used for the test. | Numeric The available units depend on the balance model. |
| Weight class | Defines the weight class according to OIML or ASTM. Alternatively, create a customized tolerance class with the parameter Own . | E1 E2 F1 F2* M1 M2 M3 ASTM000 ASTM00 ASTM0 ASTM1 ASTM2 ASTM3 ASTM4 ASTM5 ASTM6 ASTM7 Own |
| Control limit | Defines the error tolerance of a process with respect to its set value. Exceeding the value Control limit is a violation of quality requirements and therefore requires a correction of the process. If the value Control limit is exceeded: The test failed, the balance is out of specification. | Numeric The available units depend on the balance model. |
| Warning limit | Defines the upper or lower limit that, if exceeded or not reached, makes more stringent process monitoring necessary. The value Warning limit must be smaller than the value Control limit . Result if the value Warning limit is exceeded: The test is passed, but the difference is higher than expected. | Active Inactive* Numeric |

* Factory setting

Management

| Parameter | Description | Values |
|------------------|---|--------------------|
| Error management | Blocks the balance automatically when the adjustment fails. | Active Inactive* |

| | | |
|----------|---|--------------------|
| Planning | Serves to plan when a test is automatically executed. | Active Inactive* |
|----------|---|--------------------|

If the parameter **Planning** is activated, the following options are available.

| Parameter | Description | Values |
|--------------|---|--|
| Frequency | Serves to plan how often a test is automatically executed. | Daily* Weekly Biweekly Monthly Bimonthly Quarterly Twice a year Yearly |
| Time | Serves to plan at what time a test is automatically executed. | Numeric |
| Notification | Defines how early you will be notified about a planned test (in hours). | Numeric |

* Factory setting

Report

| Parameter | Description | Values |
|---------------------|--|--|
| Header and Footer | Defines the header and/or the footer to be published. | Title Date/time User Signature Empty lines |
| Balance information | Defines which information about the balance is published. | Balance type Balance ID Balance serial number Software version |
| Quality information | Defines whether the leveling status of the balance is published. | Level state |

See also

[Repeatability test](#) ▶ Page 55

6.4.3 Settings: Eccentricity test

Navigation: > Applications > Tests > Eccentricity >

The settings are divided into the following subsections:

- Specification**
- Management**
- Report**

Specification

| Parameter | Description | Values |
|------------------|-------------------------------|----------------|
| Weighing profile | Defines the weighing profile. | General* 10d |

* Factory setting

Test Point

| Parameter | Description | Values |
|----------------|--|---|
| Nominal weight | Defines the nominal value of the weight that is used for the test. | Numeric The available units depend on the balance model. |

| | | |
|---------------|---|--|
| Weight class | Defines the weight class according to OIML or ASTM. Alternatively, create a customized tolerance class with the parameter Own . | E1 E2 F1 F2* M1 M2 M3 ASTM000 ASTM00 ASTM0 ASTM1 ASTM2 ASTM3 ASTM4 ASTM5 ASTM6 ASTM7 Own |
| Control limit | Defines the error tolerance of a process with respect to its set value. Exceeding the value Control limit is a violation of quality requirements and therefore requires a correction of the process. If the value Control limit is exceeded: The test failed, the balance is out of specification. | Numeric The available units depend on the balance model. |
| Warning limit | Defines the upper or lower limit that, if exceeded or not reached, makes more stringent process monitoring necessary. The value Warning limit must be smaller than the value Control limit . Result if the value Warning limit is exceeded: The test is passed, but the difference is higher than expected. | Active Inactive* Numeric |

* Factory setting

Management

| Parameter | Description | Values |
|------------------|---|--------------------|
| Error management | Blocks the balance automatically when the adjustment fails. | Active Inactive* |
| Planning | Serves to plan when a test is automatically executed. | Active Inactive* |

If the parameter **Planning** is activated, the following options are available.

| Parameter | Description | Values |
|--------------|---|--|
| Frequency | Serves to plan how often a test is automatically executed. | Daily* Weekly Biweekly Monthly Bimonthly Quarterly Twice a year Yearly |
| Time | Serves to plan at what time a test is automatically executed. | Numeric |
| Notification | Defines how early you will be notified about a planned test (in hours). | Numeric |

* Factory setting

Report

| Parameter | Description | Values |
|---------------------|--|--|
| Header and Footer | Defines the header and/or the footer to be published. | Title Date/time User Signature Empty lines |
| Balance information | Defines which information about the balance is published. | Balance type Balance ID Balance serial number Software version |
| Quality information | Defines whether the leveling status of the balance is published. | Level state |

See also

[Eccentricity test](#) ▶ Page 55

7 Maintenance

To guarantee the functionality of the balance and the accuracy of the weighing results, a number of maintenance actions must be performed by the user.

7.1 Maintenance tasks

| Maintenance action | Recommended interval | Remarks |
|---|--|-----------------------|
| Performing an adjustment | <ul style="list-style-type: none">• Daily• After cleaning• After leveling• After changing the location | see "Adjustments" |
| Performing routine tests (eccentricity test, repeatability test, sensitivity test). METTLER TOLEDO recommends to at least perform a sensitivity test. | <ul style="list-style-type: none">• After cleaning• After assembling the balance• After a software update• Depending on your internal regulations (SOP) | see "Tests" |
| Cleaning | <ul style="list-style-type: none">• After every use• Depending on the degree of pollution• Depending on your internal regulations (SOP) | see "Cleaning" |
| Updating the software | <ul style="list-style-type: none">• Depending on your internal regulations (SOP).• After a new software release. | see "Software update" |

See also

- [Adjustments](#) ▶ Page 52
- [Tests](#) ▶ Page 54
- [Cleaning](#) ▶ Page 113
- [Software update](#) ▶ Page 118

7.2 Cleaning

7.2.1 Disassembling for cleaning

Note

Depending on the balance model, the components may look different.

Note

In most cases, it is not necessary to remove the protective covers to clean the balance.

7.2.1.1 Balances with draft shield




CAUTION

Injury due to sharp objects or broken glass


Instrument components, e.g., glass, can break and lead to injuries.

- Always proceed with focus and care.


- 1 Fully open the side door (1).
- 2 Remove the weighing pan (2).
- 3 Only for balances with a readability of 0.01 mg:
Remove the draft-protection element (3).
- 4 Remove the drip tray (4).



- 5 Lift up the QuickLock (5) and pull the side door (1) towards the back to remove it (right, left).
- NOTICE: Damage to the instrument**
Hold the side door (1) tightly when removing it.




- 6 Turn the QuickLock (6, right, left), tilt the front panel (7) towards the front and lift it upwards to remove it.




- 7 Pull the top door (8) towards the front to remove it.

Note


Optional, if required: Remove the protective covers for cleaning as described below.




- 8 Remove the protective cover (9) from the terminal (10).



- 9 Open the QuickLock (11) and remove the draft shield (12).




- 10 Remove the protective cover (13) from the platform (14).




7.2.1.2 Balances without draft shield

- 1 Remove the weighing pan (1).
- 2 Remove the weighing pan support (2) and/or the support caps (3) (if applicable).
- 3 Remove the drip tray (4).




- 4 Optional, if required: Remove the protective cover (5) from the terminal (6).
- 5 Optional, if required: Remove the protective cover (7) from the platform (8).



7.2.1.3 Balances, large

- 1 Remove the weighing pan (1).
- 2 Remove the support caps (2).
- 3 Optional, if required: Remove the protective cover (3) from the terminal (4).



7.2.2 Cleaning agents

In the following table, cleaning tools and cleaning agents recommended by METTLER TOLEDO are listed. Pay attention to the concentration of the agents specified in the table.

| | | Tools | | | | | Cleaning agents | | | | | |
|----------------------|----------------------------------|--------------|-------|------------|-------|---------|-----------------|-------------------|---------------------------|-------------------------|-----------------------|--|
| | | Paper tissue | Brush | Dishwasher | Water | Acetone | Ethanol (70%) | Isopropanol (70%) | Hydrochloric acid (3-10%) | Sodium hydroxide (1-4%) | Peracetic acid (2-3%) | |
| Around the balance | Balance housing | ✓ | ✓ | — | ✓ | — | ✓ | ✓ | ✓ | ✓ | ✓ | |
| | Feet | ✓ | ✓ | — | ✓ | — | ✓ | ✓ | ✓ | ✓ | ✓ | |
| Balance terminal | Terminal | ✓ | ✓ | — | ✓ | PR | ✓ | ✓ | ✓ | ✓ | ✓ | |
| | Display | ✓ | ✓ | — | ✓ | PR | ✓ | ✓ | ✓ | ✓ | ✓ | |
| | Terminal cover | ✓ | ✓ | — | ✓ | — | ✓ | ✓ | ✓ | PR | PR | |
| Balance draft shield | Glass panels | ✓ | ✓ | ✓ | ✓ | PR | ✓ | ✓ | ✓ | ✓ | ✓ | |
| | Non-removable handles and frames | ✓ | ✓ | — | ✓ | — | ✓ | ✓ | ✓ | ✓ | ✓ | |
| Weighing area | Weighing pan | ✓ | ✓ | ✓ | ✓ | PR | ✓ | ✓ | ✓ | ✓ | ✓ | |
| | Drip tray | ✓ | ✓ | ✓ | ✓ | PR | ✓ | ✓ | — | — | ✓ | |
| Accessories | Dust cover | ✓ | ✓ | — | ✓ | — | ✓ | ✓ | — | — | PR | |
| | Antistatic kit | ✓ | ✓ | — | — | — | — | — | — | — | — | |

Legend

- ✓ Recommended by METTLER TOLEDO; can be used without limitation.
- PR Partially recommended by METTLER TOLEDO: individual resistance to acid and alkali must be evaluated, including dependence to the time exposure.
- Not recommend. High risk for damage.

7.2.3 Cleaning the balance



NOTICE

Damage to the instrument due to inappropriate cleaning methods

If liquid enters the housing, it can damage the instrument. The surface of the instrument can be damaged by certain cleaning agents, solvents, or abrasives.

- 1 Do not spray or pour liquid on the instrument.
- 2 Only use the cleaning agents specified in the Reference Manual (RM) of the instrument or the guide "8 Steps to a Clean Balance".
- 3 Only use a lightly moistened, lint-free cloth or a tissue to clean the instrument.
- 4 Wipe off any spills immediately.



For further information on cleaning a balance, consult "8 Steps to a Clean Balance".

► www.mt.com/lab-cleaning-guide

Cleaning around the balance

- Remove any dirt or dust around the balance and avoid further contaminations.

Cleaning the terminal

- Clean the terminal with a damp cloth or a tissue and a mild cleaning agent.

Cleaning the removable parts

- Clean the removed part with a damp cloth or a tissue and a mild cleaning agent or clean in a dishwasher up to 80 °C.

Cleaning the balance

- 1 Disconnect the balance from the AC/DC adapter.
- 2 Use a lint-free cloth moistened with a mild cleaning agent to clean the surface of the balance.
- 3 Remove powder or dust with a disposable tissue first.
- 4 Remove sticky substances with a damp lint-free cloth and a mild solvent, e.g., isopropanol or ethanol 70%.

7.2.4 Putting into operation after cleaning

- 1 Reassemble the balance.
- 2 Check that the draft shield doors (top, sides) open and close normally (if applicable).
- 3 Reconnect the balance to the AC/DC adapter.
- 4 Check the level status, level the balance if necessary.
- 5 Respect the warm-up time specified in the "Technical Data".
- 6 Perform an internal adjustment.
- 7 Perform a routine test according to the internal regulations of your company. METTLER TOLEDO recommends performing a sensitivity test after cleaning the balance.
- 8 Press **→0←** to zero the balance.
 - ➔ The balance is ready for use.

See also

- 🔗 [Leveling the balance](#) ► Page 31
- 🔗 [Technical Data](#) ► Page 125
- 🔗 [Performing an internal adjustment](#) ► Page 53

7.3 Service

Regular servicing by an authorized service technician ensures reliability for years to come. Contact your METTLER TOLEDO representative for details about the available service options.

7.4 Software update

Search for software:

► www.mt.com/labweighing-software-download

Contact a METTLER TOLEDO service representative if you need support updating the software.

METTLER TOLEDO recommends saving the data on a storage device before updating the software.

7.4.1 Updating the software

☰ **Navigation:** ☰ **Balance menu** > ☰ **Maintenance** > 📄 **Software update**

The function **Software update** is only available to users with the corresponding rights.



NOTICE

Removing USB storage device during software update

Do not remove the USB storage device during the software update procedure. This can lead to an incomplete or faulty installation of the balance software.

- A USB storage device containing the software installer is connected to the balance.
- 1 Tap 📄 **Software update**.
- 2 Select **Software update**.
- 3 Tap → **Next**.
 - ➔ An update wizard opens and leads you step-by-step through the procedure.
- 4 When prompted, tap **I accept the terms in the license agreement** and confirm with ✓ **OK**.

7.4.2 Putting into operation after software update

- 1 Press ⏻ to switch on the balance.
- 2 Check the level status. Level the balance if required.
- 3 Perform an internal adjustment.
- 4 Perform a routine test according to the internal regulations of your company.
- 5 Press → **0** ← to zero the balance.
 - ➔ The balance is ready for use.

See also

- 🔗 [Leveling the balance](#) ▶ Page 31
- 🔗 [Performing an internal adjustment](#) ▶ Page 53

7.5 Resetting the balance

A reset puts the balance back to factory state. All user data are deleted.

If the setting **User management** is inactive, any user can reset the balance. If the setting **User management** is active, resetting the balance requires corresponding permission.



NOTICE

Reset causes data loss

Resetting the balance will delete user application data and set the user configuration back to factory state.

- **User management** > **General**: The user's configuration permission is activated.
- 1 Tap ☰ **Menu**.
- 2 Tap ☰ **Maintenance**.
- 3 Tap ↻ **Reset**.
- 4 Tap ↻ **Reset** to confirm.
 - ➔ The balance restarts with factory settings.

8 Troubleshooting

Possible errors with their cause and remedy are described in the following chapter. If there are errors that cannot be corrected through these instructions, contact METTLER TOLEDO.

8.1 Error messages

| Error message | Possible cause | Diagnostic | Remedy |
|--|---|--|---|
| The balance shows an error code. | Software or hardware error. | – | Restart the balance. If that does not help, perform a balance reset. If the issue persists, contact your METTLER TOLEDO service representative. |
| Date and time lost | The battery is low. The battery backup is lost. | Check the settings for date and time. | Connect the balance to the power outlet and let the battery charge for two to three days. Set date and time. If the issue persists, contact your METTLER TOLEDO service representative. |
| Communication with draft shield is not possible. Only applies to balances with a backlit draft shield. | The contact pins are dirty or do not fit together. | Check the contact pins between the balance and the draft shield. | Clean or adjust the contact pins. |
| | The draft shield is defective. | – | Contact your METTLER TOLEDO service representative. |
| Communication with weighing module is not possible. | The internal communication does not work properly. | – | Perform a balance reset. Reinstall the balance software. If the issue persists, contact your METTLER TOLEDO service representative. |
| Data memory defect. | EEPROM is corrupt. | – | Perform a balance reset. If the issue persists, contact your METTLER TOLEDO service representative. |
| Memory full. | The memory storage is full. | – | Perform a balance reset. |
| No standard adjustment. | The standard adjustment is missing or invalid. | – | Contact your METTLER TOLEDO service representative. |
| Program memory defect. | The checksum for the stored program is not correct anymore. | – | Reinstall the balance software. If the issue persists, contact your METTLER TOLEDO service representative. |

| Error message | Possible cause | Diagnostic | Remedy |
|---|--|------------|--|
| Temperature sensor defect. | The temperature sensor that measures the cell temperature is defective. | – | Contact your METTLER TOLEDO service representative. |
| Type data damaged. | The TDNR is corrupt. | – | Perform a balance reset. If the issue persists, contact your METTLER TOLEDO service representative. |
| Unexpected startup problem | A problem occurred while starting up the balance. Some data could not be read correctly from the memory. | – | Restart the balance. If the issue persists, contact your METTLER TOLEDO service representative. |
| Unknown error | General error for an unspecific issue. | – | Restart the balance. Perform a balance reset. If the issue persists, contact your METTLER TOLEDO service representative. |
| User data damaged. | The user data is damaged or its context is incorrect. | – | Perform a balance reset. If the issue persists, contact your METTLER TOLEDO service representative. |
| Wrong cell data. | The cell data is damaged or its checksum is incorrect. | – | Contact your METTLER TOLEDO service representative. |
| Wrong legally relevant authentication. Only applies to approved balances. | – | – | Contact your METTLER TOLEDO service representative. |

8.2 Error symptoms

| Error symptom | Possible cause | Diagnostic | Remedy |
|---|--|---|---|
| The balance shows no valid date and time. | The battery is low. The battery backup is lost. | Check the settings for date and time. | Connect the balance to the power outlet and let the battery charge for two to three days. Set date and time. If the issue persists, contact your METTLER TOLEDO service representative. |
| The display is dark. | The balance is on standby or in power-saving mode. | – | Switch on the balance. |
| | There is no power. | Check the connection to the AC/DC adapter and the power outlet. | Connect the balance to the power outlet. See "Connecting the balance". |

| Error symptom | Possible cause | Diagnostic | Remedy |
|---|--|--|---|
| | The wrong AC/DC adapter is connected to the balance. | Check the AC/DC adapter, see "Technical Data". | Use the correct AC/DC adapter. |
| | The AC/DC adapter is defective. | – | Replace the AC/DC adapter. |
| | The display is defective. | – | Contact your METTLER TOLEDO service representative. |
| The balance does not react to any input. | Software freeze. | – | Disconnect the power cable from the balance and reconnect it after a few seconds. Perform a balance reset. If the issue persists, contact your METTLER TOLEDO service representative. |
| The balance does not start up properly. | The balance has no power. | Check if the AC/DC adapter is plugged in. | Connect the AC/DC adapter. |
| | The AC/DC adapter is defective. | Check with another AC/DC adapter if available. | Replace the AC/DC adapter. See "Accessories". |
| The balance does not return to zero when the weight is removed. | Something is touching the weighing pan. Dirt or dust on the weighing pan. | Remove the weighing pan and check it for dirt or dust. | Clean the weighing pan. If the issue persists, contact your METTLER TOLEDO service representative. |
| Taring fails. | The weighing bench is vibrating. | Tab \rightarrow and check if the value on the display is still unstable. | Place the balance on a weighing bench free of vibrations. |
| | The weighing sample is electrostatically charged. | Place a test weight on the weighing pan. Check if the weighing result is stable. | For balances with a draft shield: place a water container into the weighing chamber to increase the humidity. Use an antistatic device. See "Accessories". |
| | The balance is exposed to drafts. | Check the location for sources of draft. | Place the balance in a location without draft. |
| The internal adjustment fails. | A weight is on the weighing pan. | – | Remove the weight from the weighing pan. |
| | Repeatability is poor. | – | Perform a repeatability test. |
| | The internal weight does not function properly. | – | Contact your METTLER TOLEDO service representative. |
| The sensitivity test fails. | A weight is on the weighing pan. | – | Remove the weight from the weighing pan. |
| | The internal weight does not work properly. | – | Contact your METTLER TOLEDO service representative. |

| Error symptom | Possible cause | Diagnostic | Remedy |
|---|---|---|--|
| The repeatability test fails. | Eccentricity is out of tolerance. | Perform an eccentricity test. | If the eccentricity test fails, contact your METTLER TOLEDO service representative. |
| | The environment is unstable. | — | Place the balance in a location with suitable environmental conditions. |
| The display shows overload or underload. | The wrong weighing pan is installed. | Slightly lift or press the weighing pan to see if the weight appears on the display. | Install a proper weighing pan. |
| | No weighing pan is installed. | — | Install a proper weighing pan. |
| | Incorrect zero point when the balance is switched on. | — | Disconnect the power cable and reconnect it after a few seconds. |
| | The balance is not adjusted. | — | Perform an internal adjustment. See "Performing an internal adjustment". |
| The value on the display oscillates. | Vibrations on the weighing bench, for example, building vibrations, foot traffic. | Place a beaker with water on the weighing bench. Vibrations cause ripples on the water surface. | Protect the weighing location against vibrations, for example, with an absorber. Find a different weighing location. |
| | Draft due to untight draft shield and/or open window. | Check the draft shield for gaps. | Fix the draft shield. Close the window. |
| | The weighing sample is electrostatically charged. | Check if the weighing result is stable when using a test weight. | Increase the air humidity in the weighing chamber. Use an ionizer. See "Accessories". |
| | The location is not suitable for weighing. | — | Follow the requirements for the location. See "Selecting the location". |
| | Something is touching the weighing pan. | Check for touching parts or dirt. | Remove touching parts. Clean the balance. |
| The value on the display is drifting towards plus or minus. | The location is not suitable for weighing. | — | Place the balance in a location with suitable environmental conditions. |
| | The weighing sample absorbs moisture or evaporates moisture. | Check if the weighing result is stable when using a test weight. | Cover the weighing sample. |
| | The weighing sample is electrostatically charged. | Use a test weight to check if the weighing result is stable. | Increase the humidity in the weighing chamber. Use an ionizer. See "Accessories". |

| Error symptom | Possible cause | Diagnostic | Remedy |
|---|---|--|--|
| | The weighing sample is warmer or colder than the air in the weighing chamber. | Check if the weighing result is stable when using an acclimatized test weight. | Bring the sample to room temperature. |
| | The balance has not yet warmed up. | – | Let the balance warm up. Adequate warm-up time is specified in the section "General data". |
| Drop to cursor: data transfer does not work properly | The lock of the number pad on the keyboard is activated. | The format of the transferred data is wrong. | Switch off the lock of the number pad. |
| | Asian IME (input method editor) is running. | The format of the transferred data is wrong. | Switch off IME. |

8.3 Saving a support file

When requesting help from your METTLER TOLEDO service representative, you may be asked to send a support file. This file is analyzed and can help to solve issues with the balance.

≡ **Navigation:** ≡ **Balance menu** > ≡ **Maintenance**

- The section ≡ **Maintenance** is open.
- A USB storage device is available.
- 1 Tap ≡ **Save support file**.
- 2 Connect a USB storage device to the balance.
- 3 Tap ✓ **OK**.
 - ➔ The support file is saved to the USB storage device.

8.4 Putting into operation after fixing an error


After troubleshooting, perform the following steps to put the balance into operation:

- Ensure that the balance is completely reassembled and cleaned.
- Reconnect the balance to the AC/DC adapter.


9 Technical Data

9.1 General data

Power supply for balances with a readability of 0.01 mg and 0.1 mg

| | |
|----------------------------|---|
| AC/DC adapter: | Input: 100 – 240 V AC \pm 10%, 50 – 60 Hz, 0.8 A, 61 – 80 VA Output: 12 V DC, 2.5 A, LPS |
| Cable for AC/DC adapter: | 3-core, with country-specific plug |
| Balance power consumption: | 12 V DC, 1.0 A |
| Polarity: |  |

Power supply for balances with a readability of 1 mg or higher

| | |
|----------------------------|---|
| AC/DC adapter: | Input: 100 – 240 V AC \pm 10%, 50 – 60 Hz, 0.5 A Output: 12 V DC, 1.5 A, LPS |
| Balance power consumption: | 12 V DC, 1.0 A |
| Polarity: |  |

Protection and standards

| | |
|--------------------------|--|
| Overvoltage category: | II |
| Degree of pollution: | 2 |
| Ingress protection code: | IP41 (balances with readability of 0.1 mg or 1 mg) IP54 (balances with readability of 0.01 g or higher) |

Note

Stated IP is only achieved when the balance is ready for operation. The protective covers must be installed, and the caps must cover the interface connections.

| | |
|-------------------------------|-----------------------------------|
| Standards for safety and EMC: | See Declaration of Conformity |
| Range of application: | Use only indoors in dry locations |

Environmental conditions

The limit values apply when the balance is used under the following environmental conditions:

| | |
|------------------------------|--|
| Height above mean sea level: | Up to 5000 m |
| Ambient temperature: | +10 – +30 °C |
| Temperature change, max.: | 5 °C/h |
| Relative humidity: | 30 – 70%, non-condensing |
| Acclimatization time: | Recommendation: Up to 4 hours for precision balances, or up to 8 hours for analytical balances. These values apply after placing the balance in the same location where it will be put into operation. |

Note

The acclimatization time depends on the readability of the balance, and on the environmental conditions.

| | |
|---------------|---|
| Warm-up time: | At least 30 minutes for precision balances, 60 minutes for analytical balances, or 120 minutes for balances with a readability of 0.01 mg. These values apply after connecting the balance to the power supply, or after exiting power-saving mode. When switched on from standby, the balance is ready for operation immediately. |
|---------------|---|

The balance can be used under the following environmental conditions. However, the weighing performances of the balance may be outside the limit values:

| | |
|----------------------|---|
| Ambient temperature: | +5 °C – +40 °C |
| Relative humidity: | 20% to max. 80% at 31 °C, decreasing linearly to 50% at 40 °C, non-condensing |

The balance can be disconnected and stored in its packaging under the following conditions:

| | |
|----------------------|--------------------------|
| Ambient temperature: | -25 – +70 °C |
| Relative humidity: | 10 – 90%, non-condensing |

9.2 Materials

| | |
|---------------------------|--|
| Housing regular balances: | Bottom housing: die-cast aluminum Top housing: die-cast aluminum, powder-coated Housing frame: POM Terminal frame: chrome-plated aluminum |
| Housing large balances: | Die-cast aluminum, powder-coated Housing frame: POM Terminal frame: chrome-plated aluminum |
| Draft shield: | POM (U-shaped top frame, rear QuickLock), PBT (bottom plate), glass (doors, front panel), powder-coated aluminum (posts), PA 12 (handles, front QuickLock) |
| Weighing pan: | Balances with readability of 0.01 mg or 0.1 mg: stainless steel X2CrNiMo17-12-2 (1.4404) Balances with readability of 1 mg: die-cast zinc, chrome-plated Balances with readability of 0.01 g: die-cast aluminum, chrome-plated Balances with readability of 0.1 g: stainless steel X2CrNiMo17-12-2 (1.4404) Large balances: stainless steel X5CrNi18-10 (1.4301) |
| Draft-protection element: | Balances with readability of 0.01 mg: stainless steel X2CrNiMo17-12-2 (1.4404) |
| Drip tray: | Die-cast aluminum, powder-coated |
| Touchscreen: | Glass |
| Protective cover: | PET |
| Feet: | TPE, stainless steel X5CrNi18-10 (1.4301) |
| Battery: | ML2032 |


9.3 Explanatory notes for the METTLER TOLEDO AC/DC adapter

The certified external AC/DC adapter complies to the requirements for Class II double insulated equipment. It is not provided with a protective earth connection but with a functional earth connection for EMC purposes. This earth connection **is not** a safety feature. Further information about the compliance of our products can be found in the "Declaration of Conformity" delivered with every product.

In case of testing with regard to the European Directive 2001/95/EC, the AC/DC adapter and the instrument have to be handled as Class II double insulated equipment.

Consequently, a grounding test is not required. It is not necessary to carry out a grounding test between the earth connector of the power plug and any exposed part of the metallic housing of the instrument.

Because the instrument is sensitive to static charges, a leakage resistor of 10 k Ω is connected between the earth connector (1) and the negative pole (2) of the AC/DC adapter. The arrangement is shown in the equivalent circuit diagram. This resistor is not part of the electrical safety arrangement and does not require testing at regular intervals.



9.4 Model-specific data

9.4.1 Analytical balances, readability 0.01 mg or 0.1 mg

| | MX105 | MX105DU | MX205DU |
|--|-------------------------------|-------------------------------|--------------------------------|
| Limit values | | | |
| Capacity | 120 g | 120 g | 220 g |
| Nominal load | 100 g | 100 g | 200 g |
| Readability | 0.01 mg | 0.1 mg | 0.1 mg |
| Capacity of fine range | – | 42 g | 82 g |
| Readability in fine range | – | 0.01 mg | 0.01 mg |
| Repeatability (at 5% load) | 0.02 mg | 0.02 mg | 0.02 mg |
| Linearity deviation | 0.1 mg | 0.2 mg | 0.2 mg |
| Eccentricity deviation (at test load) | 0.3 mg (50 g) | 0.3 mg (50 g) | 0.3 mg (100 g) |
| Sensitivity offset (at nominal load) ▲ | 0.4 mg | 0.4 mg | 0.8 mg |
| Sensitivity temperature drift | 0.0002%/°C | 0.0002%/°C | 0.0002%/°C |
| Typical values | | | |
| Repeatability (at 5% load) | 0.0125 mg | 0.0125 mg | 0.0125 mg |
| Linearity deviation | 0.06 mg | 0.06 mg | 0.06 mg |
| Eccentricity deviation (at test load) | 0.1 mg (50 g) | 0.1 mg (50 g) | 0.1 mg (100 g) |
| Sensitivity offset (at nominal load) ▲ | 0.25 mg | 0.25 mg | 0.5 mg |
| Minimum weight (USP, tolerance = 0.10%) ▼ | 25 mg | 25 mg | 25 mg |
| Minimum weight (tolerance = 1%) ▼ | 2.5 mg | 2.5 mg | 2.5 mg |
| Settling time | 2 s | 2 s | 2 s |
| Dimensions and other specifications | | | |
| Balance dimensions (W × D × H) | 222 × 379 × 353 mm | 222 × 379 × 353 mm | 222 × 379 × 353 mm |
| Weighing pan diameter | 80 mm | 80 mm | 80 mm |
| Usable height of draft shield | 235 mm | 235 mm | 235 mm |
| Balance weight | 7.5 kg | 7.5 kg | 7.5 kg |
| Weights for routine testing | | | |
| Weights (OIML class) | 100 g (F2) / 5 g (F2) | 100 g (F2) / 5 g (F2) | 200 g (F2) / 10 g (F2) |
| Weights (ASTM class) | 100 g (ASTM 1) / 5 g (ASTM 1) | 100 g (ASTM 1) / 5 g (ASTM 1) | 200 g (ASTM 1) / 10 g (ASTM 1) |

▲ after adjustment with internal weight

▼ determined at 5% load, k = 2

| | MX104 | MX204 | MX304 |
|--|-------------------------------|--------------------------------|--------------------------------|
| Limit values | | | |
| Capacity | 120 g | 220 g | 320 g |
| Nominal load | 100 g | 200 g | 300 g |
| Readability | 0.1 mg | 0.1 mg | 0.1 mg |
| Capacity of fine range | – | – | – |
| Readability in fine range | – | – | – |
| Repeatability (at 5% load) | 0.1 mg | 0.1 mg | 0.1 mg |
| Linearity deviation | 0.2 mg | 0.2 mg | 0.3 mg |
| Eccentricity deviation (at test load) | 0.4 mg (50 g) | 0.4 mg (100 g) | 0.4 mg (100 g) |
| Sensitivity offset (at nominal load) ▲ | 0.5 mg | 0.8 mg | 1 mg |
| Sensitivity temperature drift | 0.0002%/°C | 0.0002%/°C | 0.0002%/°C |
| Typical values | | | |
| Repeatability (at 5% load) | 0.05 mg | 0.05 mg | 0.05 mg |
| Linearity deviation | 0.06 mg | 0.06 mg | 0.06 mg |
| Eccentricity deviation (at test load) | 0.1 mg (50 g) | 0.1 mg (100 g) | 0.1 mg (100 g) |
| Sensitivity offset (at nominal load) ▲ | 0.3 mg | 0.5 mg | 0.6 mg |
| Minimum weight (USP, tolerance = 0.10%) ▼ | 100 mg | 100 mg | 100 mg |
| Minimum weight (tolerance = 1%) ▼ | 10 mg | 10 mg | 10 mg |
| Settling time | 2 s | 2 s | 2 s |
| Dimensions and other specifications | | | |
| Balance dimensions (W × D × H) | 222 × 379 × 353 mm | 222 × 379 × 353 mm | 222 × 379 × 353 mm |
| Weighing pan diameter | 90 mm | 90 mm | 90 mm |
| Usable height of draft shield | 239 mm | 239 mm | 239 mm |
| Balance weight | 6.2 kg | 6.2 kg | 6.2 kg |
| Weights for routine testing | | | |
| Weights (OIML class) | 100 g (F2) / 5 g (F2) | 200 g (F2) / 10 g (F2) | 200 g (F2) / 10 g (F2) |
| Weights (ASTM class) | 100 g (ASTM 1) / 5 g (ASTM 1) | 200 g (ASTM 1) / 10 g (ASTM 1) | 200 g (ASTM 1) / 10 g (ASTM 1) |

▲ after adjustment with internal weight

▼ determined at 5% load, k = 2

9.4.2 Precision balances, readability 1 mg

| | MX303 | MX603 | MX1203 |
|--|--------------------------------|--------------------------------|---------------------------------|
| Limit values | | | |
| Capacity | 320 g | 620 g | 1.22 kg |
| Nominal load | 300 g | 600 g | 1.2 kg |
| Readability | 1 mg | 1 mg | 1 mg |
| Capacity of fine range | – | – | – |
| Readability in fine range | – | – | – |
| Repeatability (at 5% load) | 0.9 mg | 0.9 mg | 0.9 mg |
| Linearity deviation | 2 mg | 2 mg | 2 mg |
| Eccentricity deviation (at test load) | 3 mg (100 g) | 3 mg (200 g) | 3 mg (500 g) |
| Sensitivity offset (at nominal load) ▲ | 6 mg | 6 mg | 7 mg |
| Sensitivity temperature drift | 0.0002%/°C | 0.0002%/°C | 0.0002%/°C |
| Typical values | | | |
| Repeatability (at 5% load) | 0.5 mg | 0.5 mg | 0.5 mg |
| Linearity deviation | 0.6 mg | 0.6 mg | 0.6 mg |
| Eccentricity deviation (at test load) | 1 mg (100 g) | 1 mg (200 g) | 1 mg (500 g) |
| Sensitivity offset (at nominal load) ▲ | 4 mg | 4 mg | 4 mg |
| Minimum weight (USP, tolerance = 0.10%) ▼ | 1 g | 1 g | 1 g |
| Minimum weight (tolerance = 1%) ▼ | 100 mg | 100 mg | 100 mg |
| Settling time | 1.5 s | 1.5 s | 1.5 s |
| Dimensions and other specifications | | | |
| Balance dimensions (W × D × H) | 222 × 379 × 353 mm | 222 × 379 × 353 mm | 222 × 379 × 353 mm |
| Weighing pan dimensions (W × D) | 127 × 127 mm | 127 × 127 mm | 127 × 127 mm |
| Usable height of draft shield | 238 mm | 238 mm | 238 mm |
| Balance weight | 6.3 kg | 6.3 kg | 6.7 kg |
| Weights for routine testing | | | |
| Weights (OIML class) | 200 g (F2) / 10 g (F2) | 500 g (F2) / 20 g (F2) | 1000 g (F2) / 50 g (F2) |
| Weights (ASTM class) | 200 g (ASTM 1) / 10 g (ASTM 1) | 500 g (ASTM 1) / 20 g (ASTM 1) | 1000 g (ASTM 1) / 50 g (ASTM 1) |

▲ after adjustment with internal weight

▼ determined at 5% load, k = 2

| | MX303N | MX603N | MX1203N |
|--|--------------------------------|--------------------------------|---------------------------------|
| Limit values | | | |
| Capacity | 320 g | 620 g | 1.22 kg |
| Nominal load | 300 g | 600 g | 1.2 kg |
| Readability | 1 mg | 1 mg | 1 mg |
| Capacity of fine range | – | – | – |
| Readability in fine range | – | – | – |
| Repeatability (at 5% load) | 0.9 mg | 0.9 mg | 0.9 mg |
| Linearity deviation | 2 mg | 2 mg | 2 mg |
| Eccentricity deviation (at test load) | 3 mg (100 g) | 3 mg (200 g) | 3 mg (500 g) |
| Sensitivity offset (at nominal load) ▲ | 6 mg | 6 mg | 7 mg |
| Sensitivity temperature drift | 0.0002%/°C | 0.0002%/°C | 0.0002%/°C |
| Typical values | | | |
| Repeatability (at 5% load) | 0.5 mg | 0.5 mg | 0.5 mg |
| Linearity deviation | 0.6 mg | 0.6 mg | 0.6 mg |
| Eccentricity deviation (at test load) | 1 mg (100 g) | 1 mg (200 g) | 1 mg (500 g) |
| Sensitivity offset (at nominal load) ▲ | 4 mg | 4 mg | 4 mg |
| Minimum weight (USP, tolerance = 0.10%) ▼ | 1 g | 1 g | 1 g |
| Minimum weight (tolerance = 1%) ▼ | 100 mg | 100 mg | 100 mg |
| Settling time | 1.5 s | 1.5 s | 1.5 s |
| Dimensions and other specifications | | | |
| Balance dimensions (W × D × H) | 194 × 379 × 100 mm | 194 × 379 × 100 mm | 194 × 379 × 100 mm |
| Weighing pan dimensions (W × D) | 127 × 127 mm | 127 × 127 mm | 127 × 127 mm |
| Usable height of draft shield | – | – | – |
| Balance weight | 4 kg | 4 kg | 4.4 kg |
| Weights for routine testing | | | |
| Weights (OIML class) | 200 g (F2) / 10 g (F2) | 500 g (F2) / 20 g (F2) | 1000 g (F2) / 50 g (F2) |
| Weights (ASTM class) | 200 g (ASTM 1) / 10 g (ASTM 1) | 500 g (ASTM 1) / 20 g (ASTM 1) | 1000 g (ASTM 1) / 50 g (ASTM 1) |

▲ after adjustment with internal weight

▼ determined at 5% load, k = 2

9.4.3 Precision balances, readability 0.01 g or 0.1 g

| | MX2002 | MX4002 | MX6002 |
|--|-------------------------------------|-------------------------------------|-------------------------------------|
| Limit values | | | |
| Capacity | 2.2 kg | 4.2 kg | 6.2 kg |
| Nominal load | 2 kg | 4 kg | 6 kg |
| Readability | 0.01 g | 0.01 g | 0.01 g |
| Capacity of fine range | – | – | – |
| Readability in fine range | – | – | – |
| Repeatability (at 5% load) | 8 mg | 8 mg | 8 mg |
| Linearity deviation | 20 mg | 20 mg | 20 mg |
| Eccentricity deviation (at test load) | 30 mg (1 kg) | 30 mg (2 kg) | 30 mg (2 kg) |
| Sensitivity offset (at nominal load) ▲ | 80 mg | 80 mg | 80 mg |
| Sensitivity temperature drift | 0.0003%/°C | 0.0003%/°C | 0.0003%/°C |
| Typical values | | | |
| Repeatability (at 5% load) | 4 mg | 4 mg | 4 mg |
| Linearity deviation | 6 mg | 6 mg | 6 mg |
| Eccentricity deviation (at test load) | 10 mg (1 kg) | 10 mg (2 kg) | 10 mg (2 kg) |
| Sensitivity offset (at nominal load) ▲ | 50 mg | 50 mg | 50 mg |
| Minimum weight (USP, tolerance = 0.10%) ▼ | 8.2 g | 8.2 g | 8.2 g |
| Minimum weight (tolerance = 1%) ▼ | 820 mg | 820 mg | 820 mg |
| Settling time | 1 s | 1 s | 1 s |
| Dimensions and other specifications | | | |
| Balance dimensions (W × D × H) | 194 × 379 × 103 mm | 194 × 379 × 103 mm | 194 × 379 × 103 mm |
| Weighing pan dimensions (W × D) | 170 × 203 mm | 170 × 203 mm | 170 × 203 mm |
| Usable height of draft shield | – | – | – |
| Balance weight | 4.9 kg | 4.9 kg | 5.5 kg |
| Weights for routine testing | | | |
| Weights (OIML class) | 2000 g (F2) / 100 g (F2) | 2000 g (F2) / 200 g (F2) | 5000 g (F2) / 200 g (F2) |
| Weights (ASTM class) | 2000 g (ASTM 1) / 100 g (ASTM 1) | 2000 g (ASTM 4) / 200 g (ASTM 4) | 5000 g (ASTM 4) / 200 g (ASTM 4) |

▲ after adjustment with internal weight

▼ determined at 5% load, k = 2

| | MX6002DR | MX12002 |
|--|----------------------------------|---------------------------------|
| Limit values | | |
| Capacity | 6.2 kg | 12.2 kg |
| Nominal load | 6 kg | 12 kg |
| Readability | 0.1 g | 0.01 g |
| Capacity of fine range | 1.2 kg | – |
| Readability in fine range | 0.01 g | – |
| Repeatability (at 5% load) | 8 mg | 8 mg |
| Linearity deviation | 30 mg | 20 mg |
| Eccentricity deviation (at test load) | 100 mg (2 kg) | 40 mg (5 kg) |
| Sensitivity offset (at nominal load) ▲ | 120 mg | 70 mg |
| Sensitivity temperature drift | 0.0003%/°C | 0.0003%/°C |
| Typical values | | |
| Repeatability (at 5% load) | 4 mg | 4 mg |
| Linearity deviation | 10 mg | 6 mg |
| Eccentricity deviation (at test load) | 32 mg (2 kg) | 12 mg (5 kg) |
| Sensitivity offset (at nominal load) ▲ | 80 mg | 40 mg |
| Minimum weight (USP, tolerance = 0.10%) ▼ | 8.2 g | 8.2 g |
| Minimum weight (tolerance = 1%) ▼ | 820 mg | 820 mg |
| Settling time | 1 s | 1 s |
| Dimensions and other specifications | | |
| Balance dimensions (W × D × H) | 194 × 379 × 103 mm | 194 × 379 × 103 mm |
| Weighing pan dimensions (W × D) | 170 × 203 mm | 170 × 203 mm |
| Usable height of draft shield | – | – |
| Balance weight | 5.5 kg | 5.6 kg |
| Weights for routine testing | | |
| Weights (OIML class) | 5000 g (F2) / 200 g (F2) | 10 kg (F2) / 500 g (F2) |
| Weights (ASTM class) | 5000 g (ASTM 4) / 200 g (ASTM 4) | 10 kg (ASTM 4) / 500 g (ASTM 4) |

▲ after adjustment with internal weight

▼ determined at 5% load, k = 2

| | MX6001 | MX8001 |
|--|----------------------------------|----------------------------------|
| Limit values | | |
| Capacity | 6.2 kg | 8.2 kg |
| Nominal load | 6 kg | 8 kg |
| Readability | 0.1 g | 0.1 g |
| Capacity of fine range | – | – |
| Readability in fine range | – | – |
| Repeatability (at 5% load) | 50 mg | 50 mg |
| Linearity deviation | 60 mg | 100 mg |
| Eccentricity deviation (at test load) | 200 mg (2 kg) | 300 mg (5 kg) |
| Sensitivity offset (at nominal load) ▲ | 240 mg | 400 mg |
| Sensitivity temperature drift | 0.0005%/°C | 0.0005%/°C |
| Typical values | | |
| Repeatability (at 5% load) | 40 mg | 40 mg |
| Linearity deviation | 20 mg | 30 mg |
| Eccentricity deviation (at test load) | 60 mg (2 kg) | 100 mg (5 kg) |
| Sensitivity offset (at nominal load) ▲ | 150 mg | 250 mg |
| Minimum weight (USP, tolerance = 0.10%) ▼ | 82 g | 82 g |
| Minimum weight (tolerance = 1%) ▼ | 8.2 g | 8.2 g |
| Settling time | 0.8 s | 1 s |
| Dimensions and other specifications | | |
| Balance dimensions (W × D × H) | 194 × 379 × 104 mm | 194 × 379 × 104 mm |
| Weighing pan dimensions (W × D) | 172 × 205 mm | 172 × 205 mm |
| Usable height of draft shield | – | – |
| Balance weight | 5.2 kg | 5.2 kg |
| Weights for routine testing | | |
| Weights (OIML class) | 5000 g (F2) / 200 g (F2) | 5000 g (F2) / 200 g (F2) |
| Weights (ASTM class) | 5000 g (ASTM 4) / 200 g (ASTM 4) | 5000 g (ASTM 4) / 200 g (ASTM 4) |

▲ after adjustment with internal weight

▼ determined at 5% load, k = 2

9.4.4 Precision balances, large

| | MX12001L | MX16001L | MX32001L | MX32000L |
|--|---------------------------------|---------------------------------|--------------------------------|--------------------------------|
| Limit values | | | | |
| Capacity | 12.2 kg | 16.2 kg | 32.2 kg | 32.2 kg |
| Nominal load | 12 kg | 16 kg | 30 kg | 30 kg |
| Readability | 0.1 g | 0.1 g | 0.1 g | 1 g |
| Capacity of fine range | – | – | – | – |
| Readability in fine range | – | – | – | – |
| Repeatability (at 5% load) | 80 mg | 80 mg | 80 mg | 600 mg |
| Linearity deviation | 200 mg | 200 mg | 250 mg | 300 mg |
| Eccentricity deviation (at test load) | 300 mg (5 kg) | 300 mg (5 kg) | 300 mg (10 kg) | 1 g (10 kg) |
| Sensitivity offset (at nominal load) ▲ | 600 mg | 600 mg | 800 mg | 1 g |
| Sensitivity temperature drift | 0.0015%/°C | 0.0015%/°C | 0.0015%/°C | 0.0015%/°C |
| Typical values | | | | |
| Repeatability (at 5% load) | 40 mg | 40 mg | 40 mg | 400 mg |
| Linearity deviation | 60 mg | 60 mg | 80 mg | 100 mg |
| Eccentricity deviation (at test load) | 100 mg (5 kg) | 100 mg (5 kg) | 100 mg (10 kg) | 300 mg (10 kg) |
| Sensitivity offset (at nominal load) ▲ | 400 mg | 400 mg | 500 mg | 600 mg |
| Minimum weight (USP, tolerance = 0.10%) ▼ | 82 g | 82 g | 82 g | 820 g |
| Minimum weight (tolerance = 1%) ▼ | 8.2 g | 8.2 g | 8.2 g | 82 g |
| Settling time | 1.5 s | 1.5 s | 1.5 s | 1.2 s |
| Dimensions and other specifications | | | | |
| Balance dimensions (W × D × H) | 354 × 380 × 126 mm | 354 × 380 × 126 mm | 354 × 380 × 126 mm | 354 × 380 × 126 mm |
| Weighing pan dimensions (W × D) | 352 × 246 mm | 352 × 246 mm | 352 × 246 mm | 352 × 246 mm |
| Usable height of draft shield | – | – | – | – |
| Balance weight | 11.7 kg | 11.7 kg | 11.7 kg | 11.7 kg |
| Weights for routine testing | | | | |
| Weights (OIML class) | 10 kg (F2) / 500 g (F2) | 10 kg (F2) / 500 g (F2) | 20 kg (F2) / 1 kg (F2) | 20 kg (F2) / 1 kg (F2) |
| Weights (ASTM class) | 10 kg (ASTM 4) / 500 g (ASTM 4) | 10 kg (ASTM 4) / 500 g (ASTM 4) | 20 kg (ASTM 4) / 1 kg (ASTM 4) | 20 kg (ASTM 4) / 1 kg (ASTM 4) |

▲ after adjustment with internal weight


▼ determined at 5% load, k = 2

9.5 Dimensions

Dimensions in mm.

9.5.1 Analytical balances, readability 0.01 mg


Balance models: MX105, MX105DU, MX205DU



| | |
|----------|------------------------------------|
| ↔ | Outer dimensions [mm] |
| ⇄ | Clear dimensions [mm] |
| 3 | Position of the weighing hook axle |

9.5.2 MX analytical balances, readability 0.1 mg


Balance models: MX104, MX204, MX304



| | |
|--|------------------------------------|
| | Outer dimensions [mm] |
| | Clear dimensions [mm] |
| | Position of the weighing hook axle |

9.5.3 MX precision balances, readability 1 mg, with draft shield


Balance models: MX303, MX603, MX1203



| | |
|--|------------------------------------|
| | Outer dimensions [mm] |
| | Clear dimensions [mm] |
| | Position of the weighing hook axle |

9.5.4 MX precision balances, readability 1 mg, without draft shield


Balance models: MX303N, MX603N, MX1203N



| | |
|----------|------------------------------------|
| ↔ | Outer dimensions [mm] |
| ↔ | Clear dimensions [mm] |
| 3 | Position of the weighing hook axle |

9.5.5 MX precision balances, readability 0.01 g


Balance models: MX2002, MX4002, MX6002, MX6002DR, MX12002



| | |
|----------|------------------------------------|
| ↔ | Outer dimensions [mm] |
| ↔ | Clear dimensions [mm] |
| J | Position of the weighing hook axle |

9.5.6 MX precision balances, readability 0.1 g


Balance models: MX6001, MX8001



| | |
|----------|------------------------------------|
| ↔ | Outer dimensions [mm] |
| ◁ ▷ | Clear dimensions [mm] |
| 3 | Position of the weighing hook axle |

9.5.7 MX large balances, readability 0.1 g / 1 g

Balance models: MX12001L, MX16001L, MX32001L, MX32000L



| | |
|----------|------------------------------------|
| ↔ | Outer dimensions [mm] |
| ↔ | Clear dimensions [mm] |
| J | Position of the weighing hook axle |

10 Accessories and Spare Parts



10.1 Accessories

Accessories are additional components that could help you in your workflow.

Weighing pans

| | | |
|---|---|-----------------|
|  | Weighing pan <ul style="list-style-type: none">• Compatible with 0.01 mg models• ø 80 mm | 30938253 |
|  | Weighing pan <ul style="list-style-type: none">• Compatible with 0.1 mg models• ø 90 mm | 30938254 |
|  | Weighing pan <ul style="list-style-type: none">• Compatible with SmartPan (pro) weighing pan 127 x 127 mm• 128 x 128 mm | 30215433 |
|  | Weighing pan <ul style="list-style-type: none">• Compatible with SmartPan (pro) weighing pan 170 x 203 mm• 172 x 205 mm | 30215056 |

Antistatic kits

| | | |
|---|--|-----------------|
|  | Antistatic kit universal <ul style="list-style-type: none">• Removes electrostatic charges from weighing samples and tare containers• Including: U-electrode large (with installation instructions), high-voltage power supply (with user manual and country-specific power cable) | 11107767 |
|  | High-voltage power supply <ul style="list-style-type: none">• Supplies up to 2 U-electrodes• Including: country-specific power cable, user manual• Compatible with: U-electrode large, U-electrode small | 11107766 |



U-electrode large**11107764**

- Removes electrostatic charges from weighing samples and tare containers
- High-voltage cable with capacitively coupled connector



U-electrode small**11140161**

- Removes electrostatic charges from weighing samples and tare containers
- High-voltage cable with capacitively coupled connector



Ionizer ASK350**30893023**

- Removes small electrostatic charges from weighing samples and tare containers

Density determination

Density kit**30706714**

- Gravimetric density determination of solids



Thermometer, calibrated**11132685**

- Including: holder, calibration certificate
- For usage in density determination

Printers

Printer USB-P25**30702998**

- Printing technology: dot matrix
-



Printer P-52RUE**30237290**

- Printing technology: dot matrix



Printing paper roll, self-adhesive, dot matrix**11600388**

- Set of 3 rolls
- Compatible with: dot matrix printers



Printing paper roll, standard, dot matrix**72456**

- Set of 5 rolls
- Compatible with: dot matrix printers



Ribbon cartridge**65975**

- Including: 2 pcs
- Compatible with: dot matrix printers

Anti-theft devices

Anti-theft cable**11600361**

Hands-free accessories

Foot switch**30312558**

- Hands-free taring, zeroing, printing
-

Barcode readers



Barcode reader 1D Gryphon GD4220

30417466

- Scans barcodes and transmits the decoded information to a connected device
- Interface: USB-A

Cables



Cable USB-A (f) – USB-C (m)

30893021

- Data transfer between instrument and USB-A peripheral
- Length: 0.16 m



USB-C (m) – USB-A (m)

30893022

- Data transfer between instrument and PC
- Length: 1.5 m



Cable USB-A (m) – USB-B (m)

30241476

- Data transfer between instrument and peripheral
- Length: 1 m



Cable RS232 (f) – USB-A (m)

30576241

- Data transfer between balance and peripheral
- Length: 1.7 m



Cable RS232 (m) – USB-A (m)

64088427

- Data transfer between balance and peripheral
- Length: 2 m

Wireless interfaces



Bluetooth adapter ADP-BT-S, single

30086494

- Creates a bluetooth connection between instrument and peripheral

**Bluetooth/Wi-Fi combi adapter LM842****30893006**

- Creates a Bluetooth/Wi-Fi connection between instrument and peripheral

**Bluetooth/Wi-Fi combi adapter LM842, US****30893005**

- Creates a Bluetooth/Wi-Fi connection between instrument and peripheral

Software

EasyDirect Balance

EasyDirect Balance, 10 licenses**30540473**

- Data management software for up to 10 balances
- Collection, analysis, storage and export of weighing data



EasyDirect Balance

EasyDirect Balance, 3 licenses**30539323**

- Data management software for up to 3 balances
- Collection, analysis, storage and export of weighing data

Adjustment weights**Weights**

- For routine testing and calibration of weighing instruments
- Available in different accuracy classes
- With calibration certificate (OIML/ASTM)

► www.mt.com/weights


Dust covers**Dust cover****30893018**



Dust cover

30893019


Draft shields



External draft shield

30706715

- Protects against air currents and dust to maintain measurement accuracy
- Doors: glass; frame: acrylic, aluminium



Draft shield

30938251

- Protects against air currents and dust to maintain measurement accuracy
- Compatible with 1 mg MX balances

Various



EasyHub USB

30468768

- Connects up to 4 peripherals
- Interface to host: USB-B



SmartPrep weighing funnel

30061260

- For weighing powdery substances
- Including: 50 pcs



Protective cover

30106207

- Protects the weighing pan
- Compatible with weighing pan 172 x 205 mm
- 172 x 205 mm



Pipetting container 50 ml

30215436


- Minimizes evaporation in pipette calibration

10.2 Spare parts

Spare parts are parts that are delivered with the original instrument but that can be replaced, if needed, without the help of a service technician.

10.2.1 MX analytical balances, readability 0.01 mg


Balance models: MX105, MX105DU, MX205DU



| | Order no. | Designation | Remarks |
|----|-----------|----------------------------------|--|
| 1 | 30706612 | Door, left | Material: glass; including: door handle |
| 2 | 30706614 | Panel, front | Material: glass |
| 3 | 30706654 | Protective cover | For platform |
| 4 | 30706652 | Protective cover | For terminal |
| 5 | 30104835 | Leveling foot | Including: 2 pcs |
| 6 | 30706724 | Cover, weighing hook | Including: 1 round cover, 1 oval cover; material: silicone |
| 7 | 30706613 | Door, right | Material: glass; including: door handle |
| 8 | 30706618 | Drip tray | – |
| 9 | 30706646 | Draft-protection element | – |
| 10 | 30706631 | Weighing pan \varnothing 80 mm | Including: pan support |
| 11 | 30706611 | Door, top | Material: glass; including: door handle |

10.2.2 MX analytical balances, readability 0.1 mg


Balance models: MX104, MX204, MX304



| | Order no. | Designation | Remarks |
|-----------|-----------|----------------------------------|--|
| 1 | 30706612 | Door, left | Material: glass; including: door handle |
| 2 | 30706614 | Panel, front | Material: glass |
| 3 | 30706655 | Protective cover | For platform |
| 4 | 30706652 | Protective cover | For terminal |
| 5 | 30104835 | Leveling foot | Including: 2 pcs |
| 6 | 30706724 | Cover, weighing hook | Including: 1 round cover, 1 oval cover; material: silicone |
| 7 | 30706613 | Door, right | Material: glass; including: door handle |
| 8 | 30706618 | Drip tray | — |
| 9 | 30706632 | Weighing pan \varnothing 90 mm | Including: pan support |
| 10 | 30706615 | Panel back | Material: blurred glass |
| 11 | 30706611 | Door, top | Material: glass; including: door handle |


10.2.3 MX precision balances, readability 1 mg, with draft shield

Balance models: MX303, MX603



| | Order no. | Designation | Remarks |
|----|-----------|------------------------|--|
| 1 | 30706612 | Door, left | Material: glass; including: door handle |
| 2 | 30706614 | Panel, front | Material: glass |
| 3 | 30706655 | Protective cover | For platform |
| 4 | 30706652 | Protective cover | For terminal |
| 5 | 30104835 | Leveling foot | Including: 2 pcs |
| 6 | 30706724 | Cover, weighing hook | Including: 1 round cover, 1 oval cover; material: silicone |
| 7 | 30706613 | Door, right | Material: glass; including: door handle |
| 8 | 30706618 | Drip tray | – |
| 9 | 30706633 | SmartPan, weighing pan | 128 × 128 mm |
| 10 | 30706615 | Panel back | Material: blurred glass |
| 11 | 30706611 | Door, top | Material: glass; including: door handle |


Balance model: MX1203



| | Order no. | Designation | Remarks |
|-----------|------------------|------------------------|--|
| 1 | 30706612 | Door, left | Material: glass; including: door handle |
| 2 | 30706614 | Panel, front | Material: glass |
| 3 | 30706655 | Protective cover | For platform |
| 4 | 30706652 | Protective cover | For terminal |
| 5 | 30104835 | Leveling foot | Including: 2 pcs |
| 6 | 30706724 | Cover, weighing hook | Including: 1 round cover, 1 oval cover; material: silicone |
| 7 | 30706613 | Door, right | Material: glass; including: door handle |
| 8 | 30706618 | Drip tray | — |
| 9 | 30706634 | SmartPan, weighing pan | 128 × 128 mm |
| 10 | 30706615 | Panel back | Material: blurred glass |
| 11 | 30706611 | Door, top | Material: glass; including: door handle |


10.2.4 MX precision balances, readability 1 mg, without draft shield

Balance models: MX303N, MX603N



| | Order no. | Designation | Remarks |
|---|-----------|------------------------|--|
| 1 | 30706633 | SmartPan, weighing pan | 128 × 128 mm |
| 2 | 30706648 | Drip tray | – |
| 3 | 30706655 | Protective cover | For platform |
| 4 | 30706652 | Protective cover | For terminal |
| 5 | 30104835 | Leveling foot | Including: 2 pcs |
| 6 | 30706724 | Cover, weighing hook | Including: 1 round cover, 1 oval cover; material: silicone |


Balance model: MX1203N



| | Order no. | Designation | Remarks |
|---|-----------|------------------------|--|
| 1 | 30706634 | SmartPan, weighing pan | 128 × 128 mm |
| 2 | 30706648 | Drip tray | – |
| 3 | 30706655 | Protective cover | For platform |
| 4 | 30706652 | Protective cover | For terminal |
| 5 | 30104835 | Leveling foot | Including: 2 pcs |
| 6 | 30706724 | Cover, weighing hook | Including: 1 round cover, 1 oval cover; material: silicone |

10.2.5 MX precision balances, readability 0.01 g


Balance models: MX2002, MX4002, MX6002, MX6002DR, MX12002



| | Order no. | Designation | Remarks |
|---|-----------|------------------------|--|
| 1 | 30706635 | SmartPan, weighing pan | 170 × 203 mm |
| 2 | 30706651 | Cap, pan support | Including: 4 pcs |
| 3 | 30706649 | Drip tray | – |
| 4 | 30706653 | Protective cover | For platform |
| 5 | 30706652 | Protective cover | For terminal |
| 6 | 30104835 | Leveling foot | Including: 2 pcs |
| 7 | 30706724 | Cover, weighing hook | Including: 1 round cover, 1 oval cover; material: silicone |

10.2.6 MX precision balances, readability 0.1 g


Balance models: MX6001, MX8001



| | Order no. | Designation | Remarks |
|----------|-----------|----------------------|--|
| 1 | 30215056 | Weighing pan | 172 × 205 mm |
| 2 | 30706645 | Pan support | – |
| 3 | 30706651 | Cap, pan support | Including: 4 pcs |
| 4 | 30706649 | Drip tray | – |
| 5 | 30706653 | Protective cover | For platform |
| 6 | 30706652 | Protective cover | For terminal |
| 7 | 30104835 | Leveling foot | Including: 2 pcs |
| 8 | 30706724 | Cover, weighing hook | Including: 1 round cover, 1 oval cover; material: silicone |

10.2.7 MX large balances, readability 0.1 g / 1 g

Balance models: MX12001L, MX16001L, MX32001L, MX32000L




| | Order no. | Designation | Remarks |
|---|-----------|------------------|------------------|
| 1 | 30849994 | Weighing pan | 246 × 252 mm |
| 2 | 30849993 | Cap, pan support | Including: 4 pcs |
| 3 | 30706652 | Protective cover | For terminal |
| 4 | 30850018 | Leveling foot | Including: 4 pcs |

10.2.8 AC/DC adapters

10.2.8.1 AC/DC adapter


Compatible with all MX balance models.



| | Order no. | Designation | Remarks |
|-----------|-----------|--------------------|---------------------|
| 1 | 11107909 | AC/DC adapter | Output: 12 V, 2.5 A |
| 2 | 88751 | Power cable AU | — |
| 3 | 30015268 | Power cable BR | — |
| 4 | 87920 | Power cable CH | — |
| 5 | 30047293 | Power cable CN | — |
| 6 | 87452 | Power cable DK | — |
| 7 | 87925 | Power cable EU | — |
| 8 | 89405 | Power cable GB | — |
| 9 | 225297 | Power cable IL | — |
| 10 | 11600569 | Power cable IN | — |
| 11 | 87457 | Power cable IT | — |
| 12 | 11107881 | Power cable JP | — |
| 13 | 11107880 | Power cable TH, PE | — |
| 14 | 88668 | Power cable US | — |
| 15 | 89728 | Power cable ZA | — |

10.2.8.2 AC/DC adapter, universal


Compatible with most MX balance models. Not compatible with the following balance models: MX105, MX105DU, MX205DU, MX104, MX204, MX304



| | Order no. | Designation | Remarks |
|----------|------------------|-------------------------|--|
| 1 | 30850039 | AC/DC adapter universal | Output: 12 V, 1.5 A; including: 6 plugs (EU, UK, US, AU, CN, KR) |


10.2.9 Packaging

10.2.9.1 Balances with draft shield




| | Order no. | Designation | Remarks |
|----------|------------------|--------------------|--|
| 1 | 30706728 | Packaging | Including: export box, inner protection material |
| 2 | 30706731 | Export box | Excluding: inner protection material |

10.2.9.2 Balances without draft shield



| | Order no. | Designation | Remarks |
|---|-----------|-------------|--|
| 1 | 30706729 | Packaging | Including: export box, inner protection material |
| 2 | 30706732 | Export box | Excluding: inner protection material |

10.2.9.3 Balances, large



| | Order no. | Designation | Remarks |
|---|-----------|-------------|--|
| 1 | 30706730 | Packaging | Including: export box, inner protection material |
| 2 | 30706733 | Export box | Excluding: inner protection material |

11 Disposal

In conformance with the European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE), this equipment may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

Please dispose of this equipment in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions, please contact the responsible authority or the distributor from which you purchased this equipment. Should this equipment be passed on to other parties, the content of this directive must also be passed on to the other party.



12 Compliance Information

National approval documents, e.g., the FCC Supplier Declaration of Conformity, are available online and/or included in the packaging.

▶ www.mt.com/ComplianceSearch

Contact METTLER TOLEDO for questions about the country-specific compliance of your instrument.

▶ www.mt.com/contact

United States of America

This equipment has been tested and found to comply with the limits for a **Class A** digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

1. This device may not cause harmful interference.
2. This device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Index

A

| | |
|-----------------------|----------|
| AC/DC adapter | 125, 126 |
| acclimatization | |
| time | 125 |
| adjustment | 52 |
| external | 53 |
| internal | 31, 53 |
| setting | 106 |
| strategy | 52 |
| altitude | 125 |
| application | |
| back-weighing | 50 |
| check weighing | 46 |
| counting | 46 |
| density | 50 |
| differential weighing | 51 |
| dynamic weighing | 47 |
| formulation | 48 |
| totaling | 49 |
| weighing | 42, 45 |
| automation | 44 |

B

| | |
|----------------------------|--------|
| back-weighing | 50 |
| balance | |
| block | 72 |
| reset | 119 |
| unblock | 72 |
| barcode reader | 60 |
| below-the-balance weighing | 35 |
| block | 72 |
| Bluetooth | 56, 83 |
| brightness | 81 |

C

| | |
|------------------------|------------|
| check weighing | 46 |
| combined test weight | 41 |
| command | |
| MT-SICS | 62, 64, 83 |
| compliance information | 163 |
| configuration | |
| report | 45 |
| weighing | 43 |
| convention | 7 |
| counting | 46 |

| | |
|------------|----|
| cover | |
| protection | 18 |

D

| | |
|-----------------------|------------|
| density | 50 |
| device | |
| barcode reader | 60 |
| edit setting | 62 |
| foot switch | 60 |
| keyboard | 61 |
| peripheral | 57, 83 |
| printer | 57 |
| differential weighing | 51 |
| dimension | 136 |
| disposal | 162 |
| door | |
| ErgoDoor handle | 18, 32 |
| handle | 17, 18, 32 |
| lock | 19 |
| draft protection | 17 |
| draft shield | 16 |
| QuickLock | 18 |
| draft shield door | |
| QuickLock | 19 |
| draft shield panel | |
| release button | 19 |
| drip tray | 17 |
| drop to cursor | 63, 66 |
| dynamic weighing | 47 |

E

| | |
|-------------------------|------------|
| EasyDirect Balance | 63, 65, 83 |
| eccentricity test | |
| perform | 55 |
| setting | 111 |
| environmental condition | 25, 125 |
| ErgoDoor | |
| handle | 18, 32 |
| Ethernet | 56, 82 |
| EULA | 31 |
| external | |
| adjustment | 53 |
| device | 57, 83 |

F

| | |
|------------------|-----|
| factory settings | 119 |
|------------------|-----|

| | |
|-------------|----|
| file server | 64 |
| foot switch | 60 |
| formulation | 48 |

H

| | |
|----------|-----|
| history | 73 |
| humidity | 125 |

I

| | |
|--------------------|--------|
| ID | 43 |
| install | |
| put into operation | 30 |
| site | 25 |
| interface | |
| Bluetooth | 56, 83 |
| Ethernet | 56, 82 |
| internal | |
| adjustment | 31, 53 |

K

| | |
|----------|----|
| keyboard | 61 |
|----------|----|

L

| | |
|---------------|----|
| level | |
| balance | 31 |
| indicator | 21 |
| leveling aid | 73 |
| leveling feet | 18 |
| leveling aid | 73 |
| leveling feet | 18 |
| light | |
| StatusLight | 81 |
| location | 25 |
| lock | |
| door | 19 |
| log in | 71 |
| log out | 71 |

M

| | |
|----------|------------|
| material | 126 |
| MT-SICS | 62, 64, 83 |

O

| | |
|------------|----|
| overview | |
| terminal | 16 |
| type label | 19 |

P

| | |
|--------------------|--------|
| pack | |
| balance | 35 |
| password | |
| reset | 72 |
| perform | |
| eccentricity test | 55 |
| repeatability test | 55 |
| sensitivity test | 54 |
| power consumption | |
| balance | 125 |
| power saving | 16, 32 |
| power supply | |
| see AC/CD adapter | 125 |
| print | |
| report | 45 |
| printer | 57 |
| profile | |
| weighing | 40 |
| protective cover | 18 |

Q

| | |
|-------------------|----|
| QuickLock | |
| draft shield | 18 |
| draft shield door | 19 |

R

| | |
|--------------------|-----|
| range | |
| tolerance | 42 |
| release button | |
| draft shield panel | 19 |
| repeatability test | |
| perform | 55 |
| setting | 110 |
| report | |
| configuration | 45 |
| print | 45 |
| reset | |
| balance | 119 |
| password | 72 |
| routine test | 54 |

S

| | |
|--------------------|----|
| safety information | 10 |
| sample ID | 43 |
| sensitivity test | |
| perform | 54 |

| | | | |
|----------------------|------------|----------------|--------|
| setting | 108 | totaling | 49 |
| series | 44 | transport | |
| service | | long distance | 34 |
| drop to cursor | 63, 66 | short distance | 34 |
| EasyDirect Balance | 63, 65, 83 | type label | |
| file server | 64 | overview | 19 |
| MT-SICS | 62, 64, 83 | U | |
| setting | | unblock | 72 |
| adjustment | 106 | unpack | |
| eccentricity test | 111 | balance | 25 |
| repeatability test | 110 | USB | |
| sensitivity test | 108 | see device | 57, 83 |
| test | 108 | W | |
| software | | warm up | |
| version | 7 | time | 125 |
| sound | | warning symbol | 10 |
| terminal | 81 | weighing | |
| standby | 16, 32 | application | 42 |
| statistics | 44 | configuration | 43 |
| StatusLight | 18, 81 | profile | 40 |
| store | | series | 44 |
| balance | 35 | weighing pan | 17 |
| strategy | | weight | |
| adjustment | 52 | target | 42 |
| switch on/off | 32 | test weight | 40 |
| symbol | 7 | Z | |
| warning | 10 | zero | 16 |
| T | | | |
| tare | 16 | | |
| target weight | 42 | | |
| temperature | 125 | | |
| acclimatization time | 125 | | |
| warm-up time | 125 | | |
| terminal | 18 | | |
| brightness | 81 | | |
| overview | 16 | | |
| sound | 81 | | |
| StatusLight | 81 | | |
| test | 54 | | |
| setting | 108 | | |
| test weight | 40, 41 | | |
| threshold | 91 | | |
| time | | | |
| acclimatization | 125 | | |
| warm up | 125 | | |
| tolerance | 42 | | |

To protect your product's future:

METTLER TOLEDO Service assures the quality, measuring accuracy and preservation of value of this product for years to come.

Please request full details about our attractive terms of service.

► www.mt.com/service

www.mt.com/MX-balances

For more information

Mettler-Toledo GmbH

Im Langacher 44
8606 Greifensee, Switzerland
www.mt.com/contact

Subject to technical changes.
© 08/2024 METTLER TOLEDO. All rights reserved.
30491834B en



30491834