
<b>English</b>	User Manual <b>Analytical and Precision Balances MX</b>
<b>Deutsch</b>	Benutzerhandbuch <b>Analysen- und Präzisionswaagen MX</b>
<b>Español</b>	Manual de usuario <b>Balanzas analíticas y de precisión MX</b>
<b>Français</b>	Guide de l'utilisateur <b>Balances d'analyse et de précision fine MX</b>
<b>简体中文</b>	用户手册 <b>分析与精密天平 MX</b>
<b>日本語</b>	ユーザマニュアル <b>分析天びん・上皿天びん MX</b>
<b>Čeština</b>	Návod k použití <b>Analytické a přesné váhy MX</b>
<b>Dansk</b>	Brugervejledning <b>Analyse- og præcisionsvægte MX</b>
<b>Hrvatski</b>	Korisnički priručnik <b>Analitičke i precizne vage MX</b>
<b>Magyar</b>	Felhasználói útmutató <b>Analitikai és precíziós mérlegek MX</b>

**METTLER TOLEDO**

► [www.mt.com/MX-RM](http://www.mt.com/MX-RM)

**en** For full information, always refer to the Reference Manual (RM).

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---

User Manual **Analytical and Precision Balances**

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Magyar

## Table of Contents

<b>1</b>	<b>Introduction</b>	<b>3</b>
1.1	Document purpose .....	3
1.2	Further documents and information .....	3
1.3	Acronyms and abbreviations .....	3
<b>2</b>	<b>Safety Information</b>	<b>4</b>
2.1	Definition of signal words and warning symbols .....	4
2.2	Product-specific safety notes .....	4
<b>3</b>	<b>Design and Function</b>	<b>5</b>
3.1	Overview analytical balances .....	5
3.2	Overview precision balances, small .....	6
3.2.1	Balances with draft shield .....	6
3.2.2	Balances without draft shield .....	7
3.3	Overview precision balances, large .....	8
3.4	Overview terminal .....	8
3.5	Overview interface connections .....	9
3.6	Components description .....	9
3.6.1	Draft shield .....	9
3.6.2	Weighing pan .....	10
3.6.3	Draft-protection element .....	10
3.6.4	Drip tray .....	10
3.6.5	Door handle .....	10
3.6.6	ErgoDoor handle .....	11
3.6.7	Leveling feet .....	11
3.6.8	Terminal .....	11
3.6.9	QuickLock for draft shield .....	11
3.6.10	QuickLock for top door and front panel .....	12
3.6.11	QuickLock for side door .....	12
3.6.12	Release button for back panel .....	12
3.7	User interface .....	13
3.7.1	Main sections at a glance .....	13
3.7.2	Main weighing screen .....	13
<b>4</b>	<b>Installation and Putting into Operation</b>	<b>14</b>
4.1	Selecting the location .....	14
4.2	Unpacking the balance .....	15
4.3	Installation .....	15
4.3.1	Balances with draft shield .....	16
4.3.2	Balances without draft shield .....	18
4.3.2.1	Assembling balances 1 mg .....	18
4.3.2.2	Assembling balances 0.01 g .....	18
4.3.2.3	Assembling balances 0.1 g .....	18
4.3.3	Balances, large .....	19
4.4	Putting into operation .....	19
4.4.1	Connecting the balance .....	19
4.4.2	Switching on the balance .....	19
4.4.3	Leveling the balance .....	20
4.4.4	Performing an internal adjustment .....	20
4.4.5	Entering / Exiting standby mode .....	20
4.4.6	Entering / Exiting power-saving mode .....	20
4.4.7	Switching off the balance .....	21
4.5	Performing a simple weighing .....	21

4.5.1	Opening and closing draft shield doors .....	21
4.5.2	Zeroing the balance .....	22
4.5.3	Taring the balance .....	22
4.5.4	Performing a weighing .....	22
4.6	Transporting, packing, and storing .....	22
4.6.1	Transporting the balance over short distances .....	22
4.6.2	Transporting the balance over long distances .....	23
4.6.3	Packing and storing .....	23
<b>5</b>	<b>Maintenance</b> .....	<b>23</b>
5.1	Maintenance tasks .....	24
5.2	Cleaning .....	24
5.2.1	Disassembling for cleaning .....	24
5.2.1.1	Balances with draft shield .....	24
5.2.1.2	Balances without draft shield .....	27
5.2.1.3	Balances, large .....	27
5.2.2	Cleaning the balance .....	28
5.2.3	Putting into operation after cleaning .....	28
5.3	Service .....	29
<b>6</b>	<b>Technical Data</b> .....	<b>29</b>
6.1	General data .....	29
<b>7</b>	<b>Disposal</b> .....	<b>30</b>
<b>8</b>	<b>Compliance Information</b> .....	<b>30</b>

## 1 Introduction

Thank you for choosing a METTLER TOLEDO balance. The balance combines high performance with ease of use.

### EULA

The software in this product is licensed under the METTLER TOLEDO End User License Agreement (EULA) for Software.

When using this product you agree to the terms of the EULA.

▶ [www.mt.com/EULA](http://www.mt.com/EULA)

### 1.1 Document purpose

This User Manual provides brief instructions about the first steps to take with the instrument. This ensures safe and efficient handling. Personnel must have carefully read and understood this manual before performing any task.

### 1.2 Further documents and information

This document is available in other languages online.

▶ [www.mt.com/MX-UM](http://www.mt.com/MX-UM)

Product page:

▶ [www.mt.com/MX-balances](http://www.mt.com/MX-balances)

Instructions for cleaning a balance, "8 Steps to a Clean Balance":

▶ [www.mt.com/lab-cleaning-guide](http://www.mt.com/lab-cleaning-guide)

Search for software:

▶ [www.mt.com/labweighing-software-download](http://www.mt.com/labweighing-software-download)

Search for documents:

▶ [www.mt.com/library](http://www.mt.com/library)

For further questions, please contact your authorized METTLER TOLEDO dealer or service representative.

▶ [www.mt.com/contact](http://www.mt.com/contact)

### 1.3 Acronyms and abbreviations

Original term	Explanation
AC	Alternating Current
ASTM	American Society for Testing and Materials
DC	Direct Current
EMC	Electromagnetic Compatibility
FCC	Federal Communications Commission
GWP	Good Weighing Practice
ID	Identification
IP	Ingress Protection
LAN	Local Area Network
LED	Light-Emitting Diode
LPS	Limited Power Source
MAC	Media Access Control
MT-SICS	METTLER TOLEDO Standard Interface Command Set

NA	Not Applicable
OIML	Organisation Internationale de Métrologie Légale (International Organization of Legal Metrology)
RM	Reference Manual
SOP	Standard Operating Procedure
TDNR	Type Definition Number
UM	User Manual
USB	Universal Serial Bus
USP	United States Pharmacopeia

## 2 Safety Information

Two documents named "User Manual" and "Reference Manual" are available for this instrument.

- The User Manual is available online in various languages.
- A printed version of the User Manual is delivered with the instrument.
- The Reference Manual is available online. This manual contains a full description of the instrument and its use.
- Keep both documents for future reference.
- Include both documents if you transfer the instrument to other parties.

Only use the instrument according to the User Manual and the Reference Manual. If you do not use the instrument according to these documents or if the instrument is modified, the safety of the instrument may be impaired and Mettler-Toledo GmbH assumes no liability.

### 2.1 Definition of signal words and warning symbols

Safety notes contain important information on safety issues. Ignoring the safety notes may lead to personal injury, damage to the instrument, malfunctions and false results. Safety notes are marked with the following signal words and warning symbols:

#### Signal words

<b>DANGER</b>	A hazardous situation with high risk, resulting in death or severe injury if not avoided.
<b>WARNING</b>	A hazardous situation with medium risk, possibly resulting in death or severe injury if not avoided.
<b>CAUTION</b>	A hazardous situation with low risk, resulting in minor or moderate injury if not avoided.
<b>NOTICE</b>	A hazardous situation with low risk, resulting in damage to the instrument, other material damage, malfunctions and erroneous results, or loss of data.

#### Warning symbols

General hazard

Notice

### 2.2 Product-specific safety notes

#### Intended use

This instrument is designed to be used by trained staff. The instrument is intended for weighing purposes.

Any other type of use and operation beyond the limits of use stated by Mettler-Toledo GmbH without consent from Mettler-Toledo GmbH is considered as not intended.

#### Responsibilities of the instrument owner

The instrument owner is the person holding the legal title to the instrument and who uses the instrument or authorizes any person to use it, or the person who is deemed by law to be the operator of the instrument. The instrument owner is responsible for the safety of all users of the instrument and third parties.

Mettler-Toledo GmbH assumes that the instrument owner trains users to safely use the instrument in their workplace and deal with potential hazards. Mettler-Toledo GmbH assumes that the instrument owner provides the necessary protective gear.

### Safety notes

#### **WARNING**

##### **Death or serious injury due to electric shock**

Contact with parts that carry a live current can lead to death or injury.

- 1 Only use the METTLER TOLEDO power cable and AC/DC adapter designed for your instrument.
- 2 Connect the power cable to a grounded power outlet.
- 3 Keep all electrical cables and connections away from liquids and moisture.
- 4 Check the cables and the power plug for damage and replace them if damaged.

#### **NOTICE**

##### **Damage to the instrument or malfunction due to the use of unsuitable parts**

- Only use parts from METTLER TOLEDO that are intended to be used with your instrument.

A list of spare parts and accessories can be found in the Reference Manual.

## 3 Design and Function


For further information, consult the Reference Manual (RM).


[www.mt.com/MX-RM](http://www.mt.com/MX-RM)

### 3.1 Overview analytical balances

0.01 mg

0.1 mg


1	StatusLight	9	QuickLock, top door/panel
2	Terminal	10	Top door, draft shield
3	Draft-protection element	11	QuickLock, side door
4	Weighing pan	12	Side door, draft shield (right/left)
5	Drip tray	13	Handle, side door
6	Back panel *, draft shield	14	ErgoDoor handle
7	Front panel, draft shield	15	Leveling feet
8	Handle, top door		

\* On balances with a readability of 0.01 mg, the back panel has a backlight.

## 3.2 Overview precision balances, small

### 3.2.1 Balances with draft shield

1 mg


1	StatusLight	8	QuickLock, top door/panel
2	Terminal	9	Top door, draft shield
3	SmartPan weighing pan	10	QuickLock, side door
4	Drip tray	11	Side door, draft shield (right/left)
5	Back panel, draft shield	12	Handle, side door
6	Front panel, draft shield	13	ErgoDoor handle
7	Handle, top door	14	Leveling feet


### 3.2.2 Balances without draft shield


1 mg

0.01 g


0.1 g

1	StatusLight	6	Weighing pan support
2	Terminal	7	Weighing pan
3	Drip tray	8	Safety feet
4	Weighing pan support cap	9	Leveling feet
5	SmartPan weighing pan		


### 3.3 Overview precision balances, large


0.1 g / 1 g


1	StatusLight	4	Weighing pan
2	Terminal	5	Leveling feet
3	Weighing pan support cap		

### 3.4 Overview terminal

	Name	Description
	<b>Standby / Power-saving mode</b>	By tapping  , the balance goes into standby mode. By tapping and holding  , the balance goes into power-saving mode. To switch the balance completely off, it must be unplugged from the power supply.  <b>Note</b> Do not disconnect the balance from the power supply unless the balance is not used for an extended period of time. After switching on the instrument, it must warm up before giving accurate results.
	Tare	Tares the balance.  This function is used when the weighing process involves containers. After taring the balance, the screen shows <b>Net</b> which indicates that all displayed weights are net.
	Zero	Zeroes the balance.  The balance must always be zeroed before starting the weighing process. After zeroing, the balance sets a new zero point.
	Home	To return from any menu level to the main weighing screen.

### 3.5 Overview interface connections


<b>1</b>	Slot for anti-theft cable	<b>4</b>	USB-C port
<b>2</b>	Ethernet port (LAN)	<b>5</b>	Socket for AC/DC adapter
<b>3</b>	USB-A port		

### 3.6 Components description


#### 3.6.1 Draft shield

The draft shield protects the weighing area against environmental impacts like drafts or moisture. The side doors and the top door can be opened manually.


### 3.6.2 Weighing pan

The weighing pan is the load receptor that serves to accommodate the weighing item.


### 3.6.3 Draft-protection element

The draft-protection element protects the weighing pan against draft. This element is only available for balances with a readability of 0.01 mg.


### 3.6.4 Drip tray

The drip tray is positioned below the weighing pan. The primary purpose of the drip tray is to ensure quick cleaning of the balance.


### 3.6.5 Door handle

The door handles are mounted on the draft shield doors. The handles are used to manually open the side doors and the top door of the draft shield.


### 3.6.6 ErgoDoor handle

The ErgoDoor handle is mounted on the platform. The ErgoDoor handle can be engaged with the handle of the side door. This allows to customize opening/closing of the side doors according to your needs.


### 3.6.7 Leveling feet

The balance stands on height-adjustable feet. These feet are used to level the balance.


### 3.6.8 Terminal

The balance terminal has a 7-inch touch-sensitive display. A StatusLight LED strip on the front side of the terminal indicates the current status of the balance. The terminal is protected by a replaceable cover.


### 3.6.9 QuickLock for draft shield

The QuickLock for the draft shield is used to secure the draft shield to the platform.


### 3.6.10 QuickLock for top door and front panel

Depending on the position, the QuickLock is used to lock/unlock the top door and the front panel of the draft shield.


### 3.6.11 QuickLock for side door

The QuickLock is used to lock/unlock the side door of the draft shield.

### 3.6.12 Release button for back panel


The release button is used to lock/unlock the back panel of the draft shield. This feature is only available for the draft shield of balances with a readability of 0.1 mg and 1 mg.

## 3.7 User interface

### 3.7.1 Main sections at a glance

The main weighing screen (1) is the central navigation point where all the menus and settings can be found. The sections **Balance menu** (2), **Main configuration** (3), and the applications section (4) open when tapping the corresponding icon or tab.

#### See also

[Main weighing screen](#) ▶ Page 13

### 3.7.2 Main weighing screen

	Name	Description
1	Weighing results	Shows the results of the current weighing process.
2	Level indicator	Indicates if the balance is leveled (green) or not (red).

	Name	Description
3	Applications	Accesses available applications: <b>Weighing, Adjustments, Tests.</b>
4	Readability and capacity	Shows the readability and the capacity of the balance.
5	Additional information	Shows additional information about the current activity. Example: current weighing value in another unit
6	Information and warnings	Shows current information, warnings, and error messages.
7	Functions area	Shows the active functions according to the settings of the current weighing application.
8	<b>Publish</b> button	Publishes the results according to the settings of the current weighing application.  Depending on the selected weighing application, the button can have different functions.
9	Action bar	Contains actions referring to the current weighing application.
10	<b>Balance menu</b>	Accesses the balance properties.
11	SmartTrac	Used as a weighing aid to define a target weight with upper and lower tolerances.
12	<b>Main configuration</b>	Accesses the configuration options for the current weighing application.

## 4 Installation and Putting into Operation

### 4.1 Selecting the location

A balance is a sensitive precision instrument. The location where it is placed will have a profound effect on the accuracy of the weighing results.

#### Requirements of the location

Place indoors on stable table


Ensure sufficient spacing

Level the instrument

Provide adequate lighting table


Avoid direct sunlight


Avoid vibrations

Avoid strong drafts

Avoid temperature fluctuations

Take into account the environmental conditions. See "Technical Data".

Sufficient spacing for balances: > 15 cm all around the instrument


## 4.2 Unpacking the balance

Check the package, the packaging elements and the delivered components for damages. If any components are damaged, please contact your METTLER TOLEDO service representative.


### Note

Depending on the balance model, the packaging elements and the components may look different.


- 1 Open the box and lift the package out using the lifting strap (1).

- 2 Open the lifting strap (1) and remove the User Manual (2).


- 3 Remove the upper part of the package and carefully unpack the platform (3).  
**NOTICE: Damage to the instrument**  
Do not touch the cone (4) protruding from the platform.

- 4 Remove the protective bag.
- 5 Keep the protective covers installed on the platform and on the terminal.

- 6 Open the lower part of the package (5).
- 7 Carefully lift out the draft shield (6), the draft shield doors (7), and all other items.
- 8 Store all parts of the packaging in a safe place for future use.  
⇒ The balance is ready for assembling.

## 4.3 Installation

### Note

Depending on the balance model, the components may look different.

### 4.3.1 Balances with draft shield


#### **CAUTION**

##### **Injury due to sharp objects or broken glass**


Instrument components, e.g., glass, can break and lead to injuries.

- Always proceed with focus and care.


- 1 Assemble the draft shield: Turn the QuickLock (1, right, left) and slide in the top door (2).

- 2 Attach the front panel (3), then turn the QuickLock (1, right, left) to hold the panel in place.


- 3 For balances without backlight: Attach the back panel (4) with the smooth side facing towards the front.

- 4 Slide in the side door (5) until the QuickLock (6) snaps in (right, left).

➔ The draft shield is assembled.


- 5 Place the draft shield (7) on top of the platform (8).

**Note**

To protect your balance, keep the protective covers installed on the platform (8) and on the terminal (9).

- 6 Secure the draft shield (7) to the platform (8) by turning the QuickLock (10).


- 7 Fully open the side door (11).

- 8 Insert the drip tray (12).

**Note**


If you cannot install the drip tray properly, make sure the QuickLock (10) is correctly locked.

- 9 Only for balances with a readability of 0.01 mg: Place the draft-protection element (13) on top of the drip tray (12).

- 10 Install the weighing pan (14).


➔ The balance is ready for use.

### 4.3.2 Balances without draft shield


#### 4.3.2.1 Assembling balances 1 mg

- 1 Place the drip tray (1) on top of the platform (2).
- Note**  
To protect your balance, keep the protective covers installed on the platform (2) and on the terminal (3).
- 2 Place the weighing pan (4) on top of the drip tray (1).  
➔ The balance is ready for use.


#### 4.3.2.2 Assembling balances 0.01 g

- 1 Place the drip tray (1) on top of the platform (2).
- Note**  
To protect your balance, keep the protective covers installed on the platform (2) and on the terminal (3).
- 2 Attach the support caps (4).
- 3 Place the weighing pan (5) on top of the support caps (4).  
➔ The balance is ready for use.

#### 4.3.2.3 Assembling balances 0.1 g

- 1 Place the drip tray (1) on top of the platform (2).
- Note**  
To protect your balance, keep the protective covers installed on the platform (2) and on the terminal (3).
- 2 Attach the support caps (4).
- 3 Place the weighing pan support (5) on top of the support caps (4).
- 4 Place the weighing pan (6) on top of the weighing pan support (5).  
➔ The balance is ready for use.


### 4.3.3 Balances, large

- 1 Remove the transport locks (1).
- 2 Place the weighing pan (2) on top of the support caps (3).

**Note**

To protect your balance, keep the protective cover installed on the terminal (4).

⇒ The balance is ready for use.

## 4.4 Putting into operation

### 4.4.1 Connecting the balance


**WARNING**

**Death or serious injury due to electric shock**

Contact with parts that carry a live current can lead to death or injury.

- 1 Only use the METTLER TOLEDO power cable and AC/DC adapter designed for your instrument.
- 2 Connect the power cable to a grounded power outlet.
- 3 Keep all electrical cables and connections away from liquids and moisture.
- 4 Check the cables and the power plug for damage and replace them if damaged.

- 1 Install the cables in such a way that they cannot be damaged or interfere with operation.
  - 2 Insert the plug of the AC/DC adapter (1) into the power socket of the instrument.
  - 3 Secure the plug by firmly tightening the knurled nut.
  - 4 Insert the plug of the power cable into a grounded power outlet that is easily accessible.
- ⇒ The balance automatically switches on.

**Note**

Do not connect the instrument to a power outlet controlled by a switch. After switching on the instrument, it must warm up before giving accurate results.

**See also**

🔗 General data ▶ Page 29

### 4.4.2 Switching on the balance

When connected to the power supply, the balance automatically switches on.

**EULA (End User License Agreement)**

When the balance is switched on the first time, the EULA (End User License Agreement) appears on the screen.

- 1 Read the conditions.
- 2 Tap **I accept the terms in the license agreement** and confirm with **✓ OK**.

- ➔ The main weighing screen appears.

### Acclimatization and warm up

Before the balance gives reliable results, it must:

- acclimatize to the room temperature
- warm up by being connected to the power supply

The acclimatization time and warm-up time for balances are available in "General data".

#### Note

When the balance is exiting standby, it is ready immediately.

#### See also

-  General data ▶ Page 29
-  Entering / Exiting standby mode ▶ Page 20
-  Entering / Exiting power-saving mode ▶ Page 20
-  Switching off the balance ▶ Page 21

### 4.4.3 Leveling the balance


Exact horizontal and stable positioning are essential for repeatable and accurate weighing results.

If the balance is out of level, the level indicator on the main screen turns red.

- 1 On the main weighing screen, tap .
  - ➔ The dialog **Leveling aid** opens.
- 2 Turn the leveling feet (1) as instructed on the display until the dot is in the center of the level indicator.

Alternative access to the dialog **Leveling aid**:

 **Navigation:**  **Balance menu** >  **Leveling aid**

### 4.4.4 Performing an internal adjustment

 **Navigation:**  >  **Applications** >  **Adjustments**

- **Adjustments** is set to **Internal**.

- 1 Option 1: On the main weighing screen, tap  **Adjust**.  
 Option 2: Open the applications section, tap  **Adjustments**, select the adjustment, and tap  **Start**.
  - ➔ The adjustment is executed.
  - ➔ The adjustment results appear.
- 2 Tap  **Finish**.
  - ➔ The balance is ready.

### 4.4.5 Entering / Exiting standby mode

- 1 To enter standby mode, short press .
  - ➔ The display is blue. A QR code for more information about the balance is shown.
- 2 To exit standby mode, short press .
  - ➔ The display is turned on.

### 4.4.6 Entering / Exiting power-saving mode

- 1 To enter power-saving mode, long press  (more than two seconds).

- ➔ The display is dark. The balance is in power-saving mode.
- 2 To exit power-saving mode, long press .
- ➔ The balance is switched on.

### Note

We recommend configuring power-saving times. When the balance exits power-saving mode automatically at the defined time, the balance is ready for use immediately.

If the power-saving mode is terminated manually, the balance must warm up before it can be used.

### See also

 General data ▶ Page 29

## 4.4.7 Switching off the balance

To completely switch off the balance, it must be disconnected from the power supply. By pressing , the balance goes only into standby mode or into power-saving mode.

### Note

When the balance has been completely switched off for some time, it must warm up before it can be used.

### See also

 Switching on the balance ▶ Page 19

 Entering / Exiting standby mode ▶ Page 20

 Entering / Exiting power-saving mode ▶ Page 20

## 4.5 Performing a simple weighing

### Note

A balance with draft shield is used to explain the procedure. For balances without a draft shield, skip the instructions steps concerning the draft shield.



### 4.5.1 Opening and closing draft shield doors

- 1 Open the door manually with the door handle (1).
- 2 Alternatively, use the ErgoDoor handle (2) to open the side door.


#### Note

The following instructions describe one use case for weighing in the sample from the right side.

- 3 Connect the ErgoDoor handle (3) to the door handle on the right side (4).

- 4 Move the ErgoDoor handle on the left side (5) to open and close the door on the right side.

#### See also

🔗 ErgoDoor handle ▶ Page 11


### 4.5.2 Zeroing the balance

- 1 Open the draft shield.
- 2 Clear the weighing pan.
- 3 Close the draft shield.
- 4 Press  $\rightarrow 0 \leftarrow$  to zero the balance.
  - ➔ The balance is zeroed.

### 4.5.3 Taring the balance

If a sample vessel is used, the balance must be tared.

- The balance is zeroed.
- 1 Place the sample vessel (1) on the weighing pan (2).
  - 2 Press  $\rightarrow T \leftarrow$  to tare the balance.
    - ➔ The balance is tared. The icon  appears.

### 4.5.4 Performing a weighing

- 1 Open the draft shield.
- 2 Place the weighing object into the sample vessel.
- 3 Close the draft shield.
  - ➔ The result is displayed.
- 4 Optional, if a printer is connected: Tap  to print the weighing result.

## 4.6 Transporting, packing, and storing

### 4.6.1 Transporting the balance over short distances

#### NOTICE

**Flickering backlight due to disconnection of the draft shield (balances with readability of 0.01 mg)**


Hold the balance by the platform. Never lift the balance by the draft shield.

- 1 Disconnect the AC/DC adapter and unplug all interface cables.
- 2 Hold the balance with both hands and carry it in horizontal position to the target location. Consider the requirements of the location.

If you want to put the balance into operation, proceed as follows:

- 1 Connect in reverse order.
- 2 Give the balance sufficient time to warm up.
- 3 Level the balance.
- 4 Perform an internal adjustment.

#### See also

- 🔗 Selecting the location ▶ Page 14
- 🔗 Switching on the balance ▶ Page 19
- 🔗 Leveling the balance ▶ Page 20
- 🔗 Performing an internal adjustment ▶ Page 20

### 4.6.2 Transporting the balance over long distances

METTLER TOLEDO recommends using the original packaging for transportation or shipment of the balance or balance components over long distances. The elements of the original packaging are developed specifically for the balance and its components and ensure maximum protection during transportation.

#### See also

- 🔗 Unpacking the balance ▶ Page 15

### 4.6.3 Packing and storing

#### Packing the balance

Store all parts of packaging in a safe place. The elements of the original packaging are developed specifically for the balance and its components, and ensure maximum protection during transportation and storage.

#### Storing the balance

Only store the balance under the following conditions:

- Indoor and in the original packaging
- According to the environmental conditions, see "Technical Data"

#### **i** Note

When storing for longer than 6 months, the rechargeable battery may become empty (only date and time get lost).

#### See also

- 🔗 Technical Data ▶ Page 29

## 5 Maintenance

To guarantee the functionality of the balance and the accuracy of the weighing results, a number of maintenance actions must be performed by the user.

For further information, consult the Reference Manual (RM).

▶ [www.mt.com/MX-RM](http://www.mt.com/MX-RM)

## 5.1 Maintenance tasks

Maintenance action	Recommended interval	Remarks
Performing an internal adjustment	<ul style="list-style-type: none"> <li>• Daily</li> <li>• After cleaning</li> <li>• After leveling</li> <li>• After changing the location</li> </ul>	see "Performing an internal adjustment"
Performing routine tests (eccentricity test, repeatability test, sensitivity test). METTLER TOLEDO recommends to at least perform a sensitivity test.	<ul style="list-style-type: none"> <li>• After cleaning</li> <li>• After assembling the balance</li> <li>• After a software update</li> <li>• Depending on your internal regulations (SOP)</li> </ul>	see "Tests" in the Reference Manual
Cleaning	<ul style="list-style-type: none"> <li>• After every use</li> <li>• Depending on the degree of pollution</li> <li>• Depending on your internal regulations (SOP)</li> </ul>	see "Cleaning"
Updating the software	<ul style="list-style-type: none"> <li>• Depending on your internal regulations (SOP).</li> <li>• After a new software release.</li> </ul>	see "Software update" in the Reference Manual

### See also

- 🔗 Performing an internal adjustment ▶ Page 20
- 🔗 Cleaning ▶ Page 24

## 5.2 Cleaning

### 5.2.1 Disassembling for cleaning

#### 📘 Note

Depending on the balance model, the components may look different.

#### 📘 Note

In most cases, it is not necessary to remove the protective covers to clean the balance.

#### 5.2.1.1 Balances with draft shield


#### ⚠️ CAUTION

##### Injury due to sharp objects or broken glass

Instrument components, e.g., glass, can break and lead to injuries.

- Always proceed with focus and care.


- 1 Fully open the side door (1).
- 2 Remove the weighing pan (2).
- 3 Only for balances with a readability of 0.01 mg: Remove the draft-protection element (3).
- 4 Remove the drip tray (4).


- 5 Lift up the QuickLock (5) and pull the side door (1) towards the back to remove it (right, left).

**NOTICE: Damage to the instrument**

Hold the side door (1) tightly when removing it.


- 6 Turn the QuickLock (6, right, left), tilt the front panel (7) towards the front and lift it upwards to remove it.


- 7 Pull the top door (8) towards the front to remove it.

**Note**


Optional, if required: Remove the protective covers for cleaning as described below.


- 8 Remove the protective cover (9) from the terminal (10).

- 9 Open the QuickLock (11) and remove the draft shield (12).


- 10 Remove the protective cover (13) from the platform (14).


### 5.2.1.2 Balances without draft shield

- 1 Remove the weighing pan (1).
- 2 Remove the weighing pan support (2) and/or the support caps (3) (if applicable).
- 3 Remove the drip tray (4).


- 4 Optional, if required: Remove the protective cover (5) from the terminal (6).
- 5 Optional, if required: Remove the protective cover (7) from the platform (8).

### 5.2.1.3 Balances, large

- 1 Remove the weighing pan (1).
- 2 Remove the support caps (2).
- 3 Optional, if required: Remove the protective cover (3) from the terminal (4).

## 5.2.2 Cleaning the balance

### NOTICE

#### Damage to the instrument due to inappropriate cleaning methods

If liquid enters the housing, it can damage the instrument. The surface of the instrument can be damaged by certain cleaning agents, solvents, or abrasives.

- 1 Do not spray or pour liquid on the instrument.
- 2 Only use the cleaning agents specified in the Reference Manual (RM) of the instrument or the guide "8 Steps to a Clean Balance".
- 3 Only use a lightly moistened, lint-free cloth or a tissue to clean the instrument.
- 4 Wipe off any spills immediately.

For further information on cleaning a balance, consult "8 Steps to a Clean Balance".

► [www.mt.com/lab-cleaning-guide](http://www.mt.com/lab-cleaning-guide)

#### Cleaning around the balance

- Remove any dirt or dust around the balance and avoid further contaminations.

#### Cleaning the terminal

- Clean the terminal with a damp cloth or a tissue and a mild cleaning agent.

#### Cleaning the removable parts

- Clean the removed part with a damp cloth or a tissue and a mild cleaning agent or clean in a dishwasher up to 80 °C.

#### Cleaning the balance

- 1 Disconnect the balance from the AC/DC adapter.
- 2 Use a lint-free cloth moistened with a mild cleaning agent to clean the surface of the balance.
- 3 Remove powder or dust with a disposable tissue first.
- 4 Remove sticky substances with a damp lint-free cloth and a mild solvent, e.g., isopropanol or ethanol 70%.

## 5.2.3 Putting into operation after cleaning

- 1 Reassemble the balance.
- 2 Check that the draft shield doors (top, sides) open and close normally (if applicable).
- 3 Reconnect the balance to the AC/DC adapter.
- 4 Check the level status, level the balance if necessary.
- 5 Respect the warm-up time specified in the "Technical Data".
- 6 Perform an internal adjustment.
- 7 Perform a routine test according to the internal regulations of your company. METTLER TOLEDO recommends performing a sensitivity test after cleaning the balance.
- 8 Press **→0←** to zero the balance.
  - ➔ The balance is ready for use.

#### See also

- 🔗 Leveling the balance ► Page 20
- 🔗 Technical Data ► Page 29
- 🔗 Performing an internal adjustment ► Page 20

## 5.3 Service

Regular servicing by an authorized service technician ensures reliability for years to come. Contact your METTLER TOLEDO representative for details about the available service options.

## 6 Technical Data

### 6.1 General data

#### Power supply for balances with a readability of 0.01 mg and 0.1 mg

AC/DC adapter:	Input: 100 – 240 V AC $\pm$ 10%, 50 – 60 Hz, 0.8 A, 61 – 80 VA Output: 12 V DC, 2.5 A, LPS
Cable for AC/DC adapter:	3-core, with country-specific plug
Balance power consumption:	12 V DC, 1.0 A
Polarity:	

#### Power supply for balances with a readability of 1 mg or higher

AC/DC adapter:	Input: 100 – 240 V AC $\pm$ 10%, 50 – 60 Hz, 0.5 A Output: 12 V DC, 1.5 A, LPS
Balance power consumption:	12 V DC, 1.0 A
Polarity:	

#### Protection and standards

Overvoltage category:	II
Degree of pollution:	2
Ingress protection code:	IP41 (balances with readability of 0.1 mg or 1 mg) IP54 (balances with readability of 0.01 g or higher)

#### Note

Stated IP is only achieved when the balance is ready for operation. The protective covers must be installed, and the caps must cover the interface connections.

Standards for safety and EMC:	See Declaration of Conformity
Range of application:	Use only indoors in dry locations

#### Environmental conditions

The limit values apply when the balance is used under the following environmental conditions:

Height above mean sea level:	Up to 5000 m
Ambient temperature:	+10 – +30 °C
Temperature change, max.:	5 °C/h
Relative humidity:	30 – 70%, non-condensing
Acclimatization time:	Recommendation: Up to <b>4 hours</b> for precision balances, or up to <b>8 hours</b> for analytical balances. These values apply after placing the balance in the same location where it will be put into operation.

#### Note

The acclimatization time depends on the readability of the balance, and on the environmental conditions.

Warm-up time:	At least <b>30 minutes</b> for precision balances, <b>60 minutes</b> for analytical balances, or <b>120 minutes</b> for balances with a readability of 0.01 mg. These values apply after connecting the balance to the power supply, or after exiting power-saving mode. When switched on from standby, the balance is ready for operation immediately.
---------------	---

The balance can be used under the following environmental conditions. However, the weighing performances of the balance may be outside the limit values:

Ambient temperature: +5 °C – +40 °C

Relative humidity: 20% to max. 80% at 31 °C, decreasing linearly to 50% at 40 °C, non-condensing

The balance can be disconnected and stored in its packaging under the following conditions:

Ambient temperature: -25 – +70 °C

Relative humidity: 10 – 90%, non-condensing

## 7 Disposal

In conformance with the European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions, please contact the responsible authority or the distributor from which you purchased this device. Should this device be passed on to other parties, the content of this regulation must also be related.

## 8 Compliance Information

National approval documents, e.g., the FCC Supplier Declaration of Conformity, are available online and/or included in the packaging.

▶ [www.mt.com/ComplianceSearch](http://www.mt.com/ComplianceSearch)

For further information, consult the Reference Manual (RM).

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► [www.mt.com/contact](http://www.mt.com/contact)

## 1.3 Akronymmer og forkortelser

Originalt udtryk	Oversat udtryk	Forklaring
AC		Alternating Current
ASTM		American Society for Testing and Materials
DC		Direct Current
EMC		Electromagnetic Compatibility
FCC		Federal Communications Commission
GWP		Good Weighing Practice
ID		Identification
IP		Ingress Protection
LAN		Local Area Network
LED		Light-Emitting Diode
LPS		Limited Power Source
MAC		Media Access Control

MT-SICS	METTLER TOLEDO Standard Interface Command Set
NA	Not Applicable
OIML	Organisation Internationale de Métrologie Légale (International Organization of Legal Metrology)
RM	Reference Manual
SOP	Standard Operating Procedure
TDNR	Type Definition Number
UM	User Manual (Brugervejledning)
USB	Universal Serial Bus
USP	United States Pharmacopeia

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© 07/2024 METTLER TOLEDO. All rights reserved.