

Determination of hydrazine in water samples

Hydrazine is determined in aqueous samples by redox titration using potassium iodate as a titrant.

Sample	Hydrazine sulphate, $(\text{NH}_2)_2\text{-H}_2\text{SO}_4$ 0.04-0.09 g	Preparation and Procedures 1.) Weigh 0.04-0.09 g of hydrazine sulphate in a glass titration beaker (polypropylene beakers will be irreversibly stained). 2.) Add 50 mL deionized water. 3.) Add 2 mL concentrated hydrochloric acid (1.16 sg) in order to initiate the desired reaction.
Substance	$(\text{NH}_2)_2\text{-H}_2\text{SO}_4$, $M = 130.12$ g/mol	
Chemicals	50 mL deion. water, 2 mL conc. hydrochloric acid, HCl (1.16 sg)	
Titrant	Potassium iodate, KIO_3 $c(\text{KIO}_3) = 0.1$ mol/L	
Standard	$\text{Na}_2\text{S}_2\text{O}_3$ / KI	
Instruments	DL50Graphix, DL5x, DL7x Analytical balance, printer	
Accessories	Glass titration beaker, ME-101446 Rondo 60 Sample changer	
Indication	DM140-SC	
Chemistry	$\text{N}_2\text{H}_4 + 2\text{H}^+ + \text{IO}_3^- + \text{Cl}^- =$ $\text{ICl} + \text{N}_2 + 3\text{H}_2\text{O}$	
Calculation	$R1 = \text{VEQ} ; \text{mL}$ $R2 = Q1 * C2 / m ; \%$ $C2 = M / (10 * z)$	
Waste disposal	Neutralize the aqueous acidic phase before final disposal as anorganic salt solution.	
Author	Patrick Hogan , MT-UK	Remarks 1) The hydrazine level is determined to test the water quality for steam turbine. Hydrazine is used as a scavenger buffer system to ensure greater longevity of those parts of the turbine in contact with the constant steam output. 2) In this application, hydrazine sulphate was directly dissolved in water. 3) According to Vogel's "Textbook of quantitative inorganic analysis" (Ed. 5, 1991), hydrazine can be titrated against potassium iodate releasing iodine monochloride as one of the side-products. The concurrent change in redox states means that a straightforward redox reaction is taking place according to the equation under "Chemistry".

Results

METTLER DL70 Titrator

METTLER TOLEDO UK Ltd
Applications Laboratory

pgh1 Hydrazine Titration measured 26-Mar-1993 11:08
26-Mar-1993 9:59 Titrator J70820
SW Version 2.0 User Patrick Hogan

RESULTS

No	Identification	Weight	Results
1/1	(NH2)2S' Aldric	0.0909 g	5.689 mL Equival' Vol 81.423 % Hydrazine
1/2	(NH2)2S' Aldric	0.0909 g	2.936 mL Equival' Vol 79.920 % Hydrazine
1/3	(NH2)2S' Aldric	0.0909 g	3.470 mL Equival' Vol 80.604 % Hydrazine
1/4	(NH2)2S' Aldric	0.0909 g	3.626 mL Equival' Vol 79.947 % Hydrazine
1/5	(NH2)2S' Aldric	0.0909 g	3.554 mL Equival' Vol 81.693 % Hydrazine


STATISTICS

Number results R2 n = 5
Mean value x = 80.718 %
Standard deviation s = 0.8205 %
Rel. standard deviation srel = 1.106 %

Table of measured values

	Volume mL	Increment mL	Signal mV	Change mV	1st deriv. mV/mL	Time mins
E1	0.0000		439.9			0:00
	0.5720	0.5720	448.1	8.2	14.4	0:34
E2	0.8580	0.2860	454.1	6.0	20.9	1:08
	1.2000	0.3420	453.9	-0.4	-3.1	1:38
	1.2000	0.2000	457.8	4.2	20.9	2:08
	1.4000	0.2000	460.2	2.3	11.6	2:37
	1.6000	0.2000	463.3	3.5	17.0	3:05
	1.8000	0.2000	466.0	2.4	11.8	3:34
	2.0000	0.2000	468.0	2.8	14.0	4:02
	2.2000	0.2000	472.8	4.8	19.9	4:28
	2.4000	0.2000	476.8	4.0	20.0	4:54
	2.6000	0.2000	481.1	4.3	21.6	5:20
	2.8000	0.2000	481.5	0.3	19.0	5:48
	3.0000	0.2000	492.7	7.8	39.1	6:14
	3.1340	0.1340	499.3	6.8	50.3	6:36
	3.2640	0.1300	509.3	9.9	75.9	6:56
	3.3420	0.0780	516.9	9.5	122.4	7:11
	3.3900	0.0480	525.9	7.0	146.1	7:20
	3.4360	0.0460	536.1	10.2	212.3	7:29
	3.4660	0.0280	541.2	5.1	183.5	7:34
3.5180	0.0520	576.2	34.9	672.0	7:42	
EGF1	3.5380	0.0200	575.4	-0.8	-37.8	7:45
	3.5580	0.0200	639.8	64.4	3221.3	7:52
	3.5780	0.0200	688.3	48.5	2424.8	8:00
	3.5980	0.0200	704.1	15.8	788.9	8:09

Titration curve



Method

```

Method      pgh1      Hydrazine titration
Version    26-Mar-1993  9:59

Title
Method ID   . . . . . pgh1
Title      . . . . . Hydrazine titration
Date/time  . . . . . 26-Mar-1993  9:59

Sample
Number samples . . . . . 1
Titration stand . . . . . Stand 1
Entry type . . . . . Wight m
  Lower limit [mL] . . . . . 0
  Upper limit [mL] . . . . . 2.0
ID 1 . . . . . (NH2)2S
Molar mass M . . . . . 130.12
Equivalent number z . . . . . 1

Stir
Speed [%] . . . . . 50
Time [s] . . . . . 120

Titration
Titrant . . . . . KIO3
Concentration [mol/L] . . . . . 0.1
Sensor . . . . . DM140-SC
Unit of meas. . . . . mV
Titration mode . . . . . EQP
  Predispensing . . . . . mL
  Volume [mL] . . . . . 1.0
  Titrant addition . . . . . DYN
  dE(set) [mV] . . . . . 8.0
  Limits dV . . . . . Absolute
    dV(min) [mL] . . . . . 0.02
    dV(max) [mL] . . . . . 0.2
  Measure mode . . . . . EQU
    dE [mV] . . . . . 0.5
    dt [s] . . . . . 1.0
    t(min) [s] . . . . . 3.0
    t(max) [s] . . . . . 30.0
  Threshold . . . . . 500.0
  Maximum volume [mL] . . . . . 10.0
  Termination after n EQPs . . . . . Yes
    n = . . . . . 1
  Evaluation procedure . . . . . Standard

Calculation
Result name . . . . . Equival' Vol
Formula . . . . . R1=VEQ1
Constant . . . . .
Result unit . . . . . mL
Decimal places . . . . . 3

Calculation
Result name . . . . . Hydrazine
Formula . . . . . R2=Q1*C2/m
Constant . . . . . C2=M/(10*z)
Result unit . . . . . %
Decimal places . . . . . 3

Record
Output unit . . . . . Printer
All results . . . . . Yes
Table of values . . . . . Yes
E - V curve . . . . . Yes
dE/dV - V curve . . . . . Yes

Statistics
Ri (i=index) . . . . . R2
Standard deviation s . . . . . Yes
Rel. standard deviation srel . . . . . Yes

Rinse
Auxiliary reagent . . . . . H2O
Volume [mL] . . . . . 30.0

Conditioning
Interval . . . . . 1
Time [s] . . . . . 120
Rinse
  Auxiliary reagent . . . . . H2O
  Volume [mL] . . . . . 15.0

Record
Output unit . . . . . Printer
All results . . . . . Yes

```