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

Reducing sugar determination in beverages according to Rebelein	
Reducing sugars are determined in beverages by titration of released iodine with sodium thiosulfate.	
Sample	- Wines and juices
Preparation and procedures	<ul style="list-style-type: none"> - 5 mL of each Fehling A and B solutions - 40 mL deionized water - Boil it for exactly 2 minutes then cool it to the room temperature - 10 mL of 10% H₂SO₄ and 10 mL of 10% KI - The solution is titrated with 0.1 mol/L Na₂S₂O₃.
Compound	<ul style="list-style-type: none"> - Sugar (Glucose), C₆H₁₂O₆ - M= 180.16 g/mol z = 1
Chemicals	<ul style="list-style-type: none"> - Deionized water - Fehling solution A, Fehling solution B - 10% H₂SO₄ - 10% KI
Titrant	<ul style="list-style-type: none"> - Sodium thiosulfate, Na₂S₂O₃ - c(Na₂S₂O₃) = 0.1 mol/L
Standard	- Glucose
Instruments, Accessories	<ul style="list-style-type: none"> - T50/T70/T90 Excellence Titrator (This method can also be used on G20). - DV1020 20 mL glass burette, - 100 mL glass beakers(ME-101446), XP205 Balance, LabX[®] pro titration software.
Indication	- DM140-SC platinum ring electrode.
Results	- Glucose : 2.3589 ± 0.0238 g/L srel = 1.009 % n = 3
Waste	- Precipitate copper by adding sodium hydroxide. The filtrate has to be classified as special waste.
Comments	<ul style="list-style-type: none"> - Use safety goggles, wear gloves and a lab coat. Note: It is recommended while working in a lab. - It may be necessary to adapt the method parameters to your specific sample. - Reducing sugars are oxidized with excess alkaline copper sulphate solution. - Standard solutions are prepared. - Clean the electrode and the stirrer after each titration with deionized water.
Literature	- Mettler-Toledo Titration Applications M340 and M341
Author	Craig Gordon Compiled by: Sohel Ansari, IMSG AnaChem, April 2012 - v1.0 Revised: V. Gärtner – MSG AnaChem, August 2012

Titration method

Sample Determination			002 Sample		
001 Title	Type	General titration	Number of IDs	1	
	Compatible with	T50/T70/T90	ID 1	--	
	ID	ReducingSugar	Entry type	Fixed Volume	
	Title	Reducing Sugar	Volume	10 mL	
002 Sample			003 Titration stand (Rondolino TTL)		
	Number of IDs	1	Type	Manual stand	
	ID 1	--	Titration stand	Manual stand 1	
	Entry type	Volume	004 Stir		
	Lower limit	2 mL	Speed	40%	
	Upper limit	20 mL	Duration	10 s	
003 Titration stand (Manual stand)			005 Titration (EQP) [1]		
Type	Manual stand		Titration		
Titration stand	Manual stand 1		Titration	Na2S2O3	
004 Stir			Concentration	0.1 mol/L	
Speed	40%		Sensor		
Duration	10 s		Type	mV	
005 Titration (EQP) [1]			Sensor	DM140-SC	
Titration			Unit	mV	
Titration	Na2S2O3				
Concentration	0.1 mol/L		Stir		
Sensor			Speed	40%	
Type	mV		Predispense		
Sensor	DM140-SC		Mode	Volume	
Unit	mV		Volume	10.000 mL	
			Wait time	0 s	
Stir			Control		
Speed	40%		Control	User	
Predispense			Titration addition	Dynamic	
Mode	Volume		dE(set value)	8.0 mV	
Volume	1.000 mL		dV(min)	0.02 mL	
Wait time	0 s		dV(max)	0.5 mL	
Control			Meas. val. acquisition	Equilibrium controlled	
Control	User		dE	0.5 mV	
Titration addition	Dynamic		dt	1 s	
dE(set value)	8.0 mV		t(min)	3 s	
dV(min)	0.02 mL		t(max)	30 s	
dV(max)	0.5 mL		Evaluation and recognition		
Meas. val. acquisition	Equilibrium controlled		Procedure	Standard	
dE	0.5 mV		Threshold	100 mV/mL	
dt	1 s		Tendency	Negative	
t(min)	3 s		Ranges	0	
t(max)	30 s		Add. EQP criteria	No	
Evaluation and recognition			Termination		
Procedure	Standard		At Vmax	20.0 mL	
Threshold	100 mV/mL		At potential	No	
Tendency	Negative		At slope	No	
Ranges	0		After number of		
Add. EQP criteria	No		recognized EQPs	Yes	
Termination			Number of EQPs	1	
At Vmax	20.0 mL		Combined termination		
At potential	No		criteria	No	
At slope	No		Accompanying stating		
After number of			Accompanying stating	No	
recognized EQPs	Yes		Condition		
Number of EQPs	1		Condition	No	
Combined termination			006 Calculation R1		
criteria	No		Result	Consumption	
Accompanying stating			Result unit	mmol	
Accompanying stating	No		Formula	R1=Q	
Condition			Constant C=	1	
Condition	No		M	M[None]	
			z	z[None]	
			Decimal places	4	
006 Calculation R1			007 Record		
Result	Content		Summary	Results	
Result unit	g/L		Results	Yes	
Formula	$R1=(B[\text{Back EQP}]-Q)*C/m$				
Constant C=	M/z				
M	M[Glucose]				
z	z[Glucose]		008 End of sample		
Decimal places	3		009 Blank		
			Name	Back EQP	
			Results	Mean[R1]	
007 Record			010 Calculation R2		
Summary	Results		Result	Mean Back value	
Results	Yes		Result unit	mmol	
			Formula	R1=Mean[R1]	
			Constant C=	1	
			M	M[None]	
			z	z[None]	
			Decimal places	4	
008 End of sample			007 Record		
			Summary	Results	
			Results	Yes	
Back Value Determination					
001 Title	Type	General titration			
	Compatible with	T50/T70/T90			
	ID	BackValue			
	Title	Back value			

Titration curve

Formol Number, Acidity and True Brix Value of Orange Juice	
Method for the determination of the formol number, acidity and true Brix value of orange juice according to NEN 2843 and NEN 2861 (Dutch standards)	
Sample	- Concentrated orange juicy, 1 - 4 g
Preparation and procedures	<ul style="list-style-type: none"> - The sample is weighed in the titration beaker and mounted at the titration head. If drinkable orange juice is used, then weigh in approximately 20 g. The sample is diluted using the second pump of the sample changer. - Firstly the acidity content is determined by performing an EP titration to pH 8.2 with NaOH. this is given as citric acid content. - Formol is dispensed by the second burette. After 1 minute reaction time, an EP titration to pH 8.2 with NaOH is started. - The formol number is defined in calculation R4 as mmol NaOH/10g sample, normalized to 65 True Brix: $FN = (mmol\ NaOH * 10 * 10 * 65) / (m * True\ Brix)$ - A blank is run and stored as auxiliary value H2. - Prior titration Brix value of the sample is measured with a refractometer. The value is stored as auxiliary value as H1 and used in calculation R3 (correction for true brix) - The molar mass M is set to 1.0 and is not used.
Compound	- Organic acids, free amino acids
Chemicals	- 50 mL deionized water
Titrant	<ul style="list-style-type: none"> - Sodium Hydroxide, NaOH - $c(NaOH) = 0.5\ mol/L$
Standard	- THAM (TRIS) Tris(hydroxymethyl)aminomethane
Instruments, Accessories	<ul style="list-style-type: none"> - Titration Excellence T50/T70/T90 - Rondo 60 sample changer - 1 additional diaphragma pump (DOSE) - XS205 Balance - 10 mL DV 1010 Burette - 1 additional burette drive DV90 - Titration beaker ME - 101974 - Printer
Indication	- DG111-SC combined pH electrode.
Results	- True Brix : 57.77 ± 0.007 % srel = 0.013 % n = 2
Waste	- Aqueous solution
Comments	<ul style="list-style-type: none"> - Use safety goggles, wear gloves and a lab coat. - Since there are different juice producers, it may be necessary to slightly adapt the method to your specific sample. - The formol value is generally given as the titrant consumption of 0.1 mol/L NaOH needed for 10 g or sample. - True Brix is the corrected measured refractive index for the citric acid concentration.
Literature	- Mettler-Toledo Titration Application M240
Author	Ko Schaap, MT-NL Compiled by : Robin Isyas MSG AnaChem, June 2012 Revised: V. Gärtner – MSG AnaChem, August 2012

Titration method

001 Title
 Type General titration
 Compatible with T50/T70/T90
 ID True Brix
 Title Formol No
 Author Administrator

002 Sample
 Number of IDs 1
 ID 1 --
 Entry type Weight
 Lower limit 0.0 g
 Upper limit 20.0 g
 Density 1.0 g/mL
 Correction factor 1.0
 Temperature 25.0 C

003 Titration stand (Manual TTL)
 Type Rondo/Tower A
 Titration stand Rondo60/1A

004 Pump
 Auxiliary reagent Water
 Volume 20 mL
 Condition No

005 Stir
 Speed [%] 70
 Duration [s] 15
 Condition No

006 Titration (EQP) [1]
Titration
 Titration NaOH
 Concentration 0.5 mol/L
Sensor
 Type pH
 Sensor DG111-SC
 Unit pH
Temperature
 Temperature measurement No
 Stir
 Speed [%] 30
Predispense
 Mode None
 Wait time[s] 5
Control
 End point type Absolute
 Tendency None
 End point value 8.1 pH
 Control band 2.0 pH
 Dosing rate (max) 10 mL/min
 Dosing rate (min) 10 uL/min

Termination
 At EP 10.0
 Termination delay 0 s
 At Vmax 20 mL
 Max.time = s
 Accompanying stating
 Accompanying No
 Condition No

007 Dispense (normal) [1]
 Titration Formol
 Concentration 37
 Volume 10 mL
 Dosing rate 60.0 mL/min
 Condition No

008 Stir
 Speed 70%
 Duration 60s
 Condition No

009 Titration (EP) [2]
 Titration NaOH
 Concentration 0.5 mol/L
Sensor
 Type pH
 Sensor DG111-SC
 Unit pH
Temperature acquisition
 Temperature acquisition No
 Stir
 Speed 30%
Predispense
 Mode None
 Wait time 0 s
Control
 End point type Absolute
 Tendency None
 End point value 8.2 pH
 Control band 2.0 pH
 Dosing rate (max) 10 mL/min
 Dosing rate (min) 10 uL/min
Termination
 At EP Yes
 Termination delay 0 s
 At Vmax 20 mL
 Max.time = s
 Accompanying stating
 Accompanying stating No
 Condition No

010 Calculation R1
 Result Citric acid
 Result unit %
 Formula $R1=Q \cdot C/m$
 Constant C- 192/m*3
 M [None]
 Z z[None]
 Decimal places 1
 Result limits No
 Record statistics Yes
 Extra statistical function No
 Send to buffer No
 Condition No

011 Calculation R2
 Result True Brix corr
 Result unit %
 Formula $R2=R1 \cdot R1 \cdot (-0.00035) + C$
 Constant C- (R1*0.191)+0.024
 M [None]
 Z z[None]
 Decimal places 2
 Result limit No
 Record statistics Yes
 Extra statistical function No
 Send to buffer No
 Condition No

012 Calculation R3
 Result True Brix
 Result unit %
 Formula $R3=H[\text{Brix value}] + R2$
 Constant C- [None]
 M [None]
 Z z[None]
 Decimal places 2
 Result limit No
 Record statistics Yes
 Extra statistical function No
 Send to buffer No
 Condition No

013 Calculation R4
 Result mmol NaOH/10g
 Result unit mmol NaOH/10g
 Formula $R4=6500 \cdot C / (m \cdot R3)$
 Constant C- Q[2]-B[Blank run]+QEX[1]
 M [None]
 Z z[None]
 Decimal places 2
 Result limit No
 Record statistics Yes
 Extra statistical function No
 Send to buffer No
 Condition No

014 Record
 Summary No
 Results Per sample
 Raw results Per sample
 Table of meas.values Lat titration function
 Sample data No
 Resource data No
 E - V Last titration function
 dE/dV - V No
 log dE/dV - V No
 d2E/dV2 - V No
 BETA - V No
 E - t No
 V - t No
 dV/dt - t No
 T - t No
 E - V & dE/dV - V No
 V - t & dV/dt - t No
 Method No
 Series data No

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- ▶ www.mt.com/gwp 用於稱量
- ▶ www.mt.com/gtp 用於滴定
- ▶ www.mt.com/gpp 用於移液
- ▶ www.mt.com/gdrp 用於測定密度和折射率

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© 09/2012 Mettler-Toledo AG
Mettler-Toledo AG