## Content Determination of Benzyl Alcohol

Method for the alkalimetric content determination of Benzyl Alcohol

<b>Sample</b>	Benzyl Alcohol C <sub>7</sub> H <sub>8</sub> O	<b>Preparation and Procedures</b> <p>Weigh a sample of 0.89 to 0.91 g benzyl alcohol, accurately to 0.1 mg, into a 150 mL round bottomed flask. Add 15.0 mL of a freshly made solution of 1 part acetic anhydride to 7 parts pyridine. An equal volume of this solution is necessary for a blank determination. Both the sample and the blank must boil for 30 minutes by refluxing in a waterbath. Finally transfer into 100 mL titration beakers, rinsing twice with 25 mL deionised water.</p>
<b>Substance</b>	Benzyl Alcohol M = 108.1 g/mol	
<b>Chemicals</b>	Acetic anhydride R Pyridine R Deionised water	
<b>Titrant</b>	c(NaOH) = 1 mol/L	
<b>Standard</b>	Potassium hydrogen phthalate (M521)	
<b>Instruments</b>	DL70ES/DL77 and DL58/55+/53+	
<b>Accessories</b>	Titration beaker ME-101974, SP250 Peristaltic pump, Rondo60 Sample changer	
<b>Indication</b>	DG 111-SC	
<b>Chemistry</b>		
<b>Calculation</b>	$R(\%) = (x - Q[2]) \cdot M \cdot /10 \cdot m$ <p>Q = Titrant consumption in mmol  x = Mean value of the blank determination in mmol  M = Molar mass Benzyl Alcohol  m = Sample size in g</p>	
<b>Waste disposal</b>	Halogen-free organic solvents waste.	
<b>Author</b>	L. Pludra	<b>Remarks</b> <p>The titrant consumption for the blank determination is greater than for the benzyl alcohol determination, since the pyridinium acetate content is higher.</p>

## Results

METTLER DL77      Titrator V3.1      Dr. Kade Pharm. Fabrik GmbH  
 Tante Elfriede      Quality Control Laboratory


Method: G057 Cont. Benzyl alcohol      Date: 09-12-1999 13:11  
 User: Salomon      Measured: 10-12-1999 10:59

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