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

O₂ Analyzer InTap portable

Content

1	Safety and Intended Use	46
2	Technical Data and Further Information	47
3	Design and Connections	48
4	Display and Operating Elements	50
5	Menu Structure	52
6	Commissioning	53
7	Calibration	54
8	Maintenance and Repair	56
9	Environmental protection	57
10	Certificates and approvals	58

1 Safety and Intended Use

The O₂ Analyzer InTap™ portable must only be connected, commissioned, and maintained by qualified specialists in full compliance with the instructions in this User Manual, the applicable norms and legal regulations.

The specialist must have read and understood this User Manual and must follow the instructions it contains.

The Reference Manual provides further information on the device and can be downloaded via Internet “www.mt.com/InTap”.

The O₂ Analyzer InTap portable should be operated only by personnel familiar with the device and who are qualified for such work.

Intended Use

The O₂ Analyzer InTap portable is a portable measurement device for dissolved oxygen measurement and for calibration and data logging of optical dissolved oxygen ISM™ sensors (InPro™ 6970 i) installed in the process. The InTap and the sensors communicate via Bluetooth. For this, the separate Bluetooth interface T100 (order number 30 432 819) has to be mounted onto each sensor.

The InTap is designed for use in the beverage industries, in non-hazardous areas. The InTap is intended to be used indoor.

The following InTap variants are available:

- O₂ Analyzer InTap portable (metric connectors): order number 30 425 550
- O₂ Analyzer InTap portable (imperial inch/US connectors): order number 30 457 912


2 Technical Data and Further Information

The most important technical data such as supply voltage are given on the name plate. The name plate is placed on the bottom of the device. Pay particular attention to the technical data on the nameplate.

Measurement parameter	Optical Dissolved Oxygen Saturation
Measuring range	0 to 2000 ppb
Accuracy ¹⁾ (sensor)	$\leq \pm (1 \% + 2 \text{ ppb})$
Response time t_{98} at +25 °C (+77 °F) (air to N ₂)	< 20 s
Power supply ¹⁾	<ul style="list-style-type: none">• External: 12 V DC, 2.5 A• Internal: Li-ion battery 45.4 Wh
Battery life time	Up to 24 hours
Data storage USB 2.0	8 GB
Operating temperature	<ul style="list-style-type: none">• -5 to +45 °C (+23 to +113 °F) with media• Up to +80 °C (+176 °F) for cleaning with 2% NaOH• With power adapter only up to +40 °C (+104 °F)
Storage temperature	-20 to 50 °C (-4 to +122 °F)
Operating pressure	0 to 6 bar (0 to 87 psi)
Design pressure	10 bar (145 psi)
Protection rating	IP 67
Hose connection	<ul style="list-style-type: none">• Variant with metric connectors: Swagelok 6 mm• Variant with imperial connectors: Swagelok 1/4"
Weight	3.5 kg (7.7 lbs)
Dimensions (height → width × depth)	280 mm × 170 mm × 150 mm (11.0" × 6.7" × 5.9")
Altitude	Max. 2000 m
Relative humidity	0 to 95 % non-condensing
Pollution degree	2

1) Measurement specification based on sensor specification

3 Design and Connections

- 1 Monochrome 4" touch display
- 2 On/Off button
- 3 Housing with measurement cell, O₂ InTap sensor, valve, battery, HW board
- 4 Connection for hose "Medium in"
- 5 Connection for hose "Medium out"
- 6 Rotary knob for adjusting medium flow rate
- 7 USB stick for data logging
- 8 USB type B connection for software updates
- 9 USB type A connection
- 10 Connection for AC/DC adapter for loading battery

Scope of delivery

- Power supply 30 W: AC/DC adapter with AC adapter set EU, US, AU, UK
- USB stick, plugged-in the USB connection

Supply voltage (Power)

You have the following possibilities to power the device:

- either via the internal battery
- or via the supplied AC/DC adapter connected to 100 to 240 V AC/50 to 60 Hz.

Hoses “Medium In” and “Medium Out”

You can connect the following hoses to “Medium In” and “Medium Out”:


O₂ Analyzer InTap portable		Hoses	
Description	Order number	Description	Order number
With metric connectors	30 425 550	6 mm	30 422 575
With imperial inch/ US connectors	30 457 912	1/4"	30 432 818

Adjusting medium flow rate

The flow rate is adjusted with the rotary knob. The flow rate depends on the medium pressure. The flow rate increases with increasing medium pressure.


1. Turn rotary knob clockwise to the end stop.
2. Turn rotary knob anti-clockwise to adjust the flow rate. The rotary knob has a click-stop position every 15°. E.g. for a flow rate of approximately 200 ml/min at 2 bar you have to turn the rotary knob to the 30° position. That means you hear 2 clicks.

4 Display and Operating Elements


- 1 Home screen
- 2 Measurement screen
- 3 Menu screen, here e.g. "Calibration menu"
- 4 Parametrization screen, here e.g. "Calibration InTap"
- 5 Navigation to Home screen
- 6 Navigation to next-higher menu level, here Measurement screen
- 7 Navigation to next-higher menu level, here Calibration menu screen

Display or operating element	Description
	USB stick is connected.
	Display of the current battery status
	The InTap is powered via the AC/DC adapter.
	Failure: No battery or battery is empty.
	Failure: Battery is low (< 5 %). The InTap is not powered via the AC/DC adapter.
	Failure: USB stick is not connected or damaged.
	Failure: Leakage, high temperature, battery faulty etc. Refer to the Reference Manual.
	Enter Messages screen.
	Enter menu screen.
	Enter start screen.

Display or operating element	Description
	Enter Trend display screen.
ISM	Enter ISM menu.
	Enter Calibration menu.
	Enter Configuration menu.
	Enter Measuring Point (MP) menu.
	Return to Menu screen.
	Enter next-lower menu level.
	Return to next-higher menu level.
	Change between pages within one menu level.

5 Menu Structure

6 Commissioning

The InTap is supplied ready for use. After switching on, you can immediately use the InTap as a measuring device or as a measuring device with data logging function.

At the factory the data logging is set to automatic mode. Use the parameter "Logging Mode" to change to manual mode (CAL \CONFIG \Measurement Mode).


Note that with factory settings, the display is switched off after some minutes, but the measurement will continue. Use the menu "Display Setup" to change the settings (CAL \CONFIG \Display Setup).

For communication with an In-Line sensor (process calibration) a separate Bluetooth interface T100 (order number 30 432 819) has to be mounted onto each sensor and the sensor has to be paired with the InTap.

Pairing InTap with In-Line Sensors

Prerequisite

- The In-Line sensor is within the Bluetooth range of the InTap.

- 1 A new Bluetooth device is within the Bluetooth range and was recognized. The sensor needs to be parameterized.
- 2 The In-Line sensor "MP_01" is selected. The sensor is within the Bluetooth range and already parameterized.
- 3 The In-Line sensor "MP_02" is not selected. The sensor is within the Bluetooth range and already parameterized.
- 4 The In-Line sensor "MP_03" is not selected. The sensor is already parameterized but not within the Bluetooth range.

1. Enter "Measuring Point (MP)" menu. See Chapter 4 on Page 50.
⇒ The screen "MP Select" with the item "New BT device recognized 00" is shown.
2. Tap "+" for the item "New BT device recognized 00".
3. Configure the new measuring point (MP).
4. Tap "OK"
⇒ The screen "MP Select" is shown.
5. Activate the radio button.
6. Return to the next-higher menu level.
7. Answer the question "Save changes?" with "Yes".
⇒ The Measurement Screen for the selected measuring point is shown. After some seconds the measurement results are shown. Tap "Save" to save the results on the USB stick. Tap "Repeat" to get new values from the sensor. Tap "Cancel" to cancel the procedure.

7 Calibration

You can also calibrate the dismantled InTap sensor or the dismantled In-Line Sensor via ISMCore™ software. See Operating Instructions of the ISMCore software.

As soon as the calibration is in progress no other calibration can be started.

Calibration Menu

After every successful calibration different options are available. If "Adjust" or "Calibrate" is chosen, the message "Calibration saved successfully!" is displayed. Press "Done".

Option	ISM sensors
Adjust	Calibration values are stored in the sensor and used for the measurement. Additionally, the calibration values are stored in the calibration history.
Calibrate	Calibration values are stored in the calibration history for documentation, but not be used for the measurement. The calibration values from the last valid adjustment are further used for the measurement.
Cancel	Calibration values are discarded.

Calibrating an In-Line Sensor

Prerequisite

- The In-Line sensor is paired with the InTap and is within the Bluetooth range of the InTap.
1. Enter "Calibrate In-Line Sensor" menu. Menu path: CAL > Calibrate In-Line Sensor
⇒ The screen "MP Select" is shown.
 2. Select the MP (In-Line sensor) you want to calibrate. Therefore activate the corresponding radio button.
 3. Return to the next-higher menu level.
⇒ The screen "Calibrate In-Line Sensor" is shown.
 4. Set the options "Cal Pressure", "Rel. Humidity" and "Stability".
 5. Tap "Cal".
 6. Confirm the message "There is a "auto logging" procedure running" with "Yes".
 7. Answer the question "First calibration after OptoCap™ change".
⇒ The current values for the InTap sensor and the In-Line sensor are shown.
 8. Tap "Next".
⇒ The values for P100 and P0 are shown.
 9. Tap "Adjust" or "Calibrate".
 10. Tap "Done" to finish calibration of the In-Line sensor.

Calibrating the InTap Sensor

Note: Use as air gas dry, clean and oil-free compressed air. Use as zero gas N₂ or CO₂ with a purity of at least 99.9%. For best results use a zero gas with a purity of 99.995%.

1. Connect the compressed air to the connection "Medium In".
2. Purge the InTap for 5 minutes with compressed air to remove the excess humidity out of the measurement cell. Ensure a steady medium flow without pressure build-up in the measurement cell. For this, perform two turns with the rotary knob. See Chapter 3 on Page 48.
3. Turn off compressed air.
4. Enter "Calibrate InTap Sensor" menu. Menu path: CAL > Calibrate InTap Sensor
5. Set the units.

6. Set the options "Cal Pressure", "Rel. Humidity" and "Stability".
7. Tap "Cal".
8. Confirm the message "There is a "auto logging" procedure running" with "Yes".
9. Answer the question "First calibration after OptoCap change" correctly.
⇒ The message "Connect the gas to InTap inlet." is shown.
10. Turn on the compressed air.
11. Ensure a steady medium flow without pressure build-up in the measurement cell.
12. Tap "Next".
⇒ The current value is shown.
⇒ If the criteria are met, the message "Please change gas." is shown.
13. Connect the second reference gas (zero gas) to the connection "Medium In".
14. Tap "Next".
⇒ The current value is shown.
⇒ If the stability criteria are met, the values for P100 and P0 are shown.
15. Tap "Adjust" or "Calibrate".
16. Tap "Done" to finish calibration of the InTap sensor.

8 Maintenance and Repair

Clean the surfaces with a soft damp cloth and dry the surfaces with a cloth carefully.


Replacing the OptoCap of the InTap Sensor (Sensor's Sensing Element)

Typically, the OptoCap has to be replaced after 6 to 12 months. This interval depends on the aggressiveness (CIP/SIP) of the application conditions.

For dismantling the sensor refer to the "Exchanging the InTap Sensor" section. For replacing the OptoCap refer to the User Manual "InPro 6000 Optical O₂ Sensors".

Exchanging the InTap Sensor

1. If necessary, exchange the sensor as shown in the following figure.
2. Tighten cover screws with a maximum torque of 5 Nm.
3. Calibrate sensor. See Chapter 7 on Page 54.

Replacing the Battery

The battery may only be replaced by METTLER TOLEDO. Loosening the screws of the inner cover will lead to the expiration of any warranty claims.

9 Environmental protection

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your Local Authority or retailer for recycling advice.

10 Certificates and approvals

European directives	The EU Declaration of conformity is part of the delivery.
Safety standards	<ul style="list-style-type: none">• CAN/CSA – C22.2 No. 61010-1-12• ANSI/UL Std. No. 61010-1 (3rd Edition)• Overvoltage Category II
FCC (USA)	<ul style="list-style-type: none">• CFR 47 FCC Part 15• ANSI C63.4-2014
IC (Canada)	<ul style="list-style-type: none">• ICES-003, Issue 6• ANSI C63.4-2014

InPro, InTap, ISMCore, OptoCap and ISM are trademarks of the METTLER TOLEDO Group.

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All rights reserved