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Quick Setup Guide for Sensor Head InSUS H60i for Single-Use O₂ Sensors

1 Introduction

Thank you for buying this optical oxygen sensor head from METTLER TOLEDO. The InSUS™ H60i is intended solely for in-line measurement of the oxygen partial pressure in combination with METTLER TOLEDO optical single-use O₂ sensors of the InSUS series, for example, InSUS 607. The O₂ single-use sensors are typically pre-sterilized and integrated in single-use process devices such as single-use bags and other similar applications. Please also refer to the documentation provided by your single-use device supplier.

2 Safety Instructions

 **Notice:** Before every startup, the sensor head must be checked for:

- Damage to the connections, fastenings, etc.
- Perfect functioning
- Authorization for use in conjunction with other plant equipment and resources.

 **Notice:** The manufacturer/supplier accepts no responsibility for any damage caused by unauthorized attachments or for the incorporation of spare parts which are not of METTLER TOLEDO provenance. The risk is borne entirely by the operator.

 **Caution:** Before the sensor head is put into operation, the operator must have already clarified that use of the sensor in conjunction with the other associated equipment and resources is fully authorized.

 **Caution:** A defective sensor head may neither be installed nor put into operation. Faulty containment or installation out of conformance with regulations and instructions can lead to the escape of medium or to pressure surges (explosion), potentially harmful both to persons and to the environment.

 **Caution:** Some components inside the sensor head are energized with voltages which can cause lethal shocks in case of contact. The instrument needs to be switched to zero potential before any operation on the wiring terminal.


 **Caution:** If a sensor head is broken it must be replaced immediately.

 **Caution:** It is essential that personnel involved in the operation or maintenance of the sensor head and any of the associated equipment (e.g. single-use bags, InSUS sensors, transmitters, etc.) is properly trained in the use and handling of the equipment.

 **Warning:** Before dismantling the sensor head or commencing any maintenance work on it, ensure that the equipment in which the sensor is installed is in a safe condition (e.g. power supplies switched off).

3 Product Description

The sensor heads are delivered fully assembled and factory tested for correct function. A quality control certificate is provided together with the sensor head:

- 1: VP8 connector
- 2: InSUS H60i
- 3: Clip
- 4: Sensor connector
- 5: Protective cap
- 6: InSUS O₂ sensor – not supplied with InSUS H60i

The sensor head must be connected to a single-use O₂ sensor prior to measurements.

3.1 Digital Connection of the InSUS H60i to a Transmitter Transmitter M400 or M800

The sensor head is connected to the transmitter via a VP8 cable which is available in different lengths.

Wiring information

		M400	M800 1-channel	M800 2/4-channel
Color	Function	TB3	TB3	TB2
gray	24 V	18	7	9
blue	GND24V	17	8	10
green/yellow	shield	13	2	12
brown	RS485B	14	3	13
pink	RS485A	15	4	14

3.2 Analog Connection of the InSUS H60i to a Transmitter or Biocontroller

The InSUS H60i is available with two different analog interfaces:

- InSUS H60i nA with analog nA output plus a simulated NTC 22 kΩ temperature signal.
- InSUS H60i mA with analog 4–20 mA output plus a simulated NTC 22 kΩ temperature signal.

VP8 cable for InSUS H60i nA or InSUS H60i mA

InSUS H60i nA		InSUS H60i mA
Color	Function	Function
black/transp	Cathode (nA)	4–20 mA+ HART
red	Anode (nA)	4–20 mA– HART
gray	24 V	24 V
blue	GND 24 V	GND 24 V
white	NTC 22 kΩ	NTC 22 kΩ
green	GND NTC 22 kΩ	GND NTC 22 kΩ
green/yellow	Shield	Shield

If the sensor head is used with a 4-wire transmitter (M400/M800), the sensor head is powered by the transmitter. In this case, no additional power supply is needed. For installation with 2-wire and BUS transmitters (M400/2(X)H; M400 FF; M400 PA) and for direct integration with a nA or mA signal, an additional power supply is required.

Power requirements:

$$U_{\min} = 19.5 \text{ VDC}$$

$$U_{\max} = \leq 25 \text{ VDC}$$

$$P_{\max} = 0.75 \text{ W (33...40 mA)}$$

4 Operation

4.1 Calibration

Calibration must be performed every time an InSUS O₂ single-use sensor is connected to the InSUS H60i. For most applications, a single point calibration should be sufficient. The calibration medium can be either air or a calibration gas with known O₂ concentration or a liquid medium with a known O₂ concentration. Before starting a calibration in gas, the correct pressure and the correct humidity have to be set in the transmitter. For detailed information please refer to the transmitter or biocontroller manual.

5 Maintenance

5.1 Visual Inspection

To check your sensor head, we recommend the following procedure:

- The contacts of the connector must be dry and clean: moisture, corrosion and dirt in the connector can lead to false readings.
- Check the cable for buckling, brittle areas or ruptures.
- For protection of the optical tip of the InSUS H60i from damage and dirt, reattach the protective cap to the sensor head whenever the InSUS H60i is not connected to an InSUS O₂ sensor.

6 Environmental Protection

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your Local Authority or retailer for recycling advice.

7 Warranty

On manufacturing defects, 12 months after delivery.

8 Storage Conditions

Please inspect the instrument immediately after receiving for any transport damage. The sensor head should be stored in a dry, clean and protected area until the time of installation. If the sensor head has been removed from the process, it should be thoroughly cleaned and dried. It should be stored in a dry, clean and protected area until time of installation.

9 EU/UK Declaration of Conformity

The complete declaration is available in the declaration of conformity.

InSUS is a trademark of the METTLER TOLEDO Group.

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本表依据SJ/T 11364的规定编制。本产品符合以下标志规范：
 Table composed in accordance with SJ/T 11364 (CN). This product is bearing the following symbol:

○：表示该有害物质在该部件所有均质材料中的含量均在GB/T 26572规定的限量要求以下
 ○：Indicates that the content of the hazardous substance in all homogeneous materials of the part is below the limit specified in GB/T 26572
 x：表示该有害物质至少在该部件的某一均质材料中的含量超出GB/T 26572规定的限量要求。
 x：Indicates that the content of the toxic substance in at least one of the homogeneous materials of the part exceeds the limits specified in GB/T 26572.

For addresses of METTLER TOLEDO
Market Organizations please go to:
www.mt.com/contacts

METTLER TOLEDO Group
Process Analytics
Local contact: www.mt.com/contacts

Subject to technical changes
© 09/2022 METTLER TOLEDO
All rights reserved
Printed in Switzerland. 30 802 066 A

Management System
certified according to
ISO 9001 / ISO 14001

www.mt.com/pro

* 3 0 8 0 2 0 6 6 A *