English Reference Manual **Compact Antistatic Kit** Excellence
Deutsch Referenzhandbuch **Kompaktes Antistatik-Kit** Excellence
Français Manuel de référence **Kit antistatique compact** Excellence
Español Manual de referencia **Kit antiestático compacto** Excellence
Italiano Manuale di riferimento **Kit antistatico compatto** Excellence

METTLER TOLEDO

Declaration of Conformity

EC Declaration of Conformity

This declaration of conformity is issued under the sole responsibility of the manufacturer:

Fraser Anti-Static Techniques Ltd
Scotts Business Park, Bampton GB-EX16 9DN

who declares that the following equipment:

Mettler Toledo Compact Ionizer Part 30090337

is in conformity with the relevant European Community harmonisation legislation:

Low Voltage Directive 2006/95/EC
EMC Directive 2004/108/EC
RoHS2 Directive 2011/65/EU

Standard IEC 601010
Standard IEC 61326-1

Signed for and on behalf of the manufacturer

A handwritten signature in black ink, appearing to read 'Robert Fraser'.

Robert Fraser, Director

5th September 2013

CE marking first affixed in September 2013.

Reference Manual **Compact Antistatic Kit**

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Italiano

1 Introduction

Your Antistatic Kit is an accessory for METTLER TOLEDO analytical balances and is used to discharge electrostatically charged objects to be weighed.

The guidelines stated in the Reference Manual of the METTLER TOLEDO balance connected to your ionizer also fully apply to your Antistatic Kit.

1.1 Explanation of conventions and symbols used

Conventions and symbols

Key and button designations are indicated by a picture or text in square brackets, e.g., [] or [Define].

Elements of instructions

In this manual, step-by-step instructions are presented as follows. The action steps are numbered and can contain prerequisites, intermediate results and results, as shown in the example. Sequences with less than two steps are not numbered.

- Prerequisites that must be fulfilled before the individual steps can be executed.

- 1 Step 1
 - ➔ Intermediate result
- 2 Step 2
 - ➔ Result

2 Safety Information

2.1 Definition of signal words and symbols

Safety notes are indicated by signal words and warning symbols and contain warnings and information about safety issues. Ignoring safety notes can lead to personal injury, damage to the instrument, malfunctions and erroneous results.

Signal words

DANGER	A hazardous situation with high risk, resulting in death or severe injury if not avoided.
WARNING	A hazardous situation with medium risk, possibly resulting in death or severe injury if not avoided.
CAUTION	for a hazardous situation with low risk, resulting in minor or moderate injury if not avoided.
NOTICE	A hazardous situation with low risk, resulting in damage to the instrument, other material damage, malfunctions and erroneous results, or loss of data.

Warning symbols

General hazard

Electrical shock

2.2 Product-specific safety notes

Intended use

Your Antistatic Kit is an accessory for METTLER TOLEDO analytical balances and is used to discharge electrostatically charged objects. Use the Antistatic Kit exclusively for this purpose. Any other type of use and operation beyond the limits of technical specifications without written consent from Mettler-Toledo GmbH, is considered as not intended.

Safety notes

⚠ WARNING

Risk of electric shock

- 1 The ionizer must be disconnected from the power supply, before cleaning or other maintenance work.
- 2 Use only the power cable from METTLER TOLEDO, if it needs replacing.
- 3 Ensure that no liquid comes into contact with the balance, the terminal, the ionizer or the AC adapter.
- 4 Never open the balance, terminal, ionizer or AC adapter – they contain no components which can be cleaned, repaired or replaced by the user.
- 5 Do not turn on power until the cleaned parts are dry.


⚠ CAUTION

Damage of device

For use only in dry interior rooms.


3 Overview

Antistatic Kit

1	Ionizer	6	Reference Manual
2	Indicator light	7	Power connection
3	Side door (marked with "Conducting Glass")	8	Control cable
4	AC adapter of balance	9	Roller Antistatic Kit
5	Power supply cable of balance		

Balance back

10	Aux 1 (connection for ionizer control)	12	Socket for AC adapter
11	Aux 2 (connection for ionizer control)		

4 Installation and Putting into Operation

The compact ionizer produces a corona of ionized air. Ionized air provides ions of the opposite polarity to neutralize static charge in the object being weighed.

4.1 Scope of delivery

- Compact ionizer
- Side door marked with "Conducting Glass"
- Reference Manual

4.2 Assembling the Antistatic Kit

⚠ WARNING

Risk of electric shock

- 1 Use only an approved METTLER TOLEDO AC adapter which has an earthed supply. The use of other SELV adapters without an earthed supply could result in dangerous operation.
- 2 The emitter pin is resistively coupled to the high voltage electronics. The maximum current available from it is less than 40 μ A.
- 3 Use only the 3-pin power cable with equipment grounding conductor which was supplied with your balance.
- 4 Only a 3-pin grounded outlet for connecting your balance must be used.
- 5 Only extension cables which meet this relevant standards and also have an equipment grounding conductor may be used.
- 6 Intentional disconnection of the equipment grounding conductor is prohibited.


The ionizer should arrive in good condition, protected by its packaging.

- 1 Inspect the ionizer and all cables for damage caused in transit.
- 2 Report any damage immediately.
- 3 Do not use if there is any damage to the cables, connectors or the ionizer body.
- 4 Guide the cables so that they cannot become damaged or interfere with the measuring process.
- 5 Remove protective cap without squeezing.

4.2.1 Installing Antistatic Kit


- Balance is switched off.
 - Balance is disconnected from the power supply.
- 1 Remove the relevant side door (1).
 - 2 Remove any of the 3 clips, for example the middle one (2).


- 3 Push the compact ionizer in from back to front until it clicks into place. If necessary, gently squeeze the strips at the back.

- 4 Insert the side door (3) marked with "Conducting Glass" (4) at an angle of approx. 30° into the two slots.
- 5 Lift the side door against the balance until it clicks into place.
- 6 Push the side door completely to the front.
- 7 Check that the side door runs easily, otherwise it is not correctly inserted.
- 8 For XPE balances, return the handle to the upright position.
For XSE balances, clip the handle together.


- 9 Plug the power connection (5) into the socket for AC adapter.
- 10 Plug the control cable (6) into Aux 1 or Aux 2.
- 11 Plug the power supply cable (7) into the ionizer.
- 12 Connect the AC adapter (8) to the power supply.

4.2.2 Installing second Antistatic Kit

For the installation of a second Antistatic Kit you will also need a second AC adapter, see [Accessories ▶ Page 11].


- 1 Remove the relevant side door (1).
- 2 Remove the any of the 3 clips, for example the middle one (2).


- 3 Push the compact ionizer in from back to front until it clicks into place. If necessary, gently squeeze the strips at the back.

- 4 Insert the side door (3) marked with "Conducting Glass" (4) at an angle of approx. 30° into the two slots.
- 5 Lift the side door against the balance until it clicks into place.
- 6 Push the side door completely to the front.
- 7 Check that the side door runs easily, otherwise it is not correctly inserted.
- 8 For XPE balances, return the handle to the upright position.
For XSE balances, clip the handle together.

- 9 Plug the control cable (5) into Aux 1 or Aux 2.
- 10 Plug the power supply cable (6) into the ionizer.
- 11 Connect the second AC adapter (7) to the power supply.

4.2.3 XPE balance settings for Antistatic Kit

The arrow buttons can be used to page forward or back to a menu page.

- 1 Switch on the balance [⏻].
- 2 Press [↩].
- 3 Choose **Smart & ErgoSens** and tap the button [Define].
 - ➔ A selection window appears.
- 4 Choose **ErgoSens 1 (Aux1)** or **ErgoSens 2 (Aux2)** and tap on the associated button.
 - ➔ A selection window appears.
- 5 Scroll and select [AntiStatic Kit].
- 6 To save the settings, tap [OK].
- 7 Choose **Ionizer Setup** and tap the button [Define].
- 8 Choose settings appropriate to your application. Recommended settings are active and 0 seconds **After door closed**.
- 9 To save the settings, tap [OK].
 - ➔ Indicator light turns on and Status-Icon [🔌] appears in the display when the ionizer is working.

4.2.4 XSE balance settings for Antistatic Kit

The arrow buttons can be used to page forward or back to a menu page.

- 1 Switch on the balance [⏻].
- 2 Press [↩].
- 3 Choose **ErgoSens** and tap the button [Define].
 - ➔ A selection window appears.
- 4 Choose **ErgoSens 1 (Aux1)** or **ErgoSens 2 (Aux2)** and tap on the associated button.
 - ➔ A selection window appears.
- 5 Tap in the scroll bar to select [AntiStatic Kit].
- 6 To save the settings, tap [OK].
- 7 Choose **Function Keys** and tap the button [Define].
 - ➔ A selection window appears.
- 8 Tap [Ionizer].
- 9 To save the settings, tap [OK].
 - ➔ The function key [Ionizer] is displayed in the application at the bottom of the display.
 - ➔ Indicator light turns on and Status-Icon [🔌] appears in the display when the Ionizer is working.

5 Operation

⚠ CAUTION

Risk of ignition or sparking

The ionizer must NOT be used inside an inert gas atmosphere or with an inert gas purging system. There is a serious risk of sparks from the ionizer if such an environment is used.

The ionizer is controlled by the balance. It should be activated during loading the balance and deactivated during weighing.

Activating with XPE

With XPE balances the ionizer is controlled by the automatic door.

Activating with XSE

- To activate the ionizer, tap function key **[Ionizer]**.
 - ⇒ The indicator light on the ionizer turns on.
- To deactivate the ionizer, tap again function key **[Ionizer]**.
 - ⇒ The indicator light turns off.

If the indicator light doesn't turn on:

- 1 Check cable connections and settings.
 - 2 If the indicator light remains off, contact your local METTLER TOLEDO service representative.
- Further information, consult your balance Reference Manual.

6 Maintenance

6.1 Cleaning

⚠ WARNING

Risk of electric shock

- 1 The ionizer must be disconnected from the power supply, before cleaning or other maintenance work.
- 2 Use only the power cable from METTLER TOLEDO, if it needs replacing.
- 3 Ensure that no liquid comes into contact with the balance, the terminal, the ionizer or the AC adapter.
- 4 Never open the balance, terminal, ionizer or AC adapter – they contain no components which can be cleaned, repaired or replaced by the user.
- 5 Do not turn on power until the cleaned parts are dry.

Clean the ionizer with a soft brush or with ethyl alcohol and a lint-free cloth.

Quick cleaning

Clean the ionizer once a month.

- 1 Dampen the lint-free cloth or foam cotton swab with a small amount of ethyl alcohol.
- 2 **NOTICE: Damage or malfunction** Do not touch the emitter pin with any hard object.
NOTICE: Malfunction due to contamination Avoid touching the emitter pin with fingers. Wear gloves.
 Wipe the inside and outside of the insulation (white PTFE) around the emitter pin (1). Very gently wipe the emitter pin.

- ➔ The PTFE insulation and emitter pin are clean.
- 3 Check that no ethyl alcohol is left on the ionizer.

Cleaning

Thoroughly clean the ionizer every three months.

- 1 Deinstall the ionizer.
- 2 Dampen the lint-free cloth or foam cotton swab with a small amount of ethyl alcohol.
- 3 **NOTICE: Damage or malfunction** Do not touch the emitter pin with any hard object.
NOTICE: Malfunction due to contamination Avoid touching the emitter pin with fingers. Wear gloves.
 Wipe the inside and outside of the insulation (white PTFE) around the emitter pin (1). Very gently wipe the emitter pin.
 ➔ The PTFE insulation and emitter pin are clean.
- 4 Clean the insulation tube (3) between the steel emitter barrel (2) and the body of the ionizer.
- 5 Check that no ethyl alcohol is left on the ionizer.
- 6 Reinstall the ionizer, see [Installing Antistatic Kit ▶ Page 6].

Contact your METTLER TOLEDO service representative for details about available service options.

7 Technical Data

7.1 General data

⚠ CAUTION

Use only with a tested AC Adapter with SELV output current.
 Ensure correct polarity

Ionizer

Input voltage:	9 V to 15 V DC
Input current (operating):	560 mA
Weight:	150 g (complete with packaging)

Protection and standards

Overvoltage category:	II
Degree of pollution:	2
Range of application:	For use only in closed interior rooms

Environmental conditions


Height above mean sea level:	up to 4000 m
Ambient temperature:	5-40 °C
Relative air humidity:	Max. 80% up to 31 °C, linearly decreasing to 50% at 40 °C, noncondensing
Atmosphere:	Do not use under inert gas atmosphere

Ozone

All HV ionizers produce a trace of ozone. The METTLER TOLEDO ionizer ozone level is considerably below the international safety limit of 0.1 ppm.

7.2 Dimensions

Dimensions in mm.

8 Disposal


In conformance with the European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions, please contact the responsible authority or the distributor from which you purchased this device. Should this device be passed on to other parties, the content of this regulation must also be related.

9 Accessories and Spare Parts

9.1 Accessories


Description

AC/DC adapter (without power cable)

- Input: 100 – 240 V AC, 50/60 Hz, 0.8 A
- Output: 12 V DC, 2.5 A

Order no.

11107909

Country-specific 3-Pin power cable with grounding conductor.

Power cable AU	00088751
Power cable BR	30015268
Power cable CH	00087920
Power cable CN	30047293
Power cable DK	00087452
Power cable EU	00087925
Power cable GB	00089405
Power cable IL	00225297
Power cable IN	11600569
Power cable IT	00087457
Power cable JP	11107881
Power cable TH, PE	11107880
Power cable US	00088668
Power cable ZA	00089728

9.2 Spare Parts

	Description	Order no.
	Side door "Conducting Glass"	11106263
	1 Roller Antistatic Kit, 2 pcs	30449902

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© Mettler-Toledo GmbH 11/2021