Thermal Analysis Excellence

Thermal Analysis

Innovative Technology

Versatile Modularity

Swiss Quality


Thermal Analysis

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METTLER TOLEDO

DSC3

Unsurpassed Sensitivity
Excellent Resolution


SmartSens terminal with One Click™ function

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
Temperature data

| | | | |
|---------------------------------------|-------------------------|---------------------------|-------------------------|
| Temperature range | air cooling | RT to 500 °C (200W) | RT to 700 °C (400W) |
| | cryostat cooling | -50 to 450 °C | -50 to 700 °C |
| | intraCooler | -100 to 450 °C | -100 to 700 °C |
| | liquid nitrogen cooling | -150 to 500 °C | -150 to 700 °C |
| Temperature accuracy | | ± 0.2 K | |
| Temperature precision | | ± 0.02 K | |
| Furnace temperature resolution | | ± 0.00006 K | |
| Heating rate RT to 700 °C | | 0.02 to 300 K/min | |
| Cooling rate | | 0.02 to 50 K/min | |
| Cooling time | air cooling | 8 min (500 to 100 °C) | 9 min (700 to 100 °C) |
| | cryostat cooling | 5 min (100 to 0 °C) | |
| | Intracooler | 5 min (100 to 0 °C) | |
| | liquid nitrogen cooling | 15 min (100 to -100 °C) | |

Colorimetric data

| | | | |
|--------------------------------------|---------------------|----------|----------|
| Sensor type | FRS 5+ | HSS 8+ | |
| Sensor material | Ceramic | | |
| Number of thermocouples | 56 | 120 | |
| Signal time constant | 1.8 s | 3.1 s | |
| Indium peak (height to width) | 17 | 6.9 | |
| TAWN | resolution | 0.12 | 0.20 |
| | sensitivity | 11.9 | 56.0 |
| Measurement range | at 100 °C | ± 350 mW | ± 160 mW |
| | at 700 °C | ± 200 mW | ± 140 mW |
| Resolution | 0.04 μW | 0.02 μW | |
| Digital resolution | 16.8 million points | | |

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
Temperature data

| | | | |
|--------------------------------|-------------------------|---------------------------|-------------------------|
| Temperature range | air cooling | RT to 500 °C (200W) | RT to 700 °C (400W) |
| | cryostat cooling | -50 to 450 °C | -50 to 700 °C |
| | intraCooler | -100 to 450 °C | -100 to 700 °C |
| | liquid nitrogen cooling | -150 to 500 °C | -150 to 700 °C |
| Temperature accuracy | | ± 0.2 K | |
| Temperature precision | | ± 0.02 K | |
| Furnace temperature resolution | | ± 0.00006 K | |
| Heating rate RT to 700 °C | | 0.02 to 300 K/min | |
| Cooling rate | | 0.02 to 50 K/min | |
| Cooling time | air cooling | 8 min (500 to 100 °C) | 9 min (700 to 100 °C) |
| | cryostat cooling | 5 min (100 to 0 °C) | |
| | Intracooler | 5 min (100 to 0 °C) | |
| | liquid nitrogen cooling | 15 min (100 to -100 °C) | |

Colorimetric data

| | | | |
|-------------------------------|--------------------------|----------|----------|
| Sensor type | FRS 6+ | HSS 9+ | |
| Sensor material | Ceramic | | |
| Number of thermocouples | 56 | 120 | |
| Signal time constant | 1.8 s | 3.1 s | |
| Indium peak (height to width) | row data | 19.5 | 6.9 |
| | mathematically corrected | >155 | 0.20 |
| TAWN | resolution | 0.12 | 56.0 |
| | sensitivity | 11.9 | ± 160 mW |
| Measurement range | at 100 °C | ± 350 mW | ± 140 mW |
| | at 700 °C | ± 200 mW | 0.02 μW |
| Resolution | 0.04 μW | | |
| Digital resolution | 16.8 million points | | |

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
Examples thermal events and processes that can be determined by DSC

- | | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <ul style="list-style-type: none"> • Melting behavior • Crystallization • Polymorphism • Liquid-crystalline transitions • Phase diagrams and composition • Glass transitions • Reactivity • Reaction kinetics | <ul style="list-style-type: none"> • Curing • Stability • Miscibility • Effects of plasticizers • Thermal history • Heat capacity and heat capacity changes • Reaction and transition enthalpies • Purity |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|


Application Area : Thermoplastics, Thermosets, Elastomers, Adhesives, Safety technology, Electronics, Aircrafts, Automotives, Foodstuffs, Pharmaceuticals, Chemicals and composites

Application

TGA2

High Performance
Already Built into the Basic Configuration

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

Sample Robot Option - Automatic weighing in (pan and sample, using sample robot)

Buoyancy Compensation

Temperature data

| | Small furnace (SF) | Large furnace (LF) |
|--------------------------|-----------------------------|-----------------------------|
| Temperature range | RT to 1100 °C | RT to 1100 °C |
| Temperature accuracy | ± 1 K | ± 1 K |
| Temperature precision | ± 0.4 K | ± 0.6 K |
| Heating rate | 0.02 to 250 K/min | 0.02 to 150 K/min |
| Cooling time | 20 min (1100 to 100 °C) | 22 min (1100 to 100°C) |
| Cooling time with helium | ≤ 10 min (1100 to 100 °C) | ≤ 11 min (1100 to 100 °C) |
| Sample volume | ≤ 100 μl | ≤ 900 μl |

Balance data

| | XP1 | XP1U | XP5 | XP5U |
|-----------------------------|------------------------------------------------------|-------------|------------|-------------|
| Measurement range | ≤ 1 g | ≤ 1 g | ≤ 5 g | ≤ 5 g |
| Resolution | 1.0 μg | 0.1 μg | 1.0 μg | 0.1 μg |
| Weighing accuracy | 0.005% | 0.005% | 0.005% | 0.005% |
| Weighing precision | 0.0025% | 0.0025% | 0.0025% | 0.0025% |
| Repeatability | < 0.001 mg | < 0.0008 mg | < 0.002 mg | < 0.0009 mg |
| Typical Minimum Weight | 0.19 mg | 0.16 mg | 0.22 mg | 0.17 mg |
| Typical Minimum Weight USP | 1.9 mg | 1.6 mg | 2.2 mg | 1.7 mg |
| Internal ring weights | 2 | | | |
| Blank curve reproducibility | better than ± 10 μg over the whole temperature range | | | |

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Examples thermal events and processes that can be determined by TGA2

- Quantitative content analysis (moisture, fillers, polymer content, materials, etc.)
- Adsorption and desorption of gases
- Kinetics of decomposition processes
- Sublimation, evaporation and vaporization
- Thermal stability
- Oxidation reactions and oxidation stability
- Identification of decomposition products, solvents and solvates
- Sorption and desorption of moisture
- Pseudopolymorphism
- Determination of Curie temperatures

Application Area : Rubbers, Plastics, Building materials, Minerals, ceramics, Chemicals, Pharmaceuticals and foodstuffs

Application

TGA/DSC3+

High Performance

Already Built into the Basic Configuration

METTLER TOLEDO Balances – Uniqueness in all Respects

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Temperature data

| | Small furnace (SF) | Large furnace (LF) | High temp. furnace (HT) |
|--------------------------------|---------------------------|---------------------------|---------------------------|
| Temperature range | RT to 1100 °C | RT to 1100 °C | RT to 1600 °C |
| Temperature accuracy | ±0.25 K | ±0.3 K | ±0.5 K |
| Temperature precision | ±0.15 K | ±0.2 K | ±0.3 K |
| Furnace temperature resolution | 0.001 K | 0.001 K | 0.002 K |
| Heating time | 5 min (RT to 1100 °C) | 10 min (RT to 1100 °C) | 10 min (RT to 1600 °C) |
| Cooling time | 20 min (1100 to 100 °C) | 22 min (1100 to 100 °C) | 27 min (1600 to 100 °C) |
| Cooling time with helium | ≤ 10 min (1100 to 100 °C) | ≤ 11 min (1100 to 100 °C) | ≤ 13 min (1600 to 100 °C) |
| Heating rate | 250 K/min | 150 K/min | 100 K/min |
| Cooling rate | -20 K/min (≥ 150 °C) | -20 K/min (≥ 150 °C) | -20 K/min (≥ 200 °C) |
| Sample volume | ≤ 100 µl | ≤ 900 µl | ≤ 900 µl |

Balance data

| | XP1 | XP1U | XP5 | XP5U |
|-----------------------------|------------------------------------------------------|-------------|------------|-------------|
| Measurement range | ≤ 1 g | ≤ 1 g | ≤ 5 g | ≤ 5 g |
| Resolution | 1.0 µg | 0.1 µg | 1.0 µg | 0.1 µg |
| Weighing accuracy | 0.005% | 0.005% | 0.005% | 0.005% |
| Weighing precision | 0.0025% | 0.0025% | 0.0025% | 0.0025% |
| Repeatability | < 0.001 mg | < 0.0008 mg | < 0.002 mg | < 0.0009 mg |
| Typical Minimum Weight | 0.19 mg | 0.16 mg | 0.22 mg | 0.17 mg |
| Typical Minimum Weight USP | 1.9 mg | 1.6 mg | 2.2 mg | 1.7 mg |
| Internal ring weights | 2 | | | |
| Blank curve reproducibility | better than ± 10 µg over the whole temperature range | | | |

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Examples thermal events and processes that can be determined by TGA/DSC 3+

TGA

- Quantitative content analysis (moisture, fillers, polymer content, materials, etc.)
- Adsorption and desorption of gases
- Kinetics of decomposition processes
- Sublimation, evaporation and vaporization
- Thermal stability
- Oxidation reactions and oxidation stability
- Identification of decomposition products, solvents and solvates
- Sorption and desorption of moisture
- Pseudopolymorphism
- Determination of Curie temperatures

DSC

- Melting behavior
- Crystallization
- Polymorphism
- Phase diagrams
- Glass transitions
- Reaction kinetics
- Heat capacity
- Reaction and transition enthalpies

Application Area : Rubbers, Plastics, Building materials, Minerals, ceramics, Chemicals, Pharmaceuticals and foodstuffs

Application

TMA/SDTA 1

Unparalleled versatility
Optimum Measurement Configurations

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

Analysis of decomposition products by EGA (Evolved Gas Analysis)

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
Examples thermal events and processes that can be determined by TGA/SDTA1

- Viscoelastic behavior (Young's modulus)
- Glass transition
- Expansion coefficient
- Expansion and shrinkage of fibers and films
- Softening
- Viscous flow
- Melting and crystallization
- Gelation
- Phase transitions
- Curing and crosslinking reaction
- Swelling behavior
- Volume expansion
- Thermal effects of pharmaceuticals and foodstuffs


Application Area : Paints, Plastics, Lacquers, Composites, Films, Fibers, Ceramics, Coatings and Adhesives

Application

DMA/SDTA861^e

Excellent Performance
Over the wide force range

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Sample Holders The various deformation Modes

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- 3-point bending
- Single cantilever bending
- Dual cantilever bending
- Tension
- Compression
- Shear

Temperature data

| | |
|----------------------|----------------|
| Range | -150 to 500 °C |
| Technical resolution | 0.003 K |
| Accuracy | 0.5 K |

Force

| | |
|----------------------|--------------------------------------------|
| Range | 0.001 to 40 N, (12, 18 or 40 N) |
| Technical resolution | 0.15 mN (0 to 5 N), 1.5 mN (0 to 50 N) |
| Sensitivity | 1 mN |

Displacement

| | |
|----------------------|----------|
| Range | ± 1.6 mm |
| Technical resolution | 0.6 nm |
| Sensitivity | 5nm |

Stiffness

| | |
|-----------|---------------|
| Range | 10 to 108 N/m |
| Precision | 0.2% |

Tan delta

| | |
|----------------------|---------------|
| Range | 0.0001 to 100 |
| Technical resolution | 0.00001 |
| Sensitivity | 0.0001 |

Frequency

| | |
|----------------------|----------------------------------------------------------------------------------------------------|
| Range | 0.001 to 1000Hz |
| Technical resolution | 0.00001 |
| Accuracy | 0.0001 |
| Modes | Logarithmic and linear scans Multi-frequency (sequentially) Multi-frequency (simultaneously) |

Measurement modes

| | |
|---------------------------------|------------------------------------------------------------------------------------------------------------------|
| Bending 3-point Dual cantilever | Length : 30 to 90 mm, Length : 20 to 80 mm Width : < 15 mm, Thickness : < 5 mm Max. sample length : 100 mm |
| Stiffness range bending | 10 to 10 ⁸ N/m |
| Shear | Diameter : ≤ 15 mm, Thickness : ≤ 6.5 mm |
| Stiffness range shear | 10 to 10 ⁸ N/m |
| Tension | Length : 19.5 mm (19.5, 10.5, 9.0, 5.5) Width : ≤ 7 mm, Thickness : ≤ 3 mm |
| Stiffness range tension | 10 to 10 ⁷ N/m |
| Compression | Diameter : ≤ 20 mm, Thickness : ≤ 9 mm |
| Stiffness range compression | 10 to 10 ⁷ N/m |

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
Examples thermal events and processes that can be determined by DMA/SDTA861e

- | | |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <ul style="list-style-type: none"> • Viscoelastic behavior • Relaxation behavior • Glass transition • Mechanical modulus • Damping behavior • Softening • Viscous material flow | <ul style="list-style-type: none"> • Crystallization and melting • Gelation • Phase transitions • Composition of blends • Curing and polymerization reactions • Material defects • Effects caused by filler materials |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

Application Area : Thermoplastics, Thermosets, Elastomers, Films, Fiber, Ceramics, Metals, Aircrafts, Automotives, Viscous materials and composites


Application

DMA1

Excellent Performance

The Perfect Solution for Materials Analysis

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
Temperature data

| | |
|----------------------|-----------------------|
| Range | -190 to 600 °C |
| Technical resolution | 0.1 K |
| Temperature Accuracy | 0.75 K |
| Heating rate | 0.1 K/min to 20 K/min |
| Cooling rate | 0.1 K/min to 30 K/min |

Force data

| | |
|----------------------|-------------------|
| Force Range | ± 0.001 to ± 10 N |
| Technical resolution | 0.25 mN |
| Sensitivity | 1 mN |

Displacement

| | |
|----------------------|--------|
| Displacement Range | ± 1 mm |
| Technical resolution | 2 nm |
| Sensitivity | 30 nm |

Stiffness

| | |
|-----------------|---------------|
| Stiffness Range | 50 to 105 N/m |
| Precision | 0.5% |

Tan delta

| | |
|----------------------|--------------|
| Tan delta Range | 0.0001 to 50 |
| Technical resolution | 0.00001 |
| Sensitivity | 0.0001 |

Frequency

| | |
|----------------------|----------------------------------------------------------------------|
| Frequency Range | 0.001 to 300 Hz |
| Technical resolution | 0.0001 |
| Accuracy | 0.001 |
| Modes | 1. Logarithmic and linear scans 2. Multi-frequency (sequentially) |

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
Examples thermal events and processes that can be determined by DMA1

- Viscoelastic behavior
- Relaxation behavior
- Glass transition
- Mechanical moduli
- Damping behavior
- Softening
- Viscous flow
- Crystallization and melting
- Gelation
- Phase transformations
- Composition of blends
- Curing and polymerization reactions
- Material defects
- Effects due to fillers

Application Area : Thermoplastics, Thermosets, Elastomers, Films, Fiber, Ceramics, Metals, Aircrafts, Automotives, Viscous materials and composites

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

- Thermal Analysis in Practice
- Thermal Analysis of Polymers
- Evolved Gas Analysis

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Quality certificate. Development, production and testing according to ISO 9001.

Environmental management system according to ISO 14001.

“European conformity”. The CE conformity mark provides you with the assurance that our products comply with the EU directives.