English

Installation Instructions **Antistatic Kit** for analytical balances

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METTLER TOLEDO

1 Target Group

These Installation Instructions are addressed to the instrument owner. To ensure proper ground connection, the high-voltage discharging power pack must be installed by a qualified electrician.


⚠ WARNING

Death or serious injury due to electric shock

The discharging power pack is a high-voltage device. Additional safety precautions must be followed when using the device.


- Make sure to read the operating instructions of the discharging power pack (manufactured by HAUG) before installing the antistatic kit.

2 Scope of Delivery

1	U-electrode	3	Power cable
2	Discharging power pack		

3 Interface of the Discharging Power Pack

1	On/Off switch	3	Sockets for electrodes (high voltage)
2	Socket for ground connector	4	Socket for power cable

4 Selecting the Location

The requirements on the location of the balance are described in the Reference Manual (RM) of the balance. In addition, the location of the discharging power pack must comply to the following requirements:

- Do not install the discharging power pack on a surface generating or radiating heat.
- Do not install the discharging power pack under direct sunlight.
- Do not install the discharging power pack in areas with potentially explosive atmospheres.
- Do not keep any inflammable materials in the vicinity of the discharging power pack or its components.

5 Unpacking

Check the package, the packaging elements, and the delivered components for damage or dirt. Clean any impurities, in particular at the high-voltage sockets. If any component is damaged, do not install, secure against inadvertent use, and contact your METTLER TOLEDO service representative.

⚠ WARNING

Death or serious injury due to electric shock or fire

Unauthorized modifications or damage to the discharging power pack may result in sparks.

- Do not open or modify the discharging power pack.

⚠ WARNING

Death or serious injury due to electric shock

Short circuits can occur as a result of soiling at the high-voltage connection sockets.

- 1 Use the provided blind plugs to protect the unused high-voltage connection sockets.
- 2 Ensure that the blind plugs are clean, dry, and free of grease.

6 Installing


⚠ WARNING

Death or serious injury due to electric shock


Contact with parts that carry a live current can lead to death or injury.

- Connect the power cable to a grounded power outlet.


- 1 Position the U-electrode (1) about 5 cm away from the left or right door of the balance.
- 2 Connect the U-electrode to the discharging power pack (2). The discharging power pack can have up to two U-electrodes connected to it.


- 3 Connect the power cable (3) to the discharging power pack (2).

- 4 **⚠ WARNING: Electric shock. Make sure that the discharging power pack is switched off.**
Connect the power cable (3) to the power outlet.

- 5 **This step has to be performed by a qualified electrician.**
Connect the ground socket (4) of the discharging power pack to ground potential complying with applicable standards. Use a ground cable with a minimum conductor cross-section of 1.5 mm².

⚠ WARNING

Death or serious injury due to electric shock

- 1 Always switch off the discharging power pack before disconnecting from power or disconnecting the U-electrode.
- 2 Always switch off the discharging power pack after use.

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GWP®

Good Weighing Practice™

GWP® is the global weighing standard, ensuring consistent accuracy of weighing processes, applicable to all equipment from any manufacturer. It helps to:

- Choose the appropriate balance or scale
- Calibrate and operate your weighing equipment with security
- Comply with quality and compliance standards in laboratory and manufacturing

 www.mt.com/GWP

www.mt.com/accessories

For more information

Mettler-Toledo GmbH

Im Langacher 44
8606 Greifensee, Switzerland
www.mt.com/contact

Subject to technical changes.
© Mettler-Toledo GmbH 04/2022
30491869C

30491869