

Pipette Calibration Order Form

800 4 Rainin (800 472 4646) RaininService.com

Step 1 Service Plans and Specifications

Product Name	Standard	Advanced	Accredited ISO 17025*	Accredited+ ISO 8655 and ISO 17025*
As-found Calibration (volumes x weighings)	n/a	2 x 4	3 x 4	3 x 10
As-returned Calibration (volumes x weighings)	2 x 4	2 x 4	3 x 10	3 x 10

Pipette Type	QTY	17800031	17800040	17800048	30829339
Single channel		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
8 channels (all)		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
12 channels (all)		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Specifications	As Found	As Returned
ISO 8655	<input type="checkbox"/>	<input type="checkbox"/>
Manufacturer	<input type="checkbox"/>	<input type="checkbox"/>
Custom**	<input type="checkbox"/>	<input type="checkbox"/>

*Includes Measurement Uncertainty, which is not Included in PASS/FAIL determination.
**Must be supplied with each order and may not exceed the capability of the pipette.

Step 2 Service Interval

Annual 6 Months 3 Months Other _____

Step 3 Optional Services (Additional fees may apply)

- DNA and RNase-off decontamination
 24 hour rush service (call 800 472 4646 to schedule)
 ISO 17025 Certificate

Step 4 Decontamination Policy

If your organization maintains a current master service agreement with Mettler-Toledo Rainin, LLC and an acceptable decontamination process is specified within such agreement, no signature is required on this form for certification of decontamination.

I certify that my pipettes are free of biological, chemical and radiological contaminants.

Signed _____ Date _____

Credit card information will be destroyed after processing.

Step 5 Payment Method

Visa Amex MasterCard Purchase Order (please attach)

Credit Card or PO No. _____

Credit Card Expiration Month: _____ Year: _____	Your Reference Number (optional)
Rainin Use Only	Promo Code/Voucher No.

Step 6 Billing

Company _____
 Cardholder Name _____
 Address 1 _____
 Address 2 _____
 Address 3 _____
 City, State, ZIP _____
 Quote No. _____

For current pricing, visit www.mt.com/RaininServicePricing

I need an estimate if additional parts/repairs are needed.

Note: We will gladly design a custom plan to meet your unique needs.
 For complete Terms & Conditions, visit www.mt.com/CaITC

Orders without payment information will not be processed.

Step 7 Shipping

Contact Name _____
 Phone Number _____
 Email _____

Shipping address same as Billing address

Organization _____
 Address 1 _____
 Address 2 _____
 Address 3 _____
 City, State, ZIP _____

By submitting this service request I agree to the METTLER TOLEDO standard Terms and Conditions of Sale, which can be accessed of www.mt.com/legal and are incorporated herein by reference. If your organization maintains a current master service agreement with Mettler-Toledo, LLC, the Terms and Conditions of that master agreement will supersede the standard terms and Conditions of Sale as described above and referenced within this pipette calibration order form.

Step 8 Service Centers

7500 Edgewater Drive, Oakland, CA 94621
 150 Wells Avenue, Newton, MA 02459
 5955 Mira Mesa Boulevard, Suite A, San Diego, CA 92121
 200 Rittenhouse Circle East, Unit 2, Bristol, PA 19007

Standard Turnaround

1-25 single channels is 48 hours (in house)
 26+ single/multi-channels is 72 hours (in house)

Special Instructions:

RaininService Convenient Mail-in Calibration Service

Rainin makes it easy to get fast reliable service and calibration, guaranteed!

- Quick turnaround, with four ISO 17025 accredited service labs across North America
- We service all brands of pipettes – most with manufacturers’ original parts (Rainin®, Eppendorf®, Gilson®, Integra®, Thermo Fisher®, Sartorius®)
- We service all pipettes – single channels, multichannels and all Rainin benchtop pipettors
- Free inbound shipping
 - 4 or more, request our free return kit
 - 3 or less, pack them yourself and print a free shipping label at [mt.com/RaininServiceLabels](https://www.mt.com/RaininServiceLabels)

Service Plans

All plans include:

- Full preventive maintenance (PM) with seal replacement or restoration
- Liquid leak test on all nozzles
- Adjustment to manufacturer’s tolerances
- Calibration certificate with detailed results
- Calibration label

	Standard Plan	Advanced Plan	Accredited Plan	Accredited+ Plan
	Labs with basic service needs Perfect for Academic R&D Research Labs	When prior performance matters Pharma Biotech Clinical Research Labs	For complete ISO 17025 documentation GLP/GMP, Regulated Labs	ISO 8655 and 17025 compliance GLP/GMP, Regulated Labs
As Found	N/A	4 Test weighings @ 10, 100 %	4 Test weighings @ 10, 50, 100 %	10 Test weighings @ 10, 50, 100 %
As Returned	4 Test weighings @ 10, 100 %	4 Test weighings @ 10, 100 %	10 Test weighings @ 10, 50, 100 %	10 Test weighings @ 10, 50, 100 %
Preventive Maintenance	Yes	Yes	Yes	Yes
Calibration Certificate	Non-accredited	Non-accredited	Accredited	Accredited
Compliance	ISO 8655-7 Annex A	ISO 8655-7 Annex A	ISO 8655-7:2022	ISO 8655-6:2022