


METTLER TOLEDO

Food Industry Solutions

Standard Connectivity For

- Precision Balances
- Bench Scales
- Digital Calipers
- pH Meters
- Moisture Analyzers
- And More



FreeWeigh.Net ^{compact}

Statistical Quality Control Systems

METTLER TOLEDO

Quality Management Starts Here Get Control. Get Compliant. Get Started.

- Compact single user weighing system with capabilities of a large control and statistics program
- Filling in compliance with US Handbook 133A and worldwide tolerance systems
- User dependent multi-language interface
- On the spot feedback and analysis reduces overfill and scrap
- Easy to create informative reports
- Standard SQL database
- Intuitive operator guidance and automatic data evaluation
- Expandable to FreeWeigh.Net multi-user integrated quality management system

Optional Software Modules

Part Number	Description
21900901	Batch Handling- provides statistics and reporting by batch
21900902	Attributes & Test Plans- provides monitoring of visual or non-measured quality data i.e: expiration date, security seal, color, taste
21900903	Adjustment Algorithms- provides machine adjustment messages
21900905	Device Integration- provides connection to calipers, moisture analyzers and other instruments
21900930	Free Reports- combine sample data, statistics, activities, alarms and graphs in fully configurable reports
21900931	SPC- control limits, CuSum, CPK
21900936	21CFR Part 11- provides compliant e-Sig and audit trail functions
21901126	Errors in Interventions- provides configurable messages and response tracking



FreeWeigh.Net Compact with central PC and Viper Balance



FreeWeigh.Net Compact with ID30, Balance and Bar Code

"Get Started" Food Manufacturing Bundle

Part Number	Description
21900954	FWN Compact - base license (requires central station PC by others)
MTASC	Telephone support during start-up
MTASC	Annual software support

Note: Bundles are also available ala carte.

METTLER TOLEDO

USA and Canada
 1900 Polaris Parkway
 Columbus, Ohio 43240
 Tel: (800) 523-5123
 (614) 438-4511
 Fax: (614) 438-4900
 indinfo@mt.com

Internet:
www.mt.com

Specifications subject to change without notice.
 © 2005 Mettler-Toledo, Inc.
 METTLER TOLEDO® is a registered trademark of Mettler-Toledo, Inc.
 2.5M 0305
 SOFTWARE
 I05-S003434.0E