



METTLER TOLEDO

Pharmaceutical Solutions


Standard Connectivity For

- Precision Balances
- Tablet Feeders
- Digital Calipers
- Thickness Gauges
- Hardness Testers
- Moisture Analyzers

FreeWeigh.Net^{compact}

Statistical Quality Control Systems

METTLER TOLEDO



Quality Management Starts Here Get Control. Get Compliant. Get Started.

- Compact single-user weighing system with the capabilities of a large control and statistics program
- Precision filling of liquids, tablets or creams in accordance with statutory limits
- On-the-spot quality and final packaging control for narrower fill quantity limits and fewer tolerance violations
- Easy-to-create informative reports and statistical evaluations
- Standard SQL database
- Traceable process and batch documentation as required by GMP
- Intuitive operator guidance and automatic data evaluation
- 21 CFR Part 11 Compliance**
- Expandable to FreeWeigh.Net multi-user integrated quality management system

Optional Software Modules

Part Number	Description
21900902	Attributes & Test Plans- provides monitoring of visual or non-measured quality data i.e: chips, cracks, color, expiration date
21900903	Adjustment Algorithms- provides machine adjustment messages
21900905	Device Integration- provides connection to calipers, hardness testers and other instruments
21900930	Free Reports- combine sample data, statistics, activities, alarms and graphs in fully configurable reports
21900931	SPC- control limits, CuSum, CPK
21900950/51	Dr. Schleuniger AT4 tablet tester interface
21901126	Errors in Interventions- provides configurable messages and response tracking



FreeWeigh.Net Compact with central PC and Viper Balance



FreeWeigh.Net Compact with ID30, Balance and Tablet Feeder

"Get Started" Pharmaceutical Bundle

Part Number	Description
21900954	FWN Compact - base license (requires central station PC by others)
21900901	Batch Handling- provides statistics and reporting by batch
21900936	21CFR Part 11- provides compliant e-Sig and audit trail functions
21900937	Validation manual 1 (DQ)
2190038	Validation manual 2 (IQ/OQ)
MTASC	Telephone support during start-up
MTASC	Annual software support

Note: Bundles are also available ala carte.

METTLER TOLEDO

USA and Canada
1900 Polaris Parkway
Columbus, Ohio 43240
Tel: (800) 523-5123
(614) 438-4511
Fax: (614) 438-4900
indinfo@mt.com

Internet:
www.mt.com

Specifications subject to change without notice.
© 2005 Mettler-Toledo, Inc.
METTLER TOLEDO® is a registered trademark of Mettler-Toledo, Inc.
2.5M 0305
SOFTWARE
I05-S003433.OE